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
- Li L, Xu Z, Zhou Y, Li T, Sun L, Chen H, Zhou R.. Analysis on Actinobacillus pleuropneumoniae LuxS regulated genes reveals pleiotropic roles of LuxS/AI-2 on biofilm formation, adhesion ability and iron metabolism. Microb Pathog. 2011 Jun;50(6):293-302.
- Zhou Y, Li L, Chen Z, Yuan H, Chen H, Zhou R. Adhesion protein ApfA of Actinobacillus pleuropneumoniae is required for pathogenesis and is a potential target for vaccine development. Clin Vaccine Immunol. 2013 Feb;20(2):287-94.
- Penghui Li Limin Li Wenhao Wang Weihong Jin Xiangmei Liua, Kelvin W.K. Yeung Paul K. Chu.Enhanced corrosion resistance and hemocompatibility of biomedical NiTi alloy by atmospheric-pressure plasma polymerized fluorine-rich coating.Applied Surface Science .2014 Apr;297(1):109 - 15
- Y Cao, Y Ma, M Zhang, H Wang, X Tu. Ultrasmall Graphene Oxide Supported Gold Nanoparticles as Adjuvants Improve Humoral and Cellular Immunity in Mice. Advanced Functional Materials . 2014 Nov;24(44):6963-71.
- Yang S, Shan T, Zhou Y, Jiang Y, Tong W, Liu F, Wen F, Zhang Q, Tong G.Molecular cloning and characterizations of porcine SAMHD1 and its roles in replication of highly pathogenic porcine reproductive and respiratory syndrome virus.Dev Comp Immunol . 2014 Dec;47(2):234-46
- Han Y, Sun S, Zhao M, Zhang Z, Gong S, Gao P, Liu J, Zhou J, Ma D, Gao Q, Wu P.CYC1 Predicts Poor Prognosis in Patients with Breast Cancer.Dis Markers . 2016;2016:3528064
- Zhang G, Xu M, Song Y, Su Z, Zhang H, Zhang C.TNF- α produced by SEC2 mutant (SAM-3)-activated human T cells induces apoptosis of HepG2 cells.APPL MICROBIOL BIOT .2016 Mar;100(6):2677 - 84.
- Tang HN, Xia Y, Xu J, Tian BM, Zhang XY, Chen FM.Assessment of cellular materials generated by co-cultured 'inflamed' and healthy periodontal ligament stem cells from patient-matched groups.Exp Cell Res . 2016 Aug 1;346(1):119-29
- Chen FY, Zhou LF, Li XY, Zhao JW, Xu SF, Huang WH, Gao LJ, Hao SJ, Ye YP, Sun HX.Stephanthraniline A suppressed CD4(+) T cell-mediated immunological hepatitis through impairing PKC θ function.Eur J Pharmacol . 2016 Oct 15;789:370-84
- Zhang G, Xu M, Zhang H, Song Y, Wang J, Zhang C.Up-regulation of granzyme B and perforin by staphylococcal enterotoxin C2 mutant induces enhanced cytotoxicity in Hepa1-6 cells.TOXICOL APPL PHARM . 2016 Dec 15;313:1-9
- Yan-Wei Yang, Lei Yang, Chao Zhang, Cai-Yun Gao, Ting Ma and Ling-Yi Kong.Physagulide Q suppresses proliferation and induces apoptosis in human hepatocellular carcinoma cells by regulating the ROS-JAK2/Src-STAT3 signaling pathway.RSC Adv . 2017; 7: 12793-12804
- Zheng T, Pu J, Chen Y, Guo Z, Pan H, Zhang L, Zhang H, Sun B, Zhang B.Exosomes Secreted from HEK293-APP Swe/Ind Cells Impair the Hippocampal Neurogenesis.Neurotox Res . 2017 Jul;32(1):82-93
- Fu X, Jiang B, Zheng B, Yan Y, Wang J, Duan Y, Li S, Yan L, Wang H, Chen B, Sang X, Ji W, Xu RH, Si W.Heterogenic transplantation of bone marrow-derived rhesus macaque mesenchymal stem cells ameliorates liver fibrosis induced by carbon tetrachloride in mouse.PeerJ . 2018 Feb 12;6:e4336

14. Yuan X,Han L,Fu P,Zeng H,Lv C,Chang W,Runyon RS,Ishii M,Han L,Liu K,Fan T,Zhang W,Liu R.Cinnamaldehyde accelerates wound healing by promoting angiogenesis via up-regulation of PI3K and MAPK signaling pathways.Lab Invest . 2018 Jun;98(6):783-798
15. Zhao H,Chen Y,Shao L,Xie M,Nie J,Qiu J,Zhao P,Ramezani H,Fu J,Ouyang H,He Y.Airflow-Assisted 3D Bioprinting of Human Heterogeneous Microspheroidal Organoids with Microfluidic Nozzle.Small . 2018 Sep;14(39):e1802630
16. Song H,Guo T,Zhao Z,Wei Y,Luo H,Weng W,Zhang R,Zhong M,Chen C,Su J,Shen W.Biocompatible PEGylated Gold nanorods function As cytokinesis inhibitors to suppress angiogenesis.Biomaterials . 2018 Sep;178:23-35
17. Yao Y,Yu H,Liu Y,Xu Q,Li X,Meng X,Huang C,Li J.PSTPIP2 Inhibits the Inflammatory Response and Proliferation of Fibroblast-Like Synoviocytes in vitro.Front Pharmacol . 2018 Dec 4;9:1432
18. Wei W,Pang S,Fu X,Tan S,Wang Q,Wang S,Sun D.The role of PERK and IRE1 signaling pathways in excessive fluoride mediated impairment of lymphocytes in rats' spleen in vivo and in vitro.Chemosphere. 2019 May;223:1-11
19. Wang K,Chen B,Yin T,Zhan Y,Lu Y,Zhang Y,Chen J,Wu W,Zhou S,Mao W,Tan Y,Du B,Liu X,Ho HI,Xiao J.N-Methylparoxetine Blocked Autophagic Flux and Induced Apoptosis by Activating ROS-MAPK Pathway in Non-Small Cell Lung Cancer Cells.Int J Mol Sci. 2019 Jul 11;20(14). pii: E3415
20. Damaneh MS,Hu JP,Huan XJ,Song SS,Tian CQ,Chen DQ,Meng T,Chen YL,Shen JK,Xiong B,Miao ZH,Wang YQ.A new BET inhibitor, 171, inhibits tumor growth through cell proliferation inhibition more than apoptosis induction.INVEST NEW DRUG. 2019 Jul 3
21. Li Y,Ruan X,Xu X,Li C,Qiang T,Zhou H,Gao J,Wang X.Shengmai Injection Suppresses Angiotensin II-Induced Cardiomyocyte Hypertrophy and Apoptosis via Activation of the AMPK Signaling Pathway Through Energy-Dependent Mechanisms.Front Pharmacol. 2019 Sep 20;10:1095
22. Lu JJ,Yang WM,Li F,Zhu W,Chen Z.Tunneling Nanotubes Mediated microRNA-155 Intercellular Transportation Promotes Bladder Cancer Cells' Invasive and Proliferative Capacity.INT J NANOMED. 2019 Dec 10;14:9731-9743
23. Kun Wang,Qing Gong,Yujuan Zhan,Bonan Chen,Ting Yin,Yuhua Lu,Yilin Zhang,Huiqi Wang,Junzi Ke,Biaoyan Du,Xiaodong Liu,Jianyong Xiao.Blockage of Autophagic Flux and Induction of Mitochondria Fragmentation by Paroxetine Hydrochloride in Lung Cancer Cells Promotes Apoptosis via the ROS-MAPK Pathway.Front Cell Dev Biol. 2020 Jan 22;7:397.;doi: 10.3389/fcell.2019.00397

Version 2016.12.26