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1. Min-Xia Su, Yu-Lian Xu, Xiao-Ming Jiang, Mu-Yang Huang, Le-Le Zhang, Luo-Wei Yuan, Xiao-Huang Xu, Qi Zhu, Jian-Li Gao, Jia-Hong Lu, Xiuping Chen, Ming-Qing Huang, Yitao Wang, Jin-Jian Lu . c-MYC-mediated TRIB3/P62+ aggresomes accumulation triggers paraptosis upon the combination of everolimus and ginsenoside Rh2 Acta Pharm Sin B. 2022 Mar;12(3):1240-1253.
2. Peng Li, Runjie Song, Yaqi Du, Huijiao Liu, Xiangdong Li . Adtrp regulates thermogenic activity of adipose tissue via mediating the secretion of S100b Cell Mol Life Sci. 2022 Jul 8;79(8):407.
3. Yongjin Sun, Xu Li, Xiaoxu Yang, Bi Chen, Wenzhi Zhang . Small Extracellular Vesicles Derived from Adipocytes Attenuate Intervertebral Disc Degeneration in Rats by Rejuvenating Senescent Nucleus Pulposus Cells and Endplate Cells by Delivering Exogenous NAMP Oxid Med Cell Longev. 2021 Aug 14;2021:9955448. doi: 10.1155/2021/9955448.
4. Yongbin Ma, Dan Zhou, Huanyan Zhang, Liming Tang, Fen Qian, Jianhua Su . Human Umbilical Cord Mesenchymal Stem Cell-Derived Extracellular Vesicles Promote the Proliferation of Schwann Cells by Regulating the PI3K/AKT Signaling Pathway via Transferring miR-21 Stem Cells Int. 2021 Sep 11;2021:1496101. doi: 10.1155/2021/1496101.
5. Yanpeng Ma, Yong Zhang, Chuan Qiu, Chunhui He, Ting He, Shuang Shi, Zhongwei Liu . Rivaroxaban Suppresses Atherosclerosis by Inhibiting FXa-Induced Macrophage M1 Polarization-Mediated Phenotypic Conversion of Vascular Smooth Muscle Cells Front Cardiovasc Med. 2021 Nov 15;8:739212. doi: 10.3389/fcvm.2021.739212.
6. Guidan Wang, Jine Dong . Network pharmacology approach to evaluate the therapeutic effects of mulberry leaf components for obesity Exp Ther Med. 2022 Jan;23(1):56.
7. Xiaoming Song, Lulu Han, Xiaowan Lin, Minghui Tian, Fenglei Sun, Bo Feng . Jian Pi Tiao Gan Yin alleviates obesity phenotypes through

mTORC1/SREBP1 signaling in vitro and in vivo *Ann Transl Med.* 2022 Mar;10(6):291.

8. Jiejie Wu, Ping Chen, Linjie Ju, Renhao Gao, Silu Li, Ziqian Huang, Yiqiu Cheng, Shuqi Gui, Zhixia Qiu, Jun Cheng, Fang Huang . *Corydalis saxicola* Bunting Total Alkaloids ameliorate diet-induced non-alcoholic steatohepatitis by regulating hepatic PI3K/Akt and TLR4/NF- $\kappa$ B pathways in mice *Biomed Pharmacother.* 2022 Jul;151:113132. doi: 10.1016/j.biopha.2022.113132.
9. Yougui Yang, Xin Ding, Fuzhong Chen, Xiaomin Wu, Yuying Chen, Qiang Zhang, Jun Cao, Junhong Wang, Yang Dai . Inhibition Effects of *Nippostrongylus brasiliensis* and Its Derivatives against Atherosclerosis in ApoE<sup>-/-</sup> Mice through Anti-Inflammatory Response *Pathogens.* 2022 Oct 20;11(10):1208.
10. Qiankun Luo, Jiayin Liu, Qiang Fu, Xu Zhang, Pengfei Yu, Pan Liu, Jiali Zhang, Huiyuan Tian, Song Chen, Hongwei Zhang, Tao Qin . Identifying cancer cell-secreted proteins that activate cancer-associated fibroblasts as prognostic factors for patients with pancreatic cancer *J Cell Mol Med.* 2022 Nov;26(22):5657-5669.
11. Jingchen Xu, Yao Lin, Mi Tian, Xinyi Li, Yuanyuan Yin, Qiming Li, Ziyu Li, Jialiang Zhou, Xiaoge Jiang, Yulin Li, Song Chen . Periodontal Ligament Stem Cell-Derived Extracellular Vesicles Enhance Tension-Induced Osteogenesis *ACS Biomater Sci Eng.* 2023 Jan 9;9(1):388-398.
12. Yiqi Yang, Kai Yuan, Yihao Liu, Qishan Wang, Yixuan Lin, Shengbing Yang, Kai Huang, Tianyou Kan, Yuxin Zhang, Mingming Xu, Zhifeng Yu, Qiming Fan, Yugang Wang, Hanjun Li, Tingting Tang . Constitutively activated AMPK $\alpha$ 1 protects against skeletal aging in mice by promoting bone-derived IGF-1 secretion *Cell Prolif.* 2023 Oct;56(10):e13476. doi: 10.1111/cpr.13476.

Version 2022.07.26