



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

青霉素-链霉素溶液(100X)

产品编号	产品名称	包装
C0222	青霉素-链霉素溶液(100X)	100ml

产品简介:

- 碧云天生产的青霉素-链霉素溶液(100X) (Penicillin-Streptomycin Solution, 100X)是专门用于细胞培养的双抗，经过过滤除菌，可以直接添加到细胞培养液内。
- 青霉素-链霉素溶液(100X)中，青霉素的含量为10kU/ml，链霉素的含量为10mg/ml。该溶液用0.9%氯化钠配制。在细胞培养液中推荐的青霉素的工作浓度为100U/ml，链霉素的工作浓度为0.1mg/ml。即按照100倍稀释使用即可。
- 一个包装即100ml青霉素-链霉素溶液(100X)可以配制10L细胞培养液。

包装清单:

产品编号	产品名称	包装
C0222	青霉素-链霉素溶液(100X)	100ml
—	说明书	1份

保存条件:

-20°C保存，一年有效。

注意事项:

- 尽量减少反复冻融的次数。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明:

青霉素-链霉素溶液(100X)可以参考如下两种方法之一使用：

1. 在无菌的细胞培养液中直接添加青霉素-链霉素溶液(100X):

按照每500ml细胞培养液添加5ml的比例加入青霉素-链霉素溶液(100X)，混匀即可使用。

2. 配制细胞培养液时加入青霉素-链霉素溶液(100X)，然后再过滤除菌:

配制细胞培养液时按照每配制1L细胞培养液加入10ml的比例加入青霉素-链霉素溶液(100X)，配制完成后过滤除菌即可使用。

使用本产品的文献:

- Wang PH, Gu ZH, Wan DH, Zhang MY, Weng SP, Yu XQ, He JG. The shrimp NF-κB pathway is activated by white spot syndrome virus (WSSV) 449 to facilitate the expression of WSSV069 (ie1), WSSV303 and WSSV371. PLoS One. 2011;6(9):e24773.
- Xia M, Zhu Y. Signaling pathways of ATP-induced PGE2 release in spinal cord astrocytes are EGFRtransactivation-dependent. Glia. 2011 Apr;59(4):664-74.
- Liu XY, Wei W, Wang CL, Yue H, Ma D, Zhu C, Ma GH, Du YG. Apoferritin-camouflaged Pt nanoparticles: surface effects on cellular uptake and cytotoxicity. J Mater Chem. 2011 Apr;20:7105-10.
- He YH, Song Y, Liao XL, Wang L, Li G, Alima, Li Y, Sun CH. The calcium-sensing receptor affects fat accumulation via effects on antilipolytic pathways in adipose tissue of rats fed low-calcium diets. J Nutr. 2011 Nov;141(11):1938-46.
- Wang MN, Liu XY, Cao CB, Shi C. Synthesis of band-gap tunable Cu-In-S ternary nanocrystals in aqueous solution. RSC Adv.2012 Feb;7:2666-70.
- Fan H, Yang L, Fu F, Xu H, Meng Q, Zhu H, Teng L, Yang M, Zhang L, Zhang Z, Liu K. Cardioprotective effects of salvianolic Acid a on myocardial ischemia-reperfusion injury in vivo and in vitro. Evid Based Complement Alternat Med. 2012;2012:508938.
- Zhang J, Tang L, Shen L, Zhou S, Duan Z, Xiao L, Cao Y, Mu X, Zha L, Wang H. High level of WAVE1 expression is associated with tumor aggressiveness and unfavorable prognosis of epithelial ovarian cancer. Gynecol Oncol. 2012 Oct;127(1):223-30.
- Wang M, Liu X, Cao C, Wang L. Highly luminescent CuInS₂-ZnS nanocrystals: achieving phase transfer and nuclear homing property simultaneously through simple TTAB modification. J. Mater. Chem. 2012 Aug;41:21979-86.
- Jiang M, Gan L, Zhu C, Dong Y, Liu J, Gan Y. Cationic core-shell liponanoparticles for ocular gene delivery. Biomaterials. 2012 Oct; 33(30):7621-30.
- Chai YY, Wang F, Li YL, Liu K, Xu H. Antioxidant Activities of Stilbenoids from Rheum emodi Wall. Evid Based Complement Alternat Med. 2012;2012:603678.
- Wang T, Jiang X. Size-dependent stability of water-solubilized CdTe quantum dots and their uptake mechanism by live HeLa cells. ACS Appl Mater Interfaces. 2013 Feb;5(4):1190-6.
- Okai I, Wang L, Gong L, Arko-Boham B, Hao L, Zhou X, Qi X, Hu J, Shao S. Overexpression of JAKMIP1 associates with Wnt/beta-catenin pathway activation and promotes cancer cellproliferation in vitro. Biomed Pharmacother. 2013 Apr;67(3):228-34.
- Jiang XY, Lu DB, Jiang YZ, Zhou LN, Cheng LQ, Chen B. PGC-1α prevents apoptosis in adipose-derived stem cells by reducing reactive oxygen species production in adiabatic microenvironment. Diabetes Res Clin Pract. 2013 Jun;100(3):368-75.
- Guo S, Sun X, Cheng J, Xu H, Dan J, Shen J, Zhou Q, Zhang Y, Meng L, Cao W, Tian Y. Apoptosis of THP-1 macrophages induced by protoporphyrin IX-mediated sonodynamic therapy. Int J Nanomedicine. 2013;8:2239-46.

15. Zhou DH, Wang X, Yang M, Shi X, Huang W, Feng Q. Combination of Low Concentration of (-)-Epigallocatechin Gallate (EGCG) and Curcumin Strongly Suppresses the Growth of Non-Small Cell Lung Cancer in Vitro and in Vivo through Causing Cell Cycle Arrest. *Int J Mol Sci.* 2013 Jun 5;14(6):12023-36.
16. Wang R, Li Z, Guo H, Shi W, Xin Y, Chang W, Huang T. Caveolin 1 knockdown inhibits the proliferation, migration and invasion of human breast cancer BT474 cells. *Mol Med Rep.* 2014 May;9(5):1723-8.
17. Ding H, Sun Y, Hou Y, Li L. Effects of INPP4B gene transfection combined with PARP inhibitor on castration therapy-resistant prostate cancer cell line, PC3. *Urol Oncol.* 2014 Jul;32(5):720-6.
18. Li M, Sun L, Li Y, Xie C, Wan D, Luo Y. Oxygen glucose deprivation/reperfusion astrocytes promotes primary neural stem/progenitor cell proliferation by releasing high-mobility group box 1. *Neurochem Res.* 2014 Aug;39(8):1440-50.
19. Chen X, Du Z, Shi W, Wang C, Yang Y, Wang F, Yao Y, He K, Hao A. 2-Bromopalmitate modulates neuronal differentiation through the regulation of histone acetylation. *Stem Cell Res.* 2014 Mar;12(2):481-91.
20. Miao CG, Qin D, Du CL, Ye H, Shi WJ, Xiong YY, Zhang XL, Yu H, Dou JF, Ma ST, Qin MS, Liu HZ, Fang YX, Zhou GL, Chen JZ, He X, Huang C, Huang Y, Zhang B, Song TW, Li J. DNMT1 activates the canonical Wnt signaling in rheumatoid arthritis model rats via a crucial functional crosstalk between miR-152 and the DNMT, MeCP2. *Int Immunopharmacol.* 2015 Sep;28(1):344-53.
21. Zhang C, Long X, Ding Y, Chen X, He J, Liu S, Geng Y, Wang Y, Liu X. Expression of DROSHA in the Uterus of Mice in Early Pregnancy and Its Potential Significance During Embryo Implantation. *Reprod Sci.* 2015 Jun;23(2):154-162.
22. Li C, Wen B, Wang L, Feng H, Xia X, Ding Z, Gao B, Zhang Y, Lan X. 99mTc-labeled single-domain antibody EG2 in targeting epidermal growth factor receptor: an in-vitro and mouse model in-vivo study. *Nucl Med Commun.* 2015 May;36(5):452-60.
23. Miao CG, Shi WJ, Xiong YY, Yu H, Zhang XL, Qin MS, Du CL, Song TW, Li J. miR-375 regulates the canonical Wnt pathway through FZD8 silencing in arthritis synovial fibroblasts. *Immunol Lett.* 2015 Mar; 164(1):1-10.
24. Wang X, Zhao C, Ji W, Xu Y, Guo H. Relationship of TLR, TLR4 and tissue remodeling in chronic rhinosinusitis. *Int J Clin Exp Pathol.* 2015 Feb 1;8(2):1199-212.
25. Zhou X, Pei J, Huang G. Reactive paper spray mass spectrometry for in situ identification of quinones. *Rapid Commun Mass Spectrom.* 2015 Jan 15;29(1):100-6.
26. Cai Y, Cao C, He X, Yang C, Tian L, Zhu R, Pan Y. Enhanced magnetic resonance imaging and staining of cancer cells using ferrimagnetic H-ferritin nanoparticles with increasing core size. *Int J Nanomedicine.* 2015 Apr 1;10:2619-34.
27. Zhou L, Lai Y, Huang W, Huang S, Xu Z, Chen J, Wu D. Biofunctionalization of microgroove titanium surfaces with antimicrobial peptide to enhance their bactericidal activity and cytocompatibility. *Colloids Surf B Biointerfaces.* 2015 Apr 1;128:552-60.
28. Xiang Q, Ma Y, Dong J, Shen R. Carnosic acid induces apoptosis associated with mitochondrial dysfunction and Akt inactivation in HepG2 cells. *Int J Food Sci Nutr.* 2015 Feb;66(1):76-84.
29. Wu C, Zhang G, Xia T, Li Z, Zhao K, Deng Z, Guo D, Peng B. Bioinspired synthesis of polydopamine/Ag nanocomposite particles with antibacterial activities. *Mater Sci Eng C Mater Biol Appl.* 2015 Oct; 55:155-65.
30. Yali Zenga, Jianjun Denga, Tingting Huob, Faqin Dongb, Limin Wang. Assessment of genetic toxicity with major inhalable mineral granules in A549 cells. *Applied Clay Science* 2016;119:175-182.
31. Zhao N, Yang Z, Li B, Meng J, Shi Z, Li P, Fu S. RGD-conjugated mesoporous silica-encapsulated gold nanorods enhance the sensitization of triple-negative breast cancer to megavoltage radiation therapy. *Int J Nanomedicine.* 2016 Oct 27;11:5595-5610.
32. Zhang C, Long X, Ding Y, Chen X, He J, Liu S, Geng Y, Wang Y, Liu X. Expression of DROSHA in the Uterus of Mice in Early Pregnancy and Its Potential Significance During Embryo Implantation. *Reprod Sci.* 2016 Feb;23(2):154-62.
33. Liu N, Liu X, Yu Q, Chen X, Ding Y, He J, Gao R, Wang Y, Liu X. SP0P Regulates Endometrial Stromal Cell Decidualization in Mice. *Reprod Sci.* 2016 Nov;23(11):1565-1574.
34. Shen Q, Ma X, Hua Y, Chen M, Wang Y, Zhou Q, Ye W, Zhu X. Aquaporin 3 Expression Induced by Salvia Miltiorrhiza via ERK1/2 Signal Pathway in the Primary Human Amnion Epithelium Cells from Isolated Oligohydramnios. *Curr Mol Med.* 2016;16(3):312-9.
35. Zhang C, Wu JM, Liao M, Wang JL, Xu CJ. The ROCK/GGTase Pathway Are Essential to the Proliferation and Differentiation of Neural Stem Cells Mediated by Simvastatin. *J Mol Neurosci.* 2016 Dec;60(4):474-485.
36. Xu H, Li Z, Yu Y, Sizdahkhani S, Ho WS, Yin F, Wang L, Zhu G, Zhang M, Jiang L, Zhuang Z, Qin J. A dynamic in vivo-like organotypic blood-brain barrier model to probe metastatic brain tumors. *Sci Rep.* 2016 Nov 10;6:36670.
37. Hou L, Gong C, Zhu Y. In vitro construction and in vivo regeneration of esophageal bilamellar muscle tissue. *J Biomater Appl.* 2016 Apr;30(9):1373-84.
38. Ji MM, Liu AQ, Sima YH, Xu SQ. 20-Hydroxyecdysone stimulates nuclear accumulation of BmNep1, a nuclear ribosome biogenesis-related protein in the silkworm, *Bombyx mori*. *Insect Mol Biol.* 2016 Oct;25(5):617-28.
39. Duan P, Hu C, Quan C, Yu T, Huang W, Chen W, Tang S, Shi Y, Martin FL, Yang K. 4-Nonylphenol induces autophagy and attenuates mTOR-p70S6K/4EBP1 signaling by modulating AMPKactivation in Sertoli cells. *Toxicol Lett.* 2017 Feb 5;267:21-31.
40. Liu Y, Sun Z, Xu D, Liu J, Li X, Wu X, Zhang Y, Wang Q, Huang C, Meng X, Li J. Hesperidin derivative-11 inhibits fibroblast-like synoviocytes proliferation by activating Secreted frizzled-related protein 2 in adjuvant arthritis rats. *Eur J Pharmacol.* 2017 Jan 5;794:173-183.
41. Liu Y, Sun Z, Xu D, Liu J, Li X, Wu X, Zhang Y, Wang Q, Huang C, Meng X, Li J. Hesperidin derivative-11 inhibits fibroblast-like synoviocytes proliferation by activating Secreted frizzled-related protein 2 in adjuvant arthritis rats. *Eur J Pharmacol.* 2017 Jan 5;794:173-183.
42. Liu Y, Sun Z, Xu D, Liu J, Li X, Wu X, Zhang Y, Wang Q, Huang C, Meng X, Li J. Hesperidin derivative-11 inhibits fibroblast-like synoviocytes proliferation by activating Secreted frizzled-related protein 2 in adjuvant arthritis rats. *Eur J Pharmacol.* 2017 Jan 5;794:173-183.
43. Duan P, Hu C, Quan C, Yu T, Huang W, Chen W, Tang S, Shi Y, Martin FL, Yang K. 4-Nonylphenol induces autophagy and attenuates mTOR-p70S6K/4EBP1 signaling by modulating AMPKactivation in Sertoli cells. *Toxicol Lett.* 2017 Feb 5;267:21-31.
44. Cao Z, Liu C, Bai Y, Dou C, Li JM, Shi DW, Dong SW, Xiang Q. Inhibitory effect of dihydroartemisinin on chondrogenic and hypertrophic differentiation of mesenchymal stem cells. *Am J Transl Res.* 2017 Jun 15;9(6):2748-2759.
45. Gong Y, Ji Y, Liu F, Li J, Cao Y. Cytotoxicity, oxidative stress and inflammation induced by ZnO nanoparticles in endothelial cells: interaction with palmitate or lipopolysaccharide. *J Appl Toxicol.* 2017 Aug;37(8):895-901.
46. He F, Du T, Jiang Q, Zhang Y. Synergistic Effect of Notch-3-Specific Inhibition and Paclitaxel in Non-Small Cell Lung Cancer (NSCLC) Cells Via Activation of The Intrinsic Apoptosis Pathway. *MED SCI MONITOR.* 2017 Aug 3;23:3760-3769.
47. Wang J, Wang S, Zhang W, Wang X, Liu X, Liu L, Li L, Liang Y, Yu J, Jeong LS, Jia L, Zhao H, Zhang Y. Targeting neddylation pathway with MLN4924 (Pevonedistat) induces NOXA-dependent apoptosis in renal cell carcinoma. *BIOCHEM BIOPH RES CO.* 2017 Sep 2;490(4):1183-1188.
48. Cao Z, Bai Y, Liu C, Dou C, Li J, Xiang J, Zhao C, Xie Z, Xiang Q, Dong S. Hypertrophic differentiation of mesenchymal stem cells is suppressed by xanthotoxin via the p38 MAPK/HDAC4 pathway. *Mol Med Rep.* 2017 Sep;16(3):2740-2746.
49. Sun ZH, Liu YH, Liu JD, Xu DD, Li XF, Meng XM, Ma TT, Huang C, Li J. MeCP2 Regulates PTCH1 Expression Through DNA Methylation in Rheumatoid Arthritis. *Inflammation.* 2017 Oct;40(5):1497-1508.
50. Yuan J, Han R, Esther A, Wu Q, Yang J, Yan W, Ji X, Liu Y, Li Y, Yao W, Ni C. Polymorphisms in autophagy related genes and the coal workers' pneumoconiosis in a Chinese population. *Gene.* 2017 Oct 20;632:36-42.
51. Lu CJ, Jiang XF, Junaid M, Ma YB, Jia PP, Wang HB, Pei DS. Graphene oxide nanosheets induce DNA damage and activate the base excision repair (BER) signaling pathway both in vitro and in vivo. *Chemosphere.* 2017 Oct;184:795-805.
52. Yang Y, Wu M, Geng Y, Liu X, Yang Y, Chen X, Ding Y, He J, Wang Y, Xie L. Roles of DEK in the endometrium of mice in early pregnancy. *Gene.* 2017 Nov 3. pii: S0378-1119(17)30939-3.
53. Liu YR, Yan X, Yu HX, Yao Y, Wang JQ, Li XF, Chen RN, Xu QQ, Ma TT, Huang C, Li J. NLRC5 promotes cell proliferation via regulating the

- NF- κ B signaling pathway in Rheumatoidarthritis.*Mol Immunol.* 2017 Nov;91:24-34
54. Cao Z, Huang S, Li J, Bai Y, Dou C, Liu C, Kang F, Gong X, Ding H, Hou T, Dong S. Long noncoding RNA expression profiles in chondrogenic and hypertrophic differentiation of mouse mesenchymal stem cells.*FUNCT INTEGR GENOMIC*. 2017 Nov;17(6):739-749
55. He T, Long J, Li J, Liu L, Cao Y. Toxicity of ZnO nanoparticles (NPs) to A549 cells and A549 epithelium in vitro: Interactions with dipalmitoyl phosphatidylcholine (DPPC).*ENVIRON TOXICOL CHEM*. 2017 Dec;36:233-240
56. Fang X, Jiang L, Gong Y, Li J, Liu L, Cao Y. The presence of oleate stabilized ZnO nanoparticles (NPs) and reduced the toxicity of aged NPs to Caco-2 and HepG2 cells.*CHEM-BIOL INTERACT*. 2017 Dec 25;278:40-47
57. Zhang MC, Zhao SG, Wang SS, Luo CC, Gao HN, Zheng N, Wang JQ. d-Glucose and amino acid deficiency inhibits casein synthesis through JAK2/STAT5 and AMPK/mTOR signaling pathways in mammary epithelial cells of dairy cows.*J Dairy Sci*. 2018 Feb;101(2):1737-1746
58. Ni Y, Gu WW, Liu ZH, Zhu YM, Rong JG, Kent TA, Li M, Qiao SG, An JZ, Zhang HL. RIP1K Contributes to Neuronal and Astrocytic Cell Death in Ischemic Stroke via Activating Autophagic-lysosomal Pathway.*Neuroscience*. 2018 Feb 10;371:60-74
59. Yang Y, Wu M, Geng Y, Liu X, Yang Y, Chen X, Ding Y, He J, Wang Y, Xie L. Roles of DEK in the endometrium of mice in early pregnancy.*Gene*. 2018 Feb 5;642:261-267
60. Cao Z, Huang S, Dou C, Xiang Q, Dong S. Cyanidin suppresses autophagic activity regulating chondrocyte hypertrophic differentiation.*J Cell Physiol*. 2018 Mar;233(3):2332-2342
61. Liu F, Pan Z, Zhang J, Ni J, Wang C, Wang Z, Gu F, Dong W, Zhou W, Liu H. Overexpression of RHEB is associated with metastasis and poor prognosis in hepatocellular carcinoma.*Oncol Lett*. 2018 Mar;15(3):3838-3845
62. Chen H, Zhou X, Han TL, Baker PN, Qi H, Zhang H. Decreased IL-33 Production Contributes to Trophoblast Cell Dysfunction in Pregnancies with Preeclampsia.*MEDIAT INFLAMM*. 2018 Mar 15;2018:9787239
63. Jin X, Liu L, Zhang Y, Xiang Y, Yin G, Lu Y, Shi L, Dong J, Shen C. Advanced Glycation End Products Enhance Murine Monocyte Proliferation in Bone Marrow and Prime Them into an Inflammatory Phenotype through MAPK Signaling.*J7 Diabetes Res*. 2018 Mar 22;2018:2527406.
64. Gao R, Ma J, Wen Z, Yang P, Zhao J, Xue M, Chen Y, Aldarouish M, Hu HM, Zhu XJ, Pan N, Wang LX. Tumor cell-released autophagosomes (TRAP) enhance apoptosis and immunosuppressive functions of neutrophils.*Oncoimmunology*. 2018 Mar 6;7(6):e1438108
65. Zhang J, Zou Z, Wang B, Xu G, Wu Q, Zhang Y, Yuan Z, Yang X, Yu C. Lysosomal deposition of copper oxide nanoparticles triggers HUVEC cells death.*Biomaterials*. 2018 Apr;161:228-239
66. Deng X, Dai P, Yu M, Chen L, Zhu C, You X, Li L, Zeng Y. Cyclophilin A is the potential receptor of the Mycoplasma genitalium adhesion protein.*Int J Med Microbiol*. 2018 Apr;308(3):405-412
67. Zhang X, Wang X, Liu T, Mo M, Ao L, Liu J, Cao J, Cui Z. ZnSO₄ rescued vimentin from collapse in DBP-exposed Sertoli cells by attenuating ER stress and apoptosis.*Toxicol In Vitro*. 2018 Apr;48:195-204
68. Zhou D, Liu W, Liang S, Sun B, Liu A, Cui Z, Han X, Yuan L. Apoptin-derived peptide reverses cisplatin resistance in gastric cancer through the PI3K-AKT signaling pathway.*CANCER MED-US*. 2018 Apr;7(4):1369-1383
69. Guo HX, Ye N, Yan P, Qiu MY, Zhang J, Shen ZG, He HY, Tian ZQ, Li HL, Li JT. Sodium chloride exacerbates dextran sulfate sodium-induced colitis by tuning proinflammatory and antiinflammatory lamina propria mononuclear cells through p38/MAPK pathway in mice.*WORLD J GASTROENTERO*. 2018 Apr 28;24(16):1779-1794
70. Chen Z, Luo S, Chen Y, Xie X, Du Z, Jiang L. High expression of NPRL2 is linked to poor prognosis in patients with prostate cancer.*Hum Pathol*. 2018 Jun;76:141-148
71. Miao C, Chang J, Zhang G, Yu H, Zhou L, Zhou G, Zhao C. CUL4B promotes the pathology of adjuvant-induced arthritis in rats through the canonical Wnt signaling.*J MOL MED*. 2018 Jun;96(6):495-511
72. Deng X, Zhu Y, Dai P, Yu M, Chen L, Zhu C, You X, Li L, Zeng Y. Three polypeptides screened from phage display random peptide library may be the receptor polypeptide of Mycoplasma genitalium adhesion protein.*MICROB PATHOGENESIS*. 2018 Jul;120:140-146
73. Wang X, Wu B, Zhong Z. Downregulation of YAP inhibits proliferation, invasion and increases cisplatin sensitivity in human hepatocellular carcinoma cells.*Oncol Lett*. 2018 Jul;16(1):585-593
74. Fang T, Cao R, Wang W, Ye H, Shen L, Li Z, Hu J, Gao Q. Alterations in necroptosis during ALDH2-mediated protection against high glucose-induced H9c2 cardiac cell injury.*Mol Med Rep*. 2018 Sep;18(3):2807-2815
75. Luo Y, Wu C, Liu L, Gong Y, Peng S, Xie Y, Cao Y. 3-Hydroxyflavone enhances the toxicity of ZnO nanoparticles in vitro.*J Appl Toxicol*. 2018 Sep;38(9):1206-1214
76. Yu G, Chen J, Li S, Pu P, Huang W, Zhao Y, Peng X, Wang R, Lei H. RING finger protein 10 prevents neointimal hyperplasia by promoting apoptosis in vitro and in vivo.*Life Sci*. 2018 Sep 1;208:325-332
77. Li Z, Yanfang W, Li J, Jiang P, Peng T, Chen K, Zhao X, Zhang Y, Zhen P, Zhu J, Li X. Tumor-released exosomal circular RNA PDE8A promotes invasive growth via the miR-338/MACC1/MET pathway in pancreatic cancer.*Cancer Lett*. 2018 Sep 28;432:237-250
78. Qin Y, He Y, Zhu YM, Li M, Ni Y, Liu J, Zhang HL. CID1067700, a late endosome GTPase Rab7 receptor antagonist, attenuates brain atrophy, improves neurologic deficits and inhibits reactive astrogliosis in rat ischemic stroke.*Acta Pharmacol Sin*. 2018 Oct 12
79. Li W, Liu J, Chen JR, Zhu YM, Gao X, Ni Y, Lin B, Li H, Qiao SG, Wang C, Zhang HL, Ao GZ. Neuroprotective Effects of DTIO, A Novel Analog of Nec-1, in Acute and Chronic Stages After Ischemic Stroke.*Neuroscience*. 2018 Oct 15;390:12-29
80. Meng Y, Zhang Y, Jia N, Qiao H, Zhu M, Meng Q, Lu Q, Zu Y. Synthesis and evaluation of a novel water-soluble high Se-enriched Astragalus polysaccharide nanoparticles.*Int J Biol Macromol*. 2018 Oct 15;118(Pt B):1438-1448
81. Dai WT, White RR, Liu JX, Liu HY. Seryl-tRNA synthetase-mediated essential amino acids regulate β -casein synthesis via cell proliferation and mammalian target of rapamycin (mTOR) signaling pathway in bovine mammary epithelial cells.*J Dairy Sci*. 2018 Nov;101(11):10456-10468
82. Junaid M, Jia PP, Tang YM, Xiong WX, Huang HY, Strauss PR, Li WG, Pei DS. Mechanistic toxicity of DEHP at environmentally relevant concentrations (ERCs) and ecological risk assessment in the Three Gorges Reservoir Area, China.*Environ Pollut*. 2018 Nov;242(Pt B):1939-1949
83. Lv Y, Li B, Han K, Xiao Y, Yu X, Ma Y, Jiao Z, Gao J. The Nedd8-activating enzyme inhibitor MLN4924 suppresses colon cancer cell growth via triggering autophagy.*KOREAN J PHYSIOL PHA*. 2018 Nov;22(6):617-625
84. Xu W, Wang L, Li J, Cai Y, Xue Y. TXNIP mediated the oxidative stress response in glomerular mesangial cells partially through AMPK pathway.*Biomed Pharmacother*. 2018 Nov;107:785-792
85. Wang JL, Wang JJ, Cai ZN, Xu CJ. The effect of curcumin on the differentiation, apoptosis and cell cycle of neural stem cells is mediated through inhibiting autophagy by the modulation of Atg7 and p62.*Int J Mol Med*. 2018 Nov;42(5):2481-2488
86. Yu D, Zhang L, Wang H, Chen F, Chen J, Zhang Z, Li J, Xing C, Li H, Li J, Cai Y. A potential role for SMAD9 in goose follicular selection through regulation of mRNA levels of luteinizing hormone receptor.*Theriogenology*. 2018 Nov 22. pii: S0093-691X(18)31081-1
87. Zhang J, Wang B, Wang H, He H, Wu Q, Qin X, Yang X, Chen L, Xu G, Yuan Z, Yi Q, Zou Z, Yu C. Disruption of the superoxide anions-mitophagy regulation axis mediates copper oxide nanoparticles-induced vascular endothelial cell death.*FREE RADICAL BIO MED*. 2018 Dec;129:268-278
88. Wang Y, Dai L, Wu H, Zhang ZR, Wang WY, Fu J, Deng R, Li F, Dai XJ, Zhan X. Novel anti-inflammatory target of geniposide: Inhibiting IgG 1/Ras-Erk1/2 signal pathway via the miRNA-124a in rheumatoid arthritis synovial fibroblasts.*Int Immunopharmacol*. 2018 Dec;65:284-294
89. Yao Y, Yu H, Liu Y, Xu Q, Li X, Meng X, Huang C, Li J. PSTPIP2 Inhibits the Inflammatory Response and Proliferation of Fibroblast-Like Synoviocytes in vitro.*Front Pharmacol*. 2018 Dec 4;9:1432
90. Zhao X, Feng P, He W, Du X, Chen C, Suo L, Liang M, Zhang N, Na A, Zhang Y. The Prevention and Inhibition Effect of Anthocyanins on Colorectal Cancer.*CURR PHARM DESIGN*. 2019;25(46):4919-4927. (IF 2.208)
91. Dong M, Ye T, Bi Y, Wang Q, Kuerban K, Li J, Feng M, Wang K, Chen Y, Ye L. A novel hybrid of 3-benzyl coumarin seco-B-ring derivative and phenylsulfonylfuroxan induces apoptosis and autophagy in

- non-small-cell lung cancer. *Phytomedicine*. 2019 Jan;52:79-88
92. Gan L,Pan S,Cui J,Bai J,Jiang P,He Y.Functional analysis of the correlation between ABCB11 gene mutation and primary intrahepatic stone.*Mol Med Rep.* 2019 Jan;19(1):195-204
 93. Liu T,Liang H,Liu L,Gong Y,Ding Y,Liao G,Cao Y.Influence of pristine and hydrophobic ZnO nanoparticles on cytotoxicity and endoplasmic reticulum (ER) stress-autophagy-apoptosis gene expression in A549-macrophage co-culture.*ECOTOX ENVIRON SAFE*. 2019 Jan 15;167:188-195. (IF 4.872)
 94. Tan M,Liu W,Liu F,Zhang W,Gao H,Cheng J,Chen Y,Wang Z,Cao Y,Ran H.Silk Fibroin-Coated Nanoagents for Acidic Lysosome Targeting by a Functional Preservation Strategy in Cancer Chemotherapy.*Theranostics*. 2019 Jan 25;9(4):961-973
 95. Gao HN,Guo HY,Zhang H,Xie XL,Wen PC,Ren FZ.Yak-milk-derived exosomes promote proliferation of intestinal epithelial cells in an hypoxic environment.*J Dairy Sci*. 2019 Feb;102(2):985-996.
 96. Yue E,Tuguzbaeva G,Chen X,Qin Y,Li A,Sun X,Dong C,Liu Y,Yu Y,Zahra SM,Shan Q,Jiang Y,Du Z,Bai Y.Anthocyanin is involved in the activation of pyroptosis in oral squamous cell carcinoma.*Phytomedicine*. 2019 Mar 15;56:286-294
 97. Jiang L,Xu D,Zhang WJ,Tang Y,Peng Y.Astrocytes induce proliferation of oligodendrocyte progenitor cells via connexin 47-mediated activation of Chi3l1 expression.*EUR REV MED PHARMACO*. 2019 Apr;23(7):3012-3020
 98. Li T,Chen L,Yu Y,Yang B,Li P,Tan XQ.Resveratrol alleviates hypoxia/reoxygenation injury-induced mitochondrial oxidative stress in cardiomyocytes.*Mol Med Rep.* 2019 Apr;19(4):2774-2780
 99. Jin Q,Zhou J,Xu X,Huang F,Xu W.Hypoxia-inducible factor-1 signaling pathway influences the sensitivity of HCC827 cells to gefitinib.*Oncol Lett*. 2019 Apr;17(4):4034-4043
 - 100.Long J,Yang CS,He JL,Liu XQ,Ding YB,Chen XM,Tong C,Peng C,Wang YX,Gao RF.FOXO3a is essential for murine endometrial decidualization through cell apoptosis during early pregnancy.*J Cell Physiol*. 2019 Apr;234(4):4154-4166
 - 101.Cui YL,Xue RQ,Xi He,Ming Zhao,Yu XJ,Liu LZ,Wu Q,Si Yang,Li DL,Zang WJ.Cholinergic drugs ameliorate endothelial dysfunction by decreasing O-GlcNAcylation via M3 AChR-AMPK-ER stress signaling.*Life Sci*. 2019 Apr 1;222:1-12
 - 102.Zhang Y,Guo J,Li Y,Jiao K,Zhang Y.let-7a suppresses liver fibrosis via TGF β /SMAD signaling transduction pathway.*Exp Ther Med*. 2019 May;17(5):3935-3942
 - 103.Zhu Y,Zhang Y,Huang X,Xie Y,Qu Y,Long H,Gu N,Jiang W.Z-Ligustilide protects vascular endothelial cells from oxidative stress and rescues high fat diet-induced atherosclerosis by activating multiple NRF2 downstream genes.*Atherosclerosis*. 2019 May;284:110-120
 - 104.Liu H,Zou M,Li P,Wang H,Lin X,Ye J.Oxymatrine-mediated maturation of dendritic cells leads to activation of FOXP3+/CD4+ Treg cells and reversal of cisplatin-resistance in lung cancer cells.*Mol Med Rep.* 2019 May;19(5):4081-4090
 - 105.Jiang S,Du Y,Liu D,He J,Huang Y,Qin K,Zhou X.Inhibitory Effect of Slit2-N on VEGF165-induced proliferation of vascular endothelia via Slit2-N-Rbo4-Akt pathway in choroidal neovascularization.*Cell Cycle*. 2019 Jun;18(11):1241-1253
 - 106.Yu Q,Pu SY,Wu H,Chen XQ,Jiang JJ,Gu KS,He YH,Kong Q.P.TICRR Contributes to Tumorigenesis Through Accelerating DNA Replication in Cancers.*Front Oncol*. 2019 Jun 18;9:516
 - 107.Luo F,Wei H,Guo H,Li Y,Feng Y,Bian Q,Wang Y.LncRNA MALAT1, an lncRNA acting via the miR-204/ZEB1 pathway, mediates the EMT induced by organic extract of PM2.5 in lung bronchial epithelial cells.*AM J PHYSIOL-LUNG C*. 2019 Jul 1;317(1):L87-L98
 - 108.Hou R,Wang X,Wei Q,Feng P,Mou X,Zhu Y,Shen Z.Biological properties of a bionic scaffold for esophageal tissue engineering research.*COLLOID SURFACE B*. 2019 Jul 1;179:208-217
 - 109.Xu J,Zhai Y,Feng L,Xie T,Yao W,Shan J,Zhang L.An optimized analytical method for cellular targeted quantification of primary metabolites in tricarboxylic acid cycle and glycolysis using gas chromatography-tandem mass spectrometry and its application in PHARMACEUT BIOMED. 2019 Jul 15;171:171-179
 - 110.Li M,Dong Y,Chen Z,Meng L,Liu X,Zhang X,Wang H,Mao W,Zhang J,Jiang Z,Huang T,Hu J,Luo P,Korner H,Ying S,Li J.MicroRNA-31 Negatively Regulates Interleukin-34 Expression In Vitro.*Immunol Invest*. 2019 Aug;48(6):597-607
 - 111.Diao H,Cheng N,Zhao Y,Xu H,Dong H,Thamm DH,Zhang D,Lin D.Ivermectin inhibits canine mammary tumor growth by regulating cell cycle progression and WNT signaling.*BMC Vet Res*. 2019 Aug 2;15(1):276
 - 112.Yang F,Chen F,Li L,Yan L,Badri T,Lv C,Yu D,Zhang M,Jang X,Li J,Yuan L,Wang G,Li H,Li J,Cai Y.Three Novel Players: PTK2B, SYK, and TNFRSF21 Were Identified to Be Involved in the Regulation of Bovine Mastitis Susceptibility via GWAS and Post-transcriptional Analysis.*Front Immunol*. 2019 Aug 6;10:1579
 - 113.Song B,Li XF,Yao Y,Xu QQ,Meng XM,Huang C,Li J.BMP9 inhibits the proliferation and migration of fibroblast-like synoviocytes in rheumatoid arthritis via the PI3K/AKT signaling pathway.*Int Immunopharmacol*. 2019 Sep;74:105685
 - 114.Xiong T,Zhang Z,Zheng R,Huang J,Guo L.N-acetyl cysteine inhibits lipopolysaccharide-induced apoptosis of human umbilical vein endothelial cells via the p38MAPK signaling pathway.*Mol Med Rep*. 2019 Sep;20(3):2945-2953
 - 115.Xiao S,Yu X,Zhang L,Zhang Y,Fan W,Sun T,Zhou C,Liu Y,Liu Y,Gong M,Zhang D.Synthesis Of PEG-Coated, Ultrasmall, Manganese-Doped Iron Oxide Nanoparticles With High Relaxivity For T1/T2 Dual-Contrast Magnetic Resonance Imaging.*INT J NANOMED*. 2019 Oct 24;14:8499-8507
 - 116.Xiang Y,Yao X,Wang X,Zhao H,Zou H,Wang L,Zhang QX.Houshiheisan promotes angiogenesis via HIF-1 α /VEGF and SDF-1/CXCR4 pathways: in vivo and in vitro.*BIOSCIENCE REP*. 2019 Oct 30;39(10). pii: BSR20191006
 - 117.Cao M,Gao D,Zhang N,Duan Y,Wang Y,Mujtaba H,Wang Y.Shp2 expression is upregulated in cervical cancer, and Shp2 contributes to cell growth and migration and reduces sensitivity to cisplatin in cervical cancer cells.*Pathol Res Pract*. 2019 Nov;215(11):152621
 - 118.Li J,Tong C,Xu P,Wang L,Han TL,Wen L,Luo X,Tan B,Zhu F,Gui S,Gao R,Qi H,Baker PN.QSOX1 regulates trophoblastic apoptosis in preeclampsia through hydrogen peroxide production.*J MATERN-FETAL NEO M*. 2019 Nov;32(22):3708-3715
 - 119.Wu L,Chen P,Ying J,Zhang Q,Liu F,Lv B,Che Z,Zhang W,Wu M,Zhang J,Yang D,Liu J.MAT2B mediates invasion and metastasis by regulating EGFR signaling pathway in hepatocellular carcinoma.*Clin Exp Med*. 2019 Nov;19(4):535-546
 - 120.Xu Q,Yin S,Yao Y,Li X,Song B,Yang Y,Liu Y,Chen R,Li J,Ma T,Meng X,Huang C,Li J.MAST3 modulates the inflammatory response and proliferation of fibroblast-like synoviocytes in rheumatoid arthritis.*Int Immunopharmacol*. 2019 Dec;77:105900
 - 121.Syeda MZ,Fasae MB,Yue E,Ishimwe AP,Jiang Y,Du Z,Yang B,Bai Y.Anthocyanidin attenuates myocardial ischemia induced injury via inhibition of ROS-JNK-Bcl-2 pathway: New mechanism of anthocyanidin action.*Phytother Res*. 2019 Dec;33(12):3129-3139
 - 122.Lin F,Yin HB,Li XY,Zhu GM,He WY,Gou X.Bladder cancer cell-secreted exosomal miR-21 activates the PI3K/AKT pathway in macrophages to promote cancer progression.*Int J Oncol*. 2020 Jan;56(1):151-164
 - 123.Meng L,Chen Z,Jiang Z,Huang T,Hu J,Luo P,Zhang H,Huang M,Huang L,Chen Y,Lu M,Xu AM,Ying S.MiR-122-5p suppresses the proliferation, migration, and invasion of gastric cancer cells by targeting LYN.*ACTA BIOCH BIOPH SIN*. 2020 Jan 2;52(1):49-57124.
 - 124.Daoju Gan,Weiyang He,Hubin Yin,Xin Gou β -elemene enhances cisplatin-induced apoptosis in bladder cancer cells through the ROS-AMPK signaling pathway.*Oncol Lett*. 2020 Jan;19(1):291-300.;doi: 10.3892/ol.2019.11103
 - 125.Jialin Jia,Lin Nie,Yi Liu.Butyrate alleviates inflammatory response and NF- κ B activation in human degenerated intervertebral disc tissues.*Int Immunopharmacol*. 2020 Jan;78:106004.;doi: 10.1016/j.intimp.2019.106004
 - 126.Gong M,Zhai X,Yu L,Li C,Ma X,Shen Q,Han Y,Yang DADSCs inhibit photoaging- and photocarcinogenesis-related inflammatory responses and extracellular matrix degradation.*J Cell Biochem*. 2020 Feb;121(2):1205-1215
 - 127.Fan G,Xu P,Tu P.MiR-1827 functions as a tumor suppressor in lung adenocarcinoma by targeting MYC and FAM83F.*J Cell Biochem*. 2020 Feb;121(2):1675-1689
 - 128.Yang L,Liu J,Shan Q,Geng G,Shao P.High glucose inhibits proliferation and differentiation of osteoblast in alveolar bone by inducing pyroptosis.*BIOCHEM BIOPH RES CO*. 2020 Feb 5;522(2):471-478
 - 129.Xiaohong Wang,Jinghong Chi,Di Huang,Li Ding,Xiaojing Zhao,Lai Jiang,Yang Yu,Feng Gao. α -synuclein promotes progression of

- Parkinson's disease by upregulating autophagy signaling pathway to activate NLRP3 inflammasomeExp Ther Med. 2020 Feb;19(2):931-938.;doi: 10.3892/etm.2019.8297
- 130.Huimin Lu,Tongguo Shi,Mingyuan Wang,Xiaomi Li,Yanzheng Gu,Xueguang Zhang,Guangbo Zhang,Weichang Chen.B7-H3 inhibits the IFN- γ -dependent cytotoxicity of V γ 9V δ 2 T cells against colon cancer cellsOncoimmunology. 2020 Apr 14;9(1):1748991.;doi: 10.1080/2162402X.2020.1748991
- 131.Jin Wang,Sam Toan,Ruibing Li,Hao Zhou.Melatonin fine-tunes intracellular calcium signals and eliminates myocardial damage through the IP3R/MCU pathways in cardiorenal syndrome type 3.Biochem Pharmacol. 2020 Apr;174:113832.;doi: 10.1016/j.bcp.2020.113832
- 132.Jin Liu,Yong-Ming Zhu,Yi Guo,Liang Lin,Zhan-Xiang Wang,Feng Gu,Xin-Yi Dong,Ming Zhou,Yi-Fan Wang,Hui-Ling Zhang.Inhibition of GSK3 β and RIP1K Attenuates Glial Scar Formation Induced by Ischemic Stroke via Reduction of Inflammatory Cytokine Production Front Pharmacol. 2020 Jun 12;11:812.;doi: 10.3389/fphar.2020.00812
- 133.Yin Zhang,Yong-Xin Huang,Dan-Lan Wang,Bing Yang,Hai-Yan Yan,Le-Hang Lin,Yun Li,Jie Chen,Li-Min Xie,Yong-Sheng Huang,Jian-You Liao,Kai-Shun Hu,Jie-Hua He,Phei Er Saw,Xiaoding Xu,Dong Yin.LncRNA DSCAM-AS1 interacts with YBX1 to promote cancer progression by forming a positive feedback loop that activates FOXA1 transcription network Theranostics. 2020 Aug 29;10(23):10823-10837.;doi: 10.7150/thno.47830
- 134.Yiyu Cai,Xinzhu Li,Zihao Pan,Yue Zhu,Ji Tuo,Qinghai Meng,Guoliang Dai,Guangming Yang,Yang Pan.Anthocyanin ameliorates hypoxia and ischemia induced inflammation and apoptosis by increasing autophagic flux in SH-SY5Y cells Eur J Pharmacol. 2020 Sep 15;883:173360.;doi: 10.1016/j.ejphar.2020.173360
- 135.Lingling Li,Dan Luo,Yating Liao,Kailan Peng,Yanhua Zeng.Mycoplasma genitalium Protein of Adhesion Induces Inflammatory Cytokines via Cyclophilin A-CD147 Activating the ERK-NF- κ B Pathway in Human Urothelial Cells Front Immunol. 2020 Sep 9;11:2052.;doi: 10.3389/fimmu.2020.02052
- 136.Yingqi Liu,Haiyan Wei,Jing Tang,Jiaming Yuan,Mingmin Wu,Chenjuan Yao,Kazuo Hosoi,Shali Yu,Xinyuan Zhao,Yu Han,Gang Chen.Dysfunction of pulmonary epithelial tight junction induced by silicon dioxide nanoparticles via the ROS/ERK pathway and protein degradation Chemosphere. 2020 Sep;255:126954.;doi: 10.1016/j.chemosphere.2020.126954
- 137.Yingqi Liu,Jing Tang,Jiaming Yuan,Chenjuan Yao,Kazuo Hosoi,Yu Han,Shali Yu,Haiyan Wei,Gang Chen.Arsenite-induced downregulation of occludin in mouse lungs and BEAS-2B cells via the ROS/ERK/ELK1/MLCK and ROS/p38 MAPK signaling pathways Toxicol Lett. 2020 Oct 10;332:146-154.;doi: 10.1016/j.toxlet.2020.07.010
- 138.Qing Yao,Yue Zhou,Yanhui Yang,Lianjun Cai,Lihui Xu,Xuebo Han,Yu Guo,P Andy Li.Activation of Sirtuin1 by lyceum barbarum polysaccharides in protection against diabetic cataract J Ethnopharmacol. 2020 Oct 28;261:113165.;doi: 10.1016/j.jep.2020.113165
- 139.Junqing Gan,Shumin Li,Yu Meng,Yuanyu Liao,Mingxia Jiang,Ling Qi,Yanjing Li,Yuxian Bai.The influence of photodynamic therapy on the Warburg effect in esophageal cancer cellsLASER MED SCI. 2020 Oct;35(8):1741-1750.;doi: 10.1007/s10103-020-02966-8
- 140.Wei Wang,Feiran Huang,Weifeng Jiang,Weiwei Wang,Jie Xiang.Brilliant blue G attenuates neuro-inflammation via regulating MAPKs and NF- κ B signaling pathways in lipopolysaccharide-induced BV2 microglia cellsExp Ther Med. 2020 Nov;20(5):116.;doi: 10.3892/etm.2020.9244
- 141.Biao Chen,Qing Chai,Shengmin Xu,Qian Li,Tao Wu,Shaopeng Chen,Lijun Wu.Silver nanoparticle-activated COX2/PGE2 axis involves alteration of lung cellular senescence in vitro and in vivo ECOTOX ENVIRON SAFE. 2020 Nov;204:111070.;doi: 10.1016/j.ecoenv.2020.111070
- 142.Guifang He,Mingfeng Li,Linna Fang,Lingfan Xu,Xiaojuan Huang,Lulu Zheng,Lijie Yang,Wenwu Luo,Yongping Cai,Wei Ma,Chaozhao Liang,Yu Yin.N-Myc induces the tumor progression of prostate cancer by regulating FSCN1Oncol Rep. 2020 Nov;44(5):2265-2274.;doi: 10.3892/or.2020.7775
- 143.Lanqing Sun,Sidi Yang,Qifeng Deng,Kedi Dong,Yuanyuan Li,Shuyan Wu,Rui Huang,Salmonella Effector SpvB Disrupts Intestinal Epithelial Barrier Integrity for Bacterial TranslocationFRONT CELL INFECT MI. 2020 Dec 17;10:606541.;doi: 10.3389/fcimb.2020.606541
- 144.Xiang-Xue Kong,Shuai Jiang,Tong Liu,Gao-Feng Liu,Mei Dong.Paclitaxel increases sensitivity of SKOV3 cells to hyperthermia by inhibiting heat shock protein 27 Biomed Pharmacother. 2020 Dec;132:110907.;doi: 10.1016/j.biopha.2020.110907
- 145.Jin Xu,Jiandong Gao,Di Huang,Pinglan Lin,Dongsheng Yao,Feng Yang,Yaheng Zhang,Xuejun Yang,Ming Wu,Chaoyang Ye.Emodin ameliorates tubulointerstitial fibrosis in obstructed kidneys by inhibiting EZH2 BIOCHEM BIOPH RES CO. 2021 Jan 1;534:279-285.;doi: 10.1016/j.bbrc.2020.11.094
- 146.Xiaolin Shi,Wei Zhang,Cheng Gu,Huangge Ren,Chen Wang,Narui Yin,Zhongmin Wang,Jiahua Yu,Fenju Liu,Haowen Zhang.NAD+ depletion radiosensitizes 2-DG-treated glioma cells by abolishing metabolic adaptation FREE RADICAL BIO MED. 2021 Jan;162:514-522.;doi: 10.1016/j.freeradbiomed.2020.11.007
- 147.Yan-Jie Wang,Bing Yang,Qiao Lai,Jun-Fang Shi,Jiang-Yun Peng,Yin Zhang,Kai-Shun Hu,Ya-Qing Li,Jing-Wen Peng,Zhi-Zhi Yang,Yao-Ting Li,Yue Pan,H Phillip Koefller,Jian-You Liao,Dong Yin.Reprogramming of m 6 A epitranscriptome is crucial for shaping of transcriptome and proteome in response to hypoxia RNA Biol. 2021 Jan;18(1):131-143.;doi: 10.1080/15476286.2020.1804697
- 148.Hui Yin,Yan-Ge Liu,Fei Li,Lun-Qing Wang,Jian-Hua Zha,Ying-Chen Xia,Ben-Tong Yu,Di-Hao Wen.Resibufogenin suppresses growth and metastasis through inducing caspase-1-dependent pyroptosis via ROS-mediated NF- κ B suppression in non-small cell lung cancer ANAT REC. 2021 Feb;304(2):302-312.;doi: 10.1002/ar.24415
- 149.Linqing Sun,Juntao Li,Wenying Yan,Zhendong Yao,Ruoqin Wang,Xiaojun Zhou,Hongya Wu,Guangbo Zhang,Tongguo Shi,Weichang Chen.H19 promotes aerobic glycolysis, proliferation, and immune escape of gastric cancer cells through the microRNA-519d-3p/lactate dehydrogenase A axis Cancer Sci. 2021 Mar 23.;doi: 10.1111/cas.14896
- 150.Shan Yu,Jingjie Zhai,Jing Yu,Qiwei Yang,Jinghui Yang.miR-98-5p protects against cerebral ischemia/reperfusion injury through anti-apoptosis and anti-oxidative stress in miceJ Biochem. 2021 Mar 5;169(2):195-206.;doi: 10.1093/jb/mvaa099
- 151.Baokun Sui,Dong Chen,Wei Liu,Bin Tian,Lei Lv,Jie Pei,Qiong Wu,Ming Zhou,Zhen F Fu,Yi Zhang,Ling Zhao.Comparison of lncRNA and mRNA expression in mouse brains infected by a wild-type and a lab-attenuated Rabies lyssavirus J Gen Virol. 2021 Mar;102(3);doi: 10.1099/jgv.0.001538
- 152.Yuhui Hao,Jiawei Huang,Yonghong Ran,Shuang Wang,Juan Li,Yazhen Zhao,Xinze Ran,Binghui Lu,Jing Liu,Rong Li.Ethylmalonic encephalopathy 1 initiates overactive autophagy in depleted uranium-induced cytotoxicity in the human embryonic kidney 293 cellsJ BIOCHEM MOL TOXIC. 2021 Mar;35(3):e22669.;doi: 10.1002/jbt.22669
- 153.Yue Zhou,Lian-Jun Cai,Li-Hui Xu,Yu Guo,Nuo Chen,Qing Yao.The role of sirt1 in the retinal ganglion cells cultured by high glucose Int Ophthalmol. 2021 Mar;41(3):845-852.;doi: 10.1007/s10792-020-01638-4
- 154.Shijiao Tian,Shichuan Tan,Wenming Jia,Juan Zhao,Xiulan Sun.Activation of Wnt/ β -catenin signaling restores insulin sensitivity in insulin resistant neurons through transcriptional regulation of IRS-1 J Neurochem. 2021 May;157(3):467-478.;doi: 10.1111/jnc.15277
- 155.Huimin Lu,Yanchao Ma,Mingyuan Wang,Jin Shen,Hongya Wu,Juntao Li,Nan Gao,Yanzheng Gu,Xueguang Zhang,Guangbo Zhang,Tongguo Shi,Weichang Chen.B7-H3 confers resistance to V γ 9V δ 2 T cell-mediated cytotoxicity in human colon cancer cells via the STAT3/ULBP2 axis CANCER IMMUNOL IMMUN. 2021 May;70(5):1213-1226.;doi: 10.1007/s00262-020-02771-w

Version 2021.09.01