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| Multiple Well Plates or Dishes | Single Well Only for Plates | | | | | |
|--------------------------------|-----------------------------|---------------------------------|-----------------------|------------------------|---------------------|-------------------------|
| | Diameter (Bottom, mm)* | Growth Area (cm ²)* | Average Cell Yield | Total Well Volume (ml) | Working Volume (ml) | Recommended Volume (ml) |
| 6 well | 34.8 | 9.5 | 9.5 × 10 ⁵ | 16.8 | 1.9-2.9 | 2 |
| 12 well | 22.1 | 3.8 | 3.8 × 10 ⁵ | 6.9 | 0.76-1.14 | 1 |
| 24 well | 15.6 | 1.9 | 1.9 × 10 ⁵ | 3.4 | 0.38-0.57 | 0.5 |
| 48 well | 11.0 | 0.95 | 9.5 × 10 ⁴ | 1.6 | 0.19-0.285 | 0.25 |
| 96 well | 6.4 | 0.32 | 3.2 × 10 ⁴ | 0.36 | 0.10-0.20 | 0.1 |
| 384 well | 2.7 | 0.056 | 5.6 × 10 ³ | 0.112 | 0.025-0.050 | 0.030 |
| 1536 well | 1.63 × 1.63** | 0.025 | 2.5 × 10 ³ | 0.0125 | 0.005-0.010 | 0.010 |
| 3.5 cm dish | 34 | 9 | 9.0 × 10 ⁵ | NA | 1.8-2.7 | 2 |
| 6 cm dish | 52 | 21 | 2.1 × 10 ⁶ | NA | 4.2-6.3 | 5 |
| 10 cm dish | 8.4 | 55 | 5.5 × 10 ⁶ | NA | 11-16.5 | 12 |
| 15cm dish | 14 | 152 | 1.5 × 10 ⁷ | NA | 30.4-45.6 | 35 |
| 24.5cm dish | 22.4 × 22.4** | 500 | 5.0 × 10 ⁷ | NA | 100-150 | 120 |

*Diameter and growth area may vary depending on the manufacturer, and the listed sizes are from Corning.

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- Chen HY, He LJ, Li SQ, Zhang YJ, Huang JH, Qin HX, Wang JL, Li QY, Yang DL. Derivate of Benzimidazole-Isoquinolinone Induces SKP2 Transcriptional Inhibition to Exert Anti-Tumor Activity in Glioblastoma Cells. *Molecules*. 24(15). pii: E2722.
- Xue S, Gong R, He F, Li Y, Wang Y, Tan T, Luo SZ. Low-complexity domain of U1-70K modulates phase separation and aggregation through distinctive basic-acidic motifs. *Sci Adv*. 5(11):eaax5349.
- Tan Y, Jin Y, Wu X, Ren Z. PSMD1 and PSMD2 regulate HepG2 cell proliferation and apoptosis via modulating cellular lipid droplet metabolism. *BMC Mol Biol*. 20(1):24.
- Guo Y, Nan X, Zhang X, Wang G, Ren Y, Wang Y, Fu Y, Hou J. Molecular characterization and functional analysis of Japanese flounder (*Paralichthys olivaceus*) thbs2 in response to lymphocystis disease virus. *FISH SHELLFISH IMMUN*. 93:183-190.
- Tan Y, Jin Y, Wang Q, Huang J, Wu X, Ren Z. Perilipin 5 Protects against Cellular Oxidative Stress by Enhancing Mitochondrial Function in HepG2 Cells. *Cells*. 8(10). pii: E1241.
- Xu Bai, Yu Wang, Bo Hu, Qi Cao, Maochen Xing, Shuliang Song, Aiguo Ji. Fucoidan Induces Apoptosis of HT-29 Cells via the Activation of DR4 and Mitochondrial Pathway. *Mar Drugs*. doi: 10.3390/md18040220.
- Zhen-Yu She, Kai-Wei Yu, Ning Zhong, Yu Xiao, Ya-Lan Wei, Yang Lin, Yue-Ling Li, Ming-Hui Lu. Kinesin-7 CENP-E regulates chromosome alignment and genome stability of spermatogenic cells. *RECENT PAT ANTI-CANC*. doi: 10.1038/s41420-020-0261-8.
- Miaomiao Chai, Ce Gu, Qihua Shen, Jiaying Liu, Yi Zhou, Ziyang Jin, Wanli Xiong, Yan Zhou, Wensong Tan. Hypoxia alleviates dexamethasone-induced inhibition of angiogenesis in cocultures of HUVECs and rBMSCs via HIF-1 α . *Stem Cell Res Ther*. doi: 10.1186/s13287-020-01853-x.
- J-J Chen, P Lei, M Zhou. hsa_circ_0121582 inhibits leukemia growth by dampening Wnt/ β -catenin signaling. *Clin Transl Oncol*. doi: 10.1007/s12094-020-02377-9.
- Wenjun Li, Qian Du, Xiaoli Li, Xiangru Zheng, Feng Lv, Xin Xi, Guili Huang, Jia Yang, Songqing Liu. Eriodictyol Inhibits Proliferation, Metastasis and Induces Apoptosis of Glioma Cells via PI3K/Akt/NF- κ B Signaling Pathway. *Front Pharmacol*. doi: 10.3389/fphar.2020.00114.
- Hao Niu, Yiwei Huang, Li Yan, Li Zhang, Mengnan Zhao, Tao Lu, Xiaodong Yang, Zhengcong Chen, Cheng Zhan, Yu Shi, Qun Wang. Knockdown of SMAD3 inhibits the growth and enhances the radiosensitivity of lung adenocarcinoma via p21 in vitro and in vivo. *Int J Biol Sci*. doi: 10.7150/ijbs.40173.
- Xiangru Zheng, Jinlai Wei, Wenjun Li, Xiaoli Li, Wuyi Wang, Jinbao Guo, Zhongxue Fu. PRDX2 removal inhibits the cell cycle and autophagy in colorectal cancer cells. *AGING-US*. doi: 10.18632/aging.103690.
- Xue Han, Jielong Guo, Yilin You, Jicheng Zhan, Weidong Huang. p-Coumaric acid prevents obesity via activating thermogenesis in brown adipose tissue mediated by mTORC1-RPS6. *FASEB J*. doi: 10.1096/fj.202000333R.
- Chuanhong Zhong, Bei Tao, Yitian Chen, Zhangchao Guo, Xiaobo Yang, Lilei Peng, Xiangguo Xia, Ligang Chen. B7-H3 Regulates Glioma Growth and Cell Invasion Through a JAK2/STAT3/Slug-Dependent Signaling Pathway. *ONCOTARGETS THER*. doi: 10.2147/OTT.S237841.
- Yang Liu, Yiling Ding, Yanting Nie, Mengyuan Yang. EMP1 Promotes the Proliferation and Invasion of Ovarian Cancer Cells Through Activating the MAPK Pathway. *ONCOTARGETS THER*. doi: 10.2147/OTT.S240028.
- Fanhua Kong, Heng Zou, Xi Liu, Jun He, Yanwen Zheng, Li Xiong, Xiongming Miao. miR-7112-3p targets PERK to regulate the endoplasmic reticulum stress pathway and apoptosis induced by photodynamic therapy in colorectal cancer CX-1 cells. *PHOTODIAGN PHOTODYN*. doi: 10.1016/j.pdpdt.2020.101663.
- Xinna Du, Xuan Shen, Li Dai, Fangjie Bi, Hu Zhang, Chunfeng Lu. PSMD12 promotes breast cancer growth via inhibiting the expression of pro-apoptotic genes. *BIOCHEM BIOPH RES CO*. doi: 10.1016/j.bbrc.2020.03.095.
- Ahui Song, Jingjing Wang, Yan Tong, Junyan Fang, Yi Zhang, Huiping Zhang, Hongqiang Ruan, Kai Wang, Yingli Liu. BKCa channels regulate the immunomodulatory properties of WJ-MSCs by affecting the exosome protein profiles during the inflammatory response. *Stem Cell Res Ther*. doi: 10.1186/s13287-020-01952-9.
- Xiaohan Zhang, Xian Lin, Chenghuang Qin, Kun Huang, Xiaomei Sun, Lianzhong Zhao, Meilin Jin. Avian Chaperonin Containing TCP1 Subunit 5 Supports Influenza A Virus Replication by Interacting With Viral Nucleoprotein, PB1, and PB2 Proteins. *Front Microbiol*. doi: 10.3389/fmicb.2020.538355.
- Tingting Cheng, Yun Feng, Xiaoyan Chen, Jian Zhou, Yuanlin Song. Lung-resident mesenchymal stem cells regulated the inflammatory responses in innate and adaptive immune cells through HVEM-BTLA pathway during ARDS. *Exp Cell Res*. doi: 10.1016/j.yexcr.2020.112155.
- Chengjie Xia, Wenyi Tang, Ping Geng, Haiyan Zhu, Wei Zhou, Hai Huang, Pei Zhou, Xunlong Shi. Baicalin down-regulating hepatitis B virus transcription depends on the liver-specific HNF4 α -HNF1 α axis. *TOXICOL APPL PHARM*. doi: 10.1016/j.taap.2020.115131.
- Qiang An, Ting Liu, Ming-Yang Wang, Yu-Jia Yang, Zhen-Dong Zhang, Zhen-Jiang Lin, Bing Yang. circKRT7-miR-29a-3p-COL1A1 Axis Promotes Ovarian Cancer Cell Progression. *ONCOTARGETS THER*. doi: 10.2147/OTT.S259033.
- Ya-Jing Liang, Jia-Hao Ni, Lan-Min Wang, Fu-Qing Tan, Wan-Xi Yang. KIF3A regulates the Wnt/ β -catenin pathway via transporting β -catenin during spermatogenesis in *Eriocheir sinensis*. *Cell Tissue Res*. doi: 10.1007/s00441-020-03220-x.
- Liuli Li, Hao Li, Jun Zhang, Xin Gao, Hao Jin, Renqi Liu, Zhen Zhang, Xuan Zhang, Xichang Wang, Peng Qu, Yuejiao Zhao, Xiaobo Lu. Bisphenol A at a human exposed level can promote epithelial-mesenchymal transition in papillary thyroid carcinoma harbouring BRAF V600E mutation. *J Cell Mol Med*. doi: 10.1111/jcmm.16279.
- Zhujun Yi, Yilin Wu, Wenfeng Zhang, Tao Wang, Jianping Gong, Yao Cheng, Chunmu Miao. Activator-Mediated Pyruvate Kinase M2 Activation Contributes to Endotoxin Tolerance by Promoting Mitochondrial Biogenesis. *Front Immunol*. doi: 10.3389/fimmu.2020.595316.
- Xiaoqun Lv, Jinguo Zhang, Jun Zhang, Wencai Guan, Weifang Ren, Yujuan Liu, Guoxiong Xu. A Negative Feedback Loop Between NAMPT and TGF- β Signaling Pathway in Colorectal Cancer Cells. *ONCOTARGETS THER*. doi: 10.2147/OTT.S282367.
- Xiuzhi Li, Lingling Sun, Guokai Yan, Xianghua Yan. PFKF facilitates ATG4B phosphorylation during amino acid deprivation-induced autophagy. *Cell Signal*.

- doi: 10.1016/j.cellsig.2021.109956.
28. Zhen-Yu She, Kai-Wei Yu, Ya-Lan Wei, Ning Zhong, Yang Lin. Kinesin-7 CENP-E regulates the formation and structural maintenance of the acrosome. *Cell Tissue Res.* doi: 10.1007/s00441-020-03341-3.
29. Hua Lin, Weifeng Zeng, Yuhang Lei, Desheng Chen, Zhen Nie. Tuftelin 1 (TUFT1) Promotes the Proliferation and Migration of Renal Cell Carcinoma via PI3K/AKT Signaling Pathway. *Pathol Oncol Res.* doi: 10.3389/pore.2021.640936.
30. Qiao-Qiao Han, Meng-Yan Deng, Hao Liu, Usman Ali, Xin-Yan Li, Yong-Xiang Wang. Cynandione A and PHA-543613 inhibit inflammation and stimulate macrophageal IL-10 expression following $\alpha 7$ nAChR activation. *Biochem Pharmacol.* doi: 10.1016/j.bcp.2021.114600.
31. Yuanlin Yu, Xiaopeng Chen, Weidong Zhang, Jun Liu. Abnormal Expression of Centromere Protein U Is Associated with Hepatocellular Cancer Progression. *Biomed Res Int.* doi: 10.1155/2021/4051192.
32. Menglin Zou, Jingfeng Zou, Xingxing Hu, Weishuai Zheng, Mingyang Zhang, Zhenshun Cheng. Latent Transforming Growth Factor- β Binding Protein-2 Regulates Lung Fibroblast-to-Myofibroblast Differentiation in Pulmonary Fibrosis via NF- κ B Signaling. *Front Pharmacol.* doi: 10.3389/fphar.2021.788714.
33. Mingfang Huang, Xuqiang Liao, Liang Li, Gao Li, Minbiao Chen. MiR-552-3p facilitated cell proliferation, migration and invasion by sponging Fibulin 5 in non-small cell lung cancer via activation of ERK/GSK3 β / β -catenin signaling pathway. *Tissue Cell.* doi: 10.1016/j.tice.2021.101672.
34. Ruifang Wang, Yinxing Zhu, Chenwei Ren, Shuaike Yang, Shan Tian, Huanchun Chen, Meilin Jin, Hongbo Zhou. Influenza A virus protein PB1-F2 impairs innate immunity by inducing mitophagy. *Autophagy.* doi: 10.1080/15548627.2020.1725375.
35. Li Jia, Jun Li, Pengfei Li, Didi Liu, Jing Li, Jiechen Shen, Bojing Zhu, Chen Ma, Ting Zhao, Rongxia Lan, Liuyi Dang, Wang Li, Shisheng Sun. Site-specific glycoproteomic analysis revealing increased core-fucosylation on FOLR1 enhances folate uptake capacity of HCC cells to promote EMT. *Theranostics.* doi: 10.7150/thno.56882.
36. Haorui Zhang, Bocen Chen, Zixiu Zou, Jian Feng, Yutao Li, Yi Wang, Xing He, Chang Xu, Haijian Wang, Shicheng Guo, Li Jin, Qiang Li, Jiucun Wang, Man Xiao, Feng Li, Junjie Wu. Associations Between CAMKK1 Polymorphism rs7214723 and the Prognosis of Patients With Lung Cancer. *Front Oncol.* doi: 10.3389/fonc.2021.757484.
37. Yang Zhou, Hui Yan, Qiang Zhou, Ruiling Feng, Penggao Wang, Fang Yang, Yaodong Zhang, Ziqiao Yuan, Bo Zhai. Beta-Lapachone Attenuates BMSC-Mediated Neuroblastoma Malignant Transformation by Inhibiting Gal-3/Gal-3BP/IL6 Axis. *Front Pharmacol.* doi: 10.3389/fphar.2021.766909.
38. Zhi-Chao Sun, Zeyin Jiang, Xiaowen Xu, Meifeng Li, Qing Zeng, Ying Zhu, Shanghong Wang, Yuanyuan Li, Xiao-Li Tian, Chengyu Hu. Fish Paralog Proteins RNASEK-a and -b Enhance Type I Interferon Secretion and Promote Apoptosis. *Front Immunol.* doi: 10.3389/fimmu.2021.762162.
39. Yan Zhang, Huiwen Ding, Xiaoxiao Wang, Xin Wang, Shengpeng Wan, Anchun Xu, Ruoyi Gan, Shou-Dong Ye. MK2 promotes Tfcp2l1 degradation via β -TrCP ubiquitin ligase to regulate mouse embryonic stem cell self-renewal. *Cell Rep.* doi: 10.1016/j.celrep.2021.109949.
40. Xiufang Ding, Siyu Gu, Song Xue, Shi-Zhong Luo. Disease-associated mutations affect TIA1 phase separation and aggregation in a proline-dependent manner. *Brain Res.* doi: 10.1016/j.brainres.2021.147589.
41. Yunyi Bian, Qihai Sui, Guoshu Bi, Yuansheng Zheng, Mengnan Zhao, Guangyu Yao, Liang Xue, Yi Zhang, Hong Fan. Identification and Validation of a Proliferation-Associated Score Model Predicting Survival in Lung Adenocarcinomas. *Dis Markers.* doi: 10.1155/2021/3219594.
42. Fanzhen Yu, Yingxin Wang, Yilong Teng, Shutong Yang, Yanming He, Zeng Zhang, Hongjie Yang, Chuan-Fan Ding, Ping Zhou. Interaction and Inhibition of a Ganoderma lucidum Proteoglycan on PTP1B Activity for Anti-diabetes. *ACS Omega.* doi: 10.1021/acsomega.1c04244.
43. Yadong Fu, Zhun Xiao, Xiaoting Tian, Wei Liu, Zhou Xu, Tao Yang, Yonghong Hu, Xiaoxi Zhou, Jing Fang, Siqi Gao, Dingqi Zhang, Yongping Mu, Hua Zhang, Yiyang Hu, Chenggang Huang, Jiamei Chen, Ping Liu. The Novel Chinese Medicine JY5 Formula Alleviates Hepatic Fibrosis by Inhibiting the Notch Signaling Pathway. *Front Pharmacol.* doi: 10.3389/fphar.2021.671152.
44. Dong-Lin Yang, Ya-Jun Zhang, Liu-Jun He, Chun-Sheng Hu, Li-Xia Gao, Jiu-Hong Huang, Yan Tang, Jie Luo, Dian-Yong Tang, Zhong-Zhu Chen. Demethylzeylasteral (T-96) initiates extrinsic apoptosis against prostate cancer cells by inducing ROS-mediated ER stress and suppressing autophagic flux. *Biol Res.* doi: 10.1186/s40659-021-00350-6.
45. Chengchen Wu, Junyi Zhao, Qiuping Wu, Qiulong Tan, Qiong Liu, Shifeng Xiao. Tau N-Terminal Inserts Regulate Tau Liquid-Liquid Phase Separation and Condensates Maturation in a Neuronal Cell Model. *Int J Mol Sci.* doi: 10.3390/ijms22189728.
46. Yaqi Bian, Yan Chen, Xiufen Wang, Guozhen Cui, Carolina Oi Lam Ung, Jia-Hong Lu, Weihong Cong, Benqin Tang, Simon Ming-Yuen Lee. Oxyphylla A ameliorates cognitive deficits and alleviates neuropathology via the Akt-GSK3 β and Nrf2-Keap1-HO-1 pathways in vitro and in vivo murine models of Alzheimer's disease. *J Adv Res.* doi: 10.1016/j.jare.2021.09.002.
47. Gangwei Jiang, Tian-Bing Ren, Elisa D'Este, Mengyi Xiong, Bin Xiong, Kai Johnsson, Xiao-Bing Zhang, Lu Wang, Lin Yuan. A synergistic strategy to develop photostable and bright dyes with long Stokes shift for nanoscopy. *Nat Commun.* doi: 10.1038/s41467-022-29547-3.
48. Meifeng Li, Jihuan Hu, Huiling Mao, Dongming Li, Zeyin Jiang, Zhichao Sun, Tingting Yu, Chengyu Hu, Xiaowen Xu. Grass Carp (Ctenopharyngodon idella) KAT8 Inhibits IFN 1 Response Through Acetylating IRF3/IRF7. *Front Immunol.* doi: 10.3389/fimmu.2021.808159.
49. Yansheng Zhai, Xiaoyan Huang, Keren Zhang, Yuchen Huang, Yanlong Jiang, Jingwei Cui, Zhe Zhang, Cookson K C Chiu, Weiye Zhong, Gang Li. Spatiotemporal-resolved protein networks profiling with photoactivation dependent proximity labeling. *Nat Commun.* doi: 10.1038/s41467-022-32689-z.
50. Zilong Zheng, Guokai Yan, Xiuzhi Li, Yuke Fei, Lingling Sun, Haonan Yu, Yaorong Niu, Weihua Gao, Qing Zhong, Xianghua Yan. Lysine crotonylation regulates leucine-deprivation-induced autophagy by a 14-3-3 ϵ -PPM1B axis. *Cell Rep.* doi: 10.1016/j.celrep.2022.111850.
51. Qinggong Wang, Qiuyuan Lu, Qiong Guo, Maikun Teng, Qingguo Gong, Xu Li, Yang Du, Zheng Liu, Yuyong Tao. Structural basis of the ligand binding and signaling mechanism of melatonin receptors. *Nat Commun.* doi: 10.1038/s41467-022-28111-3.
52. Chenyang Yi, Cong Cai, Ze Cheng, Yifan Zhao, Xu Yang, Yue Wu, Xiaoping Wang, Zehua Jin, Yaozu Xiang, Meilin Jin, Li Han, Anding Zhang. Genome-wide CRISPR-Cas9 screening identifies the CYTH2 host gene as a potential therapeutic target of influenza viral infection. *Cell Rep.* doi: 10.1016/j.celrep.2022.110559.
53. Xu Chen, Yihua Cai, Xiaoling Hu, Chuanlin Ding, Liqing He, Xiang Zhang, Fuxiang Chen, Jun Yan. Differential metabolic requirement governed by transcription factor c-Maf dictates innate γ DT17 effector functionality in mice and humans. *Sci Adv.* doi: 10.1126/sciadv.abm9120.
54. Ji Wang, Caiguangxi Deng, Ming Liu, Yihao Liu, Liubing Li, Zhangping Huang, Liru Shang, Juan Jiang, Yongyong Li, Ruohui Mo, Hui Zhang, Min Liu, Sui Peng, Haipeng Xiao. A fourth dose of the inactivated SARS-CoV-2 vaccine redistributes humoral immunity to the N-terminal domain. *Nat Commun.* doi: 10.1038/s41467-022-34633-7.
55. Zhuo-Hang Li, Xia Li, Fang-Fang Li, Zhi-Hong Wu, You-Long Xie, Shuang Zhang, Xue-Mei Chen, Ying-Xiong Wang, Yu-Bin Ding, Tai-Hang Liu. The roles of ADAMDEC1 in trophoblast differentiation during normal pregnancy and preeclampsia. *Mol Hum Reprod.* doi: 10.1093/molehr/gaac014.
56. Xiong Li, Jingwen Liu, Yan Zhou, Luman Wang, Yifan Wen, Kun Ding, Lu Zou, Xia Liu, Aiping Li, Yun Wang, Heling Fu, Min Huang, Guoxian Ding, Jianwei Zhou. Jwa participates the maintenance of intestinal epithelial homeostasis via ERK/FBXW7-mediated NOTCH1/PPAR γ /STAT5 axis and acts as a novel putative aging related gene. *Int J Biol Sci.* doi: 10.7150/ijbs.72751.

57. Xiangqi Hao, Hui Chen, Yanchao Li, Bo Chen, Weifeng Liang, Xiangyu Xiao, Pei Zhou, Shoujun Li. Molecular characterization and antiviral effects of canine interferon regulatory factor 1 (CaIRF1). *BMC Vet Res.* doi: 10.1186/s12917-022-03539-3.
58. Fanting Zeng, Yuhang Peng, Yuefeng Qin, Jianming Wang, Guoyun Jiang, Wenli Feng, Ying Yuan. Wee1 promotes cell proliferation and imatinib resistance in chronic myeloid leukemia via regulating DNA damage repair dependent on ATM-y H2AX-MDC1. *Cell Commun Signal.* doi: 10.1186/s12964-022-01021-z.
59. Adrian On-Wah Leung, Tsz-Ching Yiu, Lingxiao Liu, Hei-Yin Tam, Shen Gu, Jiajie Tu, Duanqing Pei, Hoi-Hung Cheung. Targeting G-quadruplex for rescuing impaired chondrogenesis in WRN-deficient stem cells. *Cell Biosci.* doi: 10.1186/s13578-022-00939-8.
60. Jun Wang, Beidi Wang, Biting Zhou, Jing Chen, Jia Qi, Le Shi, Shaojun Yu, Guofeng Chen, Muxing Kang, Xiaoli Jin, Lie Wang, Jinghong Xu, Linghua Zhu, Jian Chen. A novel immune-related lncRNA pair signature for prognostic prediction and immune response evaluation in gastric cancer: a bioinformatics and biological validation study. *Cancer Cell Int.* doi: 10.1186/s12935-022-02493-2.
61. Xiaohui Wu, Jianqiang Ye, Weixing Cai, Xi Yang, Qiuying Zou, Jingjing Lin, Hui Zheng, Chaoyun Wang, Liangwan Chen, Yumei Li. LDHA mediated degradation of extracellular matrix is a potential target for the treatment of aortic dissection. *Pharmacol Res.* doi: 10.1016/j.phrs.2021.106051.
62. Ning Jiao, Wan Sheng Ji, Biao Zhang, Yu Kui Shang, Yu Chen Zhang, Wei Qun Yu, Hai Long Jin, Chao Li, Cheng Ying Zhang, Cheng Yan, Wen Yue, Qing Zhang. Overexpression of Protein Phosphatase 2 Regulatory Subunit B^α Promotes Glycolysis by Regulating Hexokinase 1 in Hepatocellular Carcinoma. *Biomed Environ Sci.* doi: 10.3967/bes2022.082.
63. Haijun Cai, Hui Guo, Yixuan Deng, Jinhai Jiang, Siyuan Liu, Wenge He, Huangang Jian. RRM2 regulates osteogenesis of mouse embryo fibroblasts via the Wnt/ β -catenin signaling pathway. *Exp Ther Med.* doi: 10.3892/etm.2022.11542.
64. Yucheng Lu, Qingyun Liu, Zhiqi Huang, Xiuli Chen, Chunling Yang, Yueling Zhang, Yongzhen Zhao, Fan Wang. ATP synthase subunit e is a shrimp growth-associated breeding marker. *Genomics.* doi: 10.1016/j.ygeno.2022.110410.
65. Xiangyou Yu, Ning Jiang, Jing Li, Xiaofeng Li, Shenglin He. Upregulation of BRD7 protects podocytes against high glucose-induced apoptosis by enhancing Nrf2 in a GSK-3 β -dependent manner. *Tissue Cell.* doi: 10.1016/j.tice.2022.101813.
66. Hongchao Zhou, Yuting Zhang, Jingjing Wang, Yuchao Yan, Yi Liu, Xiaojie Shi, Qi Zhang, Xingang Xu. The CREB and AP-1-Dependent Cell Communication Network Factor 1 Regulates Porcine Epidemic Diarrhea Virus-Induced Cell Apoptosis Inhibiting Virus Replication Through the p53 Pathway. *Front Microbiol.* doi: 10.3389/fmicb.2022.831852.
67. Mengnan Zhao, Ming Li, Yuansheng Zheng, Zhengyang Hu, Jiaqi Liang, Guoshu Bi, Yunyi Bian, Qihai Sui, Cheng Zhan, Miao Lin, Qun Wang. Identification and analysis of a prognostic ferroptosis and iron-metabolism signature for esophageal squamous cell carcinoma. *J Cancer.* doi: 10.7150/jca.68568.
68. Jing Yang, Sheng Lin, Honglu Sun, Zimin Chen, Fanli Yang, Xi Lin, Liyan Guo, Lingling Wang, Ao Wen, Xindan Zhang, Yushan Dai, Bin He, Yu Cao, Haohao Dong, Xianbo Liu, Bo Chen, Jian Li, Qi Zhao, Guangwen Lu. A Potent Neutralizing Nanobody Targeting the Spike Receptor-Binding Domain of SARS-CoV-2 and the Structural Basis of Its Intimate Binding. *Front Immunol.* doi: 10.3389/fimmu.2022.820336.
69. Liying Dong, Meng Wang, Xiaolei Gao, Xuan Zheng, Yixin Zhang, Liangjie Sun, Na Zhao, Chong Ding, Zeyun Ma, Yixiang Wang. miR-9-5p promotes myogenic differentiation via the Dlx3/Myf5 axis. *PeerJ.* doi: 10.7717/peerj.13360.
70. Mengchan Zhu, Ling Ye, Guiping Zhu, Yingying Zeng, Chengyu Yang, Hui Cai, Yuqing Mo, Xixi Song, Xin Gao, Wenjun Peng, Jian Wang, Meiling Jin. ROS-Responsive miR-150-5p Downregulation Contributes to Cigarette Smoke-Induced COPD via Targeting IRE1 α . *Oxid Med Cell Longev.* doi: 10.1155/2022/5695005.
71. Xingang Xu, Lixiang Wang, Yi Liu, Xiaojie Shi, Yuchao Yan, Shuxia Zhang, Qi Zhang. TRIM56 overexpression restricts porcine epidemic diarrhoea virus replication in Marc-145 cells by enhancing TLR3-TRAF3-mediated IFN- β antiviral response. *J Gen Virol.* doi: 10.1099/jgv.0.001748.
72. Wang Bingxiang, Zhao Ningning, Wang Yakun, Sun Lingyu, Li Hongmei, Xiao Zhang, Xiaomin Zhao. Characterization of the Eimeria tenella rhopty protein with a nuclear localization sequence (EtROP30). *Parasitol Res.* doi: 10.1007/s00436-022-07499-3.
73. Zhuo Zhang, Xu Zhang, Xin Wu, Yan Zhang, Jie Lu, Dan Li. Sirt1 attenuates astrocyte activation via modulating Dnajb1 and chaperone-mediated autophagy after closed head injury. *Cereb Cortex.* doi: 10.1093/cercor/bhac007.
74. Xiaopeng Hu, Jiahua Cui, Jun Chen, Shujuan Du, Xinyu Wang, Yabin Zhang, Jiajun Qian, Haifeng Chen, Fang Wei, Qiliang Cai, Jinping Jia, Ji Wu. Identification of hACE2-interacting sites in SARS-CoV-2 spike receptor binding domain for antiviral drugs screening. *Virus Res.* doi: 10.1016/j.virusres.2022.198915.
75. Hui Jiang, Feng Chen, DianZe Song, Xiaoqin Zhou, Long Ren, Mei Zeng. Dynamin-Related Protein 1 Is Involved in Mitochondrial Damage, Defective Mitophagy, and NLRP3 Inflammasome Activation Induced by MSU Crystals. *Oxid Med Cell Longev.* doi: 10.1155/2022/5064494.
76. Cui Shengyu, Luo Yinhua, Li Yuanhong, Zhao Jinbo, Fang Can, Xia Hao, Zhang Changjiang. Selenium alleviates heart remodeling through Sirt1/AKT/GSK-3 β pathway. *Int Immunopharmacol.* doi: 10.1016/j.intimp.2022.109158.
77. Xiaoyan Ren, Jiangang Yu, Lili Guo, Zaili Zhang. CDGSH iron sulfur domain 2 mitigates apoptosis, oxidative stress and inflammatory response caused by oxygen-glucose deprivation/reoxygenation in HT22 hippocampal neurons by Akt-Nrf2-activated pathway. *Metab Brain Dis.* doi: 10.1007/s11011-022-01043-z.
78. Chunming Gu, Wenhao Zhang, Enze Yang, Congyou Gu, Zhaoyang Zhang, Jing Ke, Xiong Wang, Shengying Wu, Shan Li, Fuyun Wu. Blockage of Orai1-Nucleolin interaction mediated calcium influx attenuates breast cancer cells growth. *Oncogenesis.* doi: 10.1038/s41389-022-00429-z.
79. Qiyang Yao, Xin Zou, Shihe Liu, Haowen Wu, Qiyang Shen, Jihong Kang. Oxidative Stress as a Contributor to Insulin Resistance in the Skeletal Muscles of Mice with Polycystic Ovary Syndrome. *Int J Mol Sci.* doi: 10.3390/ijms231911384.
80. Yang Liang, Yuefeng Qin, Guoyun Jiang, Wenli Feng, Ying Yuan. Targeting MDC1 promotes apoptosis and sensitizes Imatinib resistance in CML cells by mainly disrupting non-homologous end-joining repair. *Med Oncol.* doi: 10.1007/s12032-022-01821-w.
81. Wenpeng Zhao, Maolin Xu, Herman W Barkema, Xiaochen Xie, Yushan Lin, Sohrab Khan, John P Kastelic, Dong Wang, Zhaoju Deng, Bo Han. Prototheca bovis induces autophagy in bovine mammary epithelial cells via the HIF-1 α and AMPK α /ULK1 pathway. *Front Immunol.* doi: 10.3389/fimmu.2022.934819.
82. Hao-Zheng Li, Ya-Fen Wang, Xiao-Xia Shao, Ya-Li Liu, Zeng-Guang Xu, Shi-Long Wang, Zhan-Yun Guo. FAM237A, rather than peptide PEN and proCCK56-63, binds to and activates the orphan receptor GPR83. *FEBS J.* doi: 10.1111/febs.16765.
83. Yihui Xu, Bao Ren, Min Wang. HIF-1 α contributes to metastasis in choriocarcinoma by regulating DECI1 expression. *Clin Transl Oncol.* doi: 10.1007/s12094-022-03055-8.

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