

细胞周期与细胞凋亡检测试剂盒

产品编号	产品名称	包装
C1052	细胞周期与细胞凋亡检测试剂盒	50次

产品简介:

- 碧云天生产的细胞周期与细胞凋亡检测试剂盒(Cell Cycle and Apoptosis Analysis Kit)是一种采用经典的碘化丙啶染色(Propidium staining, 即PI staining)方法进行细胞周期与细胞凋亡分析的检测试剂盒。
- 碘化丙啶(Propidium, 简称PI)是一种双链DNA的荧光染料。碘化丙啶和双链DNA结合后可以产生荧光, 并且荧光强度和双链DNA的含量成正比。细胞内的DNA被碘化丙啶染色后, 可以用流式细胞仪对细胞进行DNA含量测定, 然后根据DNA含量的分布情况, 可以进行细胞周期和细胞凋亡分析。
- 碘化丙啶染色后, 假设G₀/G₁期细胞的荧光强度为1, 那么含有双份基因组DNA的G₂/M期细胞的荧光强度的理论值为2, 正在进行DNA复制的S期细胞的荧光强度为1-2之间。凋亡细胞由于细胞核发生浓缩以及发生DNA片段化(DNA fragmentation)导致部分基因组DNA片段在染色过程中丢失, 因此凋亡细胞碘化丙啶染色后呈现明显的弱染, 即荧光强度小于1, 在流式检测的荧光图上出现所谓的sub-G₁峰, 即凋亡细胞峰。
- 细胞发生凋亡时, 由于胞浆和染色质浓缩、核碎裂, 产生凋亡小体, 使细胞的光散射性质发生变化。在细胞凋亡的早期, 细胞对前向角光散射的能力显著降低, 对侧向光散射的能力增加或没有变化。在细胞凋亡的晚期, 前向和侧向光散射的信号均降低。因此可通过流式细胞仪测定细胞光散射的变化观察细胞凋亡情况。使用本试剂盒检测经凋亡诱导的HeLa细胞的效果请参考图1。




图1. 碧云天细胞周期与细胞凋亡检测试剂盒(C1052)检测经凋亡诱导的HeLa细胞的效果图。HeLa细胞经适当凋亡诱导后, 用70%乙醇固定, 并经本试剂盒染色后进行流式检测以分析细胞周期和细胞凋亡。图中可见凋亡细胞的Sub-G₁峰。实际检测效果会因实验条件、检测仪器等的不同而存在差异, 图中效果仅供参考。

- 本试剂盒通常应用于培养的贴壁或悬浮细胞的细胞周期与细胞凋亡检测。如果用于组织的细胞周期与细胞凋亡检测, 则必须把组织消化成单细胞状态, 才可以进行检测。
- 本试剂盒足够检测50个样品, 每个样品的细胞数量可以为10-100万。

包装清单:

产品编号	产品名称	包装
C1052-1	染色缓冲液	25ml
C1052-2	碘化丙啶染色液(20X)	1.25ml
C1052-3	RNase A (50X)	0.5ml
—	说明书	1份

保存条件:

-20°C保存, 一年有效。C1052-2碘化丙啶染色液(20X)需避光保存。本试剂盒可4°C保存, 一个月内有效。

注意事项:

- 本试剂盒需使用流式细胞仪进行检测。

- 需自备PBS和70%乙醇，PBS (C0221A)可以向碧云天订购。
- 荧光染料均存在淬灭问题，保存和使用过程中请尽量注意避光，以减缓荧光淬灭。
- 碘化丙啶对人体有刺激性，操作时请小心，并注意适当防护以避免直接接触人体或吸入体内。
- 本产品仅限于专业人员的科学研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用说明：

1. 细胞样品的准备：

- 对于贴壁细胞：小心收集细胞培养液到一离心管内备用。用胰酶消化细胞，至细胞可以被轻轻用移液管或枪头吹打下来时，加入前面收集的细胞培养液，吹打下所有的贴壁细胞，并轻轻吹散细胞。再次收集到离心管内。1000g左右离心3-5分钟，沉淀细胞。对于特定的细胞，如果细胞沉淀不充分，可以适当延长离心时间或稍稍加大离心力。小心吸除上清，可以残留约50微升左右的培养液，以避免吸走细胞。加入约1毫升冰浴预冷的PBS，重悬细胞，并转移到1.5毫升离心管内。再次离心沉淀细胞，小心吸除上清，可以残留约50微升左右的PBS，以避免吸走细胞。轻轻弹击离心管底以适当分散细胞，避免细胞成团。
 - 对于悬浮细胞：1000g左右离心3-5分钟，沉淀细胞。对于特定的细胞，如果细胞沉淀不充分，可以适当延长离心时间或稍稍加大离心力。小心吸除上清，可以残留约50微升左右的培养液，以避免吸走细胞。加入约1毫升冰浴预冷的PBS，重悬细胞，并转移到1.5毫升离心管内。再次离心沉淀细胞，小心吸除上清，可以残留约50微升左右的PBS，以避免吸走细胞。轻轻弹击离心管底以适当分散细胞，避免细胞成团。
- 细胞固定：加入1毫升冰浴预冷70%乙醇中，轻轻吹打混匀，4°C固定30分钟或更长时间。通常固定2小时或以上更能保证染色效果，固定12-24小时可能效果更佳。1000g左右离心3-5分钟，沉淀细胞。对于特定的细胞，如果细胞沉淀不充分，可以适当延长离心时间或稍稍加大离心力。小心吸除上清，可以残留约50微升左右的70%乙醇，以避免吸走细胞。加入约1毫升冰浴预冷的PBS，重悬细胞。再次离心沉淀细胞，小心吸除上清，可以残留约50微升左右的PBS，以避免吸走细胞。轻轻弹击离心管底以适当分散细胞，避免细胞成团。

3. 碘化丙啶染色液的配制：参考下表，根据待检测样品的数量配制适量的碘化丙啶染色液：

	1个样品	6个样品	12个样品
染色缓冲液	0.5ml	3ml	6ml
碘化丙啶染色液(20X)	25μl	150μl	300μl
RNase A (50X)	10μl	60μl	120μl
Final volume	0.535ml	3.21ml	6.42ml

注：配制好的碘化丙啶染色液短时间内可以4°C保存，宜当日使用。

- 染色：每管细胞样品中加入0.5毫升碘化丙啶染色液，缓慢并充分重悬细胞沉淀，37°C避光温浴30分钟。随后可以4°C或冰浴避光存放。染色完成后宜在24小时内完成流式检测，最好能在当日完成流式检测。
- 流式检测和分析：用流式细胞仪在激发波长488nm波长处检测红色荧光，同时检测光散射情况。采用适当分析软件进行细胞DNA含量分析和光散射分析。

使用本产品的文献：

- Jiang X, Sha X, Xin H, Chen L, Gao X, Wang X, Law K, Gu J, Chen Y, Jiang Y, Ren X, Ren Q, Fang X. Self-aggregated pegylated poly (trimethylene carbonate) nanoparticles decorated with c(RGDyK)peptide for targeted paclitaxel delivery to integrin-rich tumors. *Biomaterials*. 2011 Dec;32(35):9457-69.
- He Y, Huang C, Sun X, Long XR, Lv XW, Li J. MicroRNA-146a modulates TGF-beta1-induced hepatic stellate cell proliferation by targeting SMAD4. *Cell Signal*. 2012 Oct;24(10):1923-30.
- Lin LY, Bao YL, Chen Y, Sun LG, Yang XG, Liu B, Lin ZX, Zhang YW, Yu CL, Wu Y, Li YX. N-Benzoyl-12-nitrodehydroabietylamine-7-one, a novel dehydroabietylamine derivative, induces apoptosis and inhibits proliferation in HepG2 cells. *CHEM-BIOL INTERACT*. 2012 Aug 30;199(2):63-73.
- Yao Y, Zhang YW, Sun LG, Liu B, Bao YL, Lin H, Zhang Y, Zheng LH, Sun Y, Yu CL, Wu Y, Wang GN, Li YX. Juglanthraquinone C, a novel natural compound derived from Juglans mandshurica Maxim, induces S phase arrest and apoptosis in HepG2 cells. *Apoptosis*. 2012 Aug;17(8):832-41.
- Hu L, Cao D, Li Y, He Y, Guo K. Resveratrol sensitized leukemia stem cell-like KG-1a cells to cytokine-induced killer cells-mediated cytotoxicity through NKG2D ligands and TRAIL receptors. *Cancer Biol Ther*. 2012 May;13(7):516-26.
- Zhao J, Qiu H, Chen DL, Zhang WX, Zhang DC, Li M. Development of nanofibrous scaffolds for vascular tissue engineering. *Int J Biol Macromol*. 2013 May;56:106-13.
- Zhu K, Guo X, Guo X, Peng W, Zhou H. Protective effects of wheat germ protein isolate hydrolysates (WGPIH) against hydrogen peroxide-induced oxidative stress in PC12 cells. *Food Research International*. 2013 Aug;53(1):297-303.
- Shen XH, Xu SJ, Jin CY, Ding F, Zhou YC, Fu GS. Interleukin-8 prevents oxidative stress-induced human endothelial cell senescence via telomerase activation. *Int Immunopharmacol*. 2013 Jun;16(2):261-7.
- Jiang Z, Wu M, Miao J, Duan H, Zhang S, Chen M, Sun L, Wang Y, Zhang X, Zhu X, Zhang L. Deoxydophyllotoxin exerts both anti-angiogenic and vascular disrupting effects. *INT J BIOCHEM CELL B*. 2013 May 20;45(8):1710-1719.
- Zhang H, Yao D, Yao M, Huang H, Wang L, Yan M, Yan X, Gu X, Wu W, Lu S. Annexin A2 silencing inhibits invasion, migration, and tumorigenic potential of hepatoma cells. *WORLD J GASTROENTERO*. 2013 June 28; 19(24): 3792-801.
- Yu HY, Zhang XQ, Li X, Zeng FB, Ruan HL. 2-Methoxyjuglone Induces Apoptosis in HepG2 Human Hepatocellular Carcinoma Cells and Exhibits in Vivo Antitumor Activity in a H22 Mouse Hepatocellular Carcinoma Model. *J Nat Prod*. 2013 May 24;76(5):889-95.
- Wang Y, Liu P, Qiu L, Sun Y, Zhu M, Gu L, Di W, Duan Y. Toxicity and therapy of cisplatin-loaded EGF modified mPEG-PLGA-PLL nanoparticles for SKOV3 cancer in mice. *Biomaterials*. 2013 May;34(16):4068-77.
- Wei Y, Xu Y, Han X, Qi Y, Xu L, Xu Y, Yin L, Sun H, Liu K, Peng J. Anti-cancer effects of dioscin on three kinds of human lung cancer cell lines through inducing DNA damage and activating mitochondrial signal pathway. *Food Chem Toxicol*. 2013 Jun 11;59C:118-128.
- Lei Y, Huang T, Su M, Luo J, Korteweg C, Li J, Chen Z, Qiu Y, Liu X, Yan M, Wang Y, Gu J. Expression and distribution of immunoglobulin G in the normal liver,

- hepatocarcinoma and postpartial hepatectomy liver. *Lab Invest.* 2014 Nov;94(11):1283-95.
15. Yuan W, Kuai R, Ran R, Fu L, Yang Y, Qin Y, Liu Y, Tang J, Fu H, Zhang Q, Yuan M, Zhang Z, Gao F, He Q. Increased delivery of doxorubicin into tumor cells using extracellularly activated TAT functionalized liposomes: in vitro and in vivo study. *J BIOMED NANOTECHNOL.* 2014 Aug;10(8):1563-73.
 16. Yao K, Xing H, Yang W, Liao A, Wu B, Li Y, Zhang R, Liu Z. Knockdown of RLIP76 expression by RNA interference inhibits proliferation, enhances apoptosis, and increases chemosensitivity to daunorubicin in U937 leukemia cells. *TUMOR BIOL.* 2014 Aug;35(8):8023-31.
 17. Wang L, Li MD, Cao PP, Zhang CF, Huang F, Xu XH, Liu BL, Zhang M, Astin B, a cyclic pentapeptide from *Aster tataricus*, induces apoptosis and autophagy in human hepatic L-02 cells. *CHEM-BIOL INTERACT.* 2014 Sep 16;223C:1-9.
 18. He Y, Nie Y, Xie L, Song H, Gu Z. p53 mediated apoptosis by reduction sensitive shielding ternary complexes based on disulfide linked PEI ternary complexes. *Biomaterials.* 2014 Feb;35(5):1657-66.
 19. Long XR, He Y, Huang C, Li J. MicroRNA-148a is silenced by hypermethylation and interacts with DNA methyltransferase 1 in hepatocellular carcinogenesis. *Int J Oncol.* 2014 Jun;44(6):1915-22.
 20. He Y, Wu YT, Huang C, Meng XM, Ma TT, Wu BM, Xu FY, Zhang L, Lv XW, Li J. Inhibitory effects of long noncoding RNA MEG3 on hepatic stellate cells activation and liver fibrogenesis. *BBA-BIOMEMBRANES.* 2014 Nov;1842(11):2204-15.
 21. Huang J, Zhang H, Yu Y, Chen Y, Wang D, Zhang G, Zhou G, Liu J, Sun Z, Sun D, Lu Y, Zhong Y. Biodegradable self-assembled nanoparticles of poly (D,L-lactide-co-glycolide)/hyaluronic acid block copolymers for target delivery of docetaxel to breast cancer. *Biomaterials.* 2014 Jan;35(1):550-66.
 22. Wu Y, Liu X, Zhou Q, Huang C, Meng X, Xu F, Li J. Silent information regulator 1 (SIRT1) ameliorates liver fibrosis via promoting activated stellate cell apoptosis and reversion. *TOXICOL APPL PHARM.* 2015 Dec 1;289(2):163-76.
 23. Liang Z, Li S, Xu X, Xu X, Wang X, Wu J, Zhu Y, Hu Z, Lin Y, Mao Y, Chen H, Luo J, Liu B, Zheng X, Xie L. MicroRNA-576-3p inhibits proliferation in bladder cancer cells by targeting cyclin D1. *Mol Cells.* 2015;38(2):130-7.
 24. Ma S, Yang LL, Niu T, Cheng C, Zhong L, Zheng MW, Xiong Y, Li LL, Xiang R, Chen LJ, Zhou Q, Wei YQ, Yang SY. SKLB-67, an FLT3 and Wnt/ β -catenin signaling inhibitor, displays potent activity in models of FLT3-driven AML. *SCI REP-UK.* 2015 Oct 26;5:15646.
 25. Zhou B, Shu B, Xi T, Su N, Liu J. Dub3 expression correlates with tumor progression and poor prognosis in human epithelial ovarian cancer. *Biomed Pharmacother.* 2015 Mar;70:84-9.
 26. Gao C, Yu Z, Liu S, Xin H, Li X. Overexpression of CUGBP1 is associated with the progression of non-small cell lung cancer. *TUMOR BIOL.* 2015 Jun;36(6):4583-9.
 27. Yu F, Guo Y, Chen B, Dong P, Zheng J. MicroRNA-17-5p activates hepatic stellate cells through targeting of Smad7. *Lab Invest.* 2015 Jul;95(7):781-9.
 28. Meng H, Wang K, Chen X, Guan X, Hu L, Xiong G, Li J, Bai Y. MicroRNA-330-3p functions as an oncogene in human esophageal cancer by targeting programmed cell death 4. *Am J Cancer Res.* 2015 Feb 15;5(3):1062-75.
 29. Qian W, Kang L, Zhang T, Meng M, Wang Y, Li Z, Xia Q, Cheng D. Ecdysone receptor (EcR) is involved in the transcription of cell cycle genes in the silkworm. *Int J Mol Sci.* 2015 Feb 3;16(2):3335-49.
 30. Wang AT, Liang DS, Liu YJ, Qi XR. Roles of ligand and TPGS of micelles in regulating internalization, penetration and accumulation against sensitive or resistant tumor and therapy for multidrug resistant tumors. *Biomaterials.* 2015 Jun;53:160-72.
 31. Gao L, Wang Y, Xu Z, Li X, Wu J, Liu S, Chu P, Sun Z, Sun B, Lin Y, Peng J, Han G, Wang S, Tang Z. SZC017, a novel oleanolic acid derivative, induces apoptosis and autophagy in human breast cancer cells. *Apoptosis.* 2015 Dec;20(12):1636-50.
 32. Xu T, Ni MM, Xing-Li, Li XF, Meng XM, Huang C, Li J. NLRCS regulates TGF- β 1-induced proliferation and activation of hepatic stellate cells during hepatic fibrosis. *INT J BIOCHEM CELL B.* 2016 Jan;70:92-104.
 33. Liu J, Wang H, Wang B, Chen T, Wang X, Huang P, Xu L, Guo Z. Microcystin-LR promotes proliferation by activating Akt/S6K1 pathway and disordering apoptosis and cell cycle associated proteins phosphorylation in HL7702 cells. *Toxicol Lett.* 2016 Jan 5;240(1):214-25.
 34. Li X, Li T, Chen D, Zhang P, Song Y, Zhu H, Xiao Y, Xing Y. Overexpression of lysine-specific demethylase 1 promotes androgen-independent transition of human prostate cancer LNCaP cells through activation of the AR signaling pathway and suppression of the p53 signaling pathway. *Oncol Rep.* 2016 Jan;35(1):584-92.
 35. Yu LT, Yang MQ, Liu JL, Alfred MO, Li X, Zhang XQ, Zhang J, Wu MY, Wang M, Luo C. Recombinant Reg3 α protein protects against experimental acute pancreatitis in mice. *Mol Cell Endocrinol.* 2016 Feb;422(15):150-9.
 36. Li X, Li T, Chen D, Zhang P, Song Y, Zhu H, Xiao Y, Xing Y. Overexpression of lysine-specific demethylase 1 promotes androgen-independent transition of human prostate cancer LNCaP cells through activation of the AR signaling pathway and suppression of the p53 signaling pathway. *Oncol Rep.* 2016 Jan;35(1):584-92.
 37. Wang D, Zhao Y, Wang Y, Rong Y, Qin H, Bao Y, Song Z, Yu C, Sun L, Li Y. 25-methoxyl-dammarane-3 β , 12 β , 20-triol and artemisinin synergistically inhibit MDA-MB-231 cell proliferation through downregulation of testes-specific protease 50 (TSP50) expression. *TUMOR BIOL.* 2016 Sep;37(9):11805-11813.
 38. Liu J, Wang H, Wang B, Chen T, Wang X, Huang P, Xu L, Guo Z. Microcystin-LR promotes proliferation by activating Akt/S6K1 pathway and disordering apoptosis and cell cycle associated proteins phosphorylation in HL7702 cells. *Toxicol Lett.* 2016 Jan 5;240(1):214-25.
 39. Mi XG, Song ZB, Sun LG, Bao YL, Yu CL, Wu Y, Li YX. Cardamonin inhibited cell viability and tumorigenesis partially through blockade of testes-specific protease 50-mediated nuclear factor-kappaB signaling pathway activation. *INT J BIOCHEM CELL B.* 2016 Apr;73:63-71.
 40. Jiao Y, Ma S, Wang Y, Li J, Shan L, Sun J, Chen J. Methacryloxyethyl Cetyl Ammonium Chloride Induces DNA Damage and Apoptosis in Human Dental Pulp Cells via Generation of Oxidative Stress. *Int J Biol Sci.* 2016 Mar 21;12(5):580-93.
 41. Wang Q, Wang Y, Xing Y, Yan Y, Guo P, Zhuang J, Qin F, Zhang J. Physcion 8-O- β -glucopyranoside induces apoptosis, suppresses invasion and inhibits epithelial to mesenchymal transition of hepatocellular carcinoma HepG2 cells. *Biomed Pharmacother.* 2016 Oct;83:372-380.
 42. Pang MJ, Yang Z, Zhang XL, Liu ZF, Fan J, Zhang HY. Physcion, a naturally occurring anthraquinone derivative, induces apoptosis and autophagy in human nasopharyngeal carcinoma. *Acta Pharmacol Sin.* 2016 Dec;37(12):1623-1640.
 43. Yang JB, Khan M, He YY, Yao M, Li YM, Gao HW, Ma TH. Tubeimoside-1 induces oxidative stress-mediated apoptosis and G0/G1 phase arrest in human prostate carcinoma cells in vitro. *Acta Pharmacol Sin.* 2016 Jul;37(7):950-62.
 44. Wang L, Wu X, Huang P, Lv Z, Qi Y, Wei X, Yang P, Zhang F, JQ. A small molecule inhibitor of BRD, suppresses cell growth and invasion in oral squamous cell carcinoma. *Oncol Rep.* 2016 Oct;36(4):1989-96.
 45. Liang YX, Lu JM, Mo RJ, He HC, Xie J, Jiang FN, Lin ZY, Chen YR, Wu YD, Luo HW, Luo Z, Zhong WD. E2F1 promotes tumor cell invasion and migration through regulating CD147 in prostate cancer. *Int J Oncol.* 2016 Apr;48(4):1650-8.
 46. Chen PJ, Cai SP, Yang Y, Li WX, Huang C, Meng XM, Li J. PTP1B confers liver fibrosis by regulating the activation of hepatic stellate cells. *TOXICOL APPL PHARM.* 2016 Feb 1;292:8-18.
 47. Zhou DD, Wang X, Wang Y, Xiang XJ, Liang ZC, Zhou Y, Xu A, Bi CH, Zhang L. MicroRNA-145 inhibits hepatic stellate cell activation and proliferation by targeting ZEB2 through Wnt/ β -catenin pathway. *Mol Immunol.* 2016 Jul;75:151-60.
 48. Zhang X, Xu J, Jiang T, Liu G, Wang D, Lu Y. MicroRNA-195 suppresses colorectal cancer cells proliferation via targeting FGF2 and regulating Wnt/ β -catenin pathway. *Am J Cancer Res.* 2016 Nov 1;6(11):2631-2640.
 49. Xu T, Ni MM, Xing-Li, Li XF, Meng XM, Huang C, Li J. NLRCS regulates TGF- β 1-induced proliferation and activation of hepatic stellate cells during

- hepaticfibrosis. *INT J BIOCHEM CELL B*. 2016 Jan;70:92-104.
50. He YH, Li MF, Zhang XY, Meng XM, Huang C, Li J. NLR5 promotes cell proliferation via regulating the AKT/VEGF-A signaling pathway in hepatocellularcarcinoma. *Toxicology*. 2016 Jun 1;359-360:47-57.
51. Zhang D, Li Y, Wang R, Li Y, Shi P, Kan Z, Pang X. Inhibition of REST Suppresses Proliferation and Migration in Glioblastoma Cells. *Int J Mol Sci*. 2016 May 3;17(5).
52. Li X, Wu XQ, Xu T, Li XF, Yang Y, Li WX, Huang C, Meng XM, Li J. Role of histone deacetylases(HDACs) in progression and reversal of liver fibrosis. *TOXICOL APPL PHARM*. 2016 Sep 1;306:58-68.
53. Chang W, Liu M, Xu J, Fu H, Zhou B, Yuan T, Chen P. MiR-377 inhibits the proliferation of pancreatic cancer by targeting Pim-3. *TUMOR BIOL*. 2016 Nov;37(11):14813-14824.
54. Ren K, Li Y, Lu H, Li Z, Li Z, Wu K, Li Z, Han X. Long Noncoding RNA HOTAIR Controls Cell Cycle by Functioning as a Competing Endogenous RNA inEsophageal Squamous Cell Carcinoma. *Transl Oncol*. 2016 Dec;9(6):489-497.
55. Gu F, Hu C, Tai Z, Yao C, Tian J, Zhang L, Xia Q, Gong C, Gao Y, Gao S. Tumour microenvironment-responsive lipoic acid nanoparticles for targeted delivery of docetaxel to lungcancer. *SCI REP-UK*. 2016 Nov 2;6:36281.
56. Mao H, Tao T, Song D, Liu M, Wang X, Liu X, Shi D. Zedoaronidol Inhibits Platelet-Derived Growth Factor-Induced Vascular Smooth Muscle CellsProliferation via Regulating AMP-Activated Protein Kinase Signaling Pathway. *CELL PHYSIOL BIOCHEM*. 2016;40(6):1506-1520.
57. Zhu JJ, Zhang XX, Miao YQ, He SF, Tian DM, Yao XS, Tang JS, Gan Y. Delivery of acetylthetevin B, an antitumor cardiac glycoside, using polymeric micelles for enhancedtherapeutic efficacy against lung cancer cells. *Acta Pharmacol Sin*. 2017 Feb;38(2):290-300.
58. Zu G, Yao J, Ji A, Ning S, Luo F, Li Z, Feng D, Rui Y, Li Y, Wang G, Tian X. Nurrl promotes intestinal regeneration after ischemia/reperfusion injury by inhibiting the expression ofp21 (Waf1/Cip1). *J MOL MED*. 2017 Jan;95(1):83-95.
59. Xiong N, Li S, Tang K, Bai H, Peng Y, Yang H, Wu C, Liu Y. Involvement of caveolin-1 in low shear stress-induced breast cancer cell motility and adhesion: Roles ofFAK/Src and ROCK/p-MLC pathways. *BBA-BIOMEMBRANES*. 2017 Jan;1864(1):12-22.
60. Qiao L, Liu X, Tang Y, Zhao Z, Zhang J, Feng Y. Down regulation ofthe long non-coding RNA PCAT-1 induced growth arrest and apoptosis of colorectal cancer cells. *Life Sci*. 2017 Nov 1;188:37-44.
61. Mei Y, Peng C, Liu YB, Wang J, Zhou HH. Silencing RIF1 decreases cell growth, migration and increases cisplatin sensitivity of humancervical cancer cells. *ONCOTARGET*. 2017 Nov 6;8(63):107044-107051.
62. Wu Y, Bu F, Yu H, Li W, Huang C, Meng X, Zhang L, Ma T, Li J. Methylation of Septin9 mediated by DNMT3a enhances hepatic stellate cells activation and liverfibrogenesis. *TOXICOL APPL PHARM*. 2017 Jan 15;315:35-49.
63. Xia S, Wang C, Ni X, Ni Z, Dong Y, Zhan W. NONHSAT076754 aids ultrasonography in predicting lymph node metastasis and promotesmigration and invasion of papillary thyroid cancer cells. *ONCOTARGET*. 2017 Jan 10;8(2):2293-2306.
64. Zu G, Yao J, Ji A, Ning S, Luo F, Li Z, Feng D, Rui Y, Li Y, Wang G, Tian X. Nurrl promotes intestinal regeneration after ischemia/reperfusion injury by inhibiting the expression of p21 (Waf1/Cip1). *J MOL MED*. 2017 Jan;95(1):83-95.
65. Yang Y, Du J, Liu F, Wang X, Li X, Li Y. Role of caspase-3/E-cadherin in helicobacter pylori-induced apoptosis of gastric epithelial cells. *ONCOTARGET*. 2017 Jul 22;8(35):59204-59216.
66. Mu W, Wang Z, Ma C, Jiang Y, Zhang N, Hu K, Li L, Wang Z. Metformin promotes the proliferation and differentiation of murine preosteoblast by regulating the expression of sirt6 and oct4. *Pharmacol Res*. 2017 Nov 21. pii: S1043-6618(17)30545-5.
67. Yang Y, Xia S, Ni X, Ni Z, Zhang L, Wang W, Kong Y, Wang Y, Ye L, Zhan W. MiR-324-5p assists ultrasonography in predicting lymph node metastasis of unifocal papillarythyroid microcarcinoma without extracapsular spread. *ONCOTARGET*. 2017 Jul 31;8(48):83802-83816.
68. Zou M, Li Y, Xia S, Chu Q, Xiao X, Qiu H, Chen Y, Zheng Z, Liu F, Zhuang L, Yu S. Knockdown of CAVEOLIN-1 Sensitizes Human Basal-Like Triple-Negative Breast Cancer Cells to Radiation. *CELL PHYSIOL BIOCHEM*. 2017;44(2):778-791.
69. Zhao W, Peng F, Shu M, Liu H, Hou X, Wang X, Ye J, Zhao B, Wang K, Zhong C, Xue L, Gao M, Liu Y, Zhao S. Isogambogenic Acid Inhibits the Growth of Glioma Through Activation of the AMPK-mTORPathway. *CELL PHYSIOL BIOCHEM*. 2017;44(4):1381-1395.
70. Li L, Wang P, Hu K, Wang X, Cai W, Ai C, Liu S, Wang Z. PFMG1 promotes osteoblast differentiation and prevents osteoporotic bone loss. *FASEB J*. 2017 Oct 11. pii: fj.201700422R.
71. Zhao Y, Jia Z, Huang S, Wu Y, Liu L, Lin L, Wang D, He Q, Ruan D. Age-Related Changes in Nucleus Pulposus Mesenchymal Stem Cells: An In Vitro Study in Rats. *Stem Cells Int*. 2017;2017:6761572.
72. Yang X, Sun J, Xia D, Can X, Liu L, Zhang J, Xu H, Du N, Liu W, Shen F, Zhang Z, Sun Y, Xi X. Capn4 Enhances Osteopontin Expression through Activation ofthe Wnt/ β -Catenin Pathway to Promote Epithelial Ovarian Carcinoma Metastasis. *CELL PHYSIOL BIOCHEM*. 2017;42(1):185-197.
73. Yao Q, Cai G, Yu Q, Shen J, Gu Z, Chen J, Shi W, Shi J. IDH1 mutation diminishes aggressive phenotype in glioma stem cells. *Int J Oncol*. 2018 Jan;52(1):270-278.
74. Yang D, Luo W, Wang J, Zheng M, Liao XH, Zhang N, Lu W, Wang L, Chen AZ, Wu WG, Liu H, Wang SB, Zhou XZ, Lu KP. A novel controlled release formulation of the Pin1 inhibitor ATRA to improve liver cancer therapy by simultaneously blocking multiple cancer pathways. *J Control Release*. 2018 Jan 10;269:405-422.
75. Feng W, Wang C, Liang C, Yang H, Chen D, Yu X, Zhao W, Geng D, Li S, Chen Z, Sun M The Dysregulated Expression of KCNQ10T1 and Its Interaction with Downstream Factors miR-145/CNNE2 in Breast Cancer Cells. *CELL PHYSIOL BIOCHEM*. 2018;49(2):432-446.
76. Wang Y, Li J, Kuang D, Wang X, Zhu Y, Xu S, Chen Y, Cheng H, Zhao Q, Duan Y, Wang G miR-148b-3p functions as a tumor suppressor in GISTs by directly targeting KIT. *Cell Commun Signal*. 2018 Apr 16;16(1):16.
77. Li L, Wang Y, Zhang X, Song G, Guo Q, Zhang Z, Diao Y, Yin H, Liu H, Jiang G Deubiquitinase USP48 promotes ATRA-induced granulocytic differentiation of acute promyelocytic leukemia cells. *Int J Oncol*. 2018 Aug;53(2):895-903.
78. Chen Q, Hang Y, Zhang T, Tan L, Li S, Jin Y USP10 promotes proliferation and migration and inhibits apoptosis of endometrial stromal cells in endometriosis through activating the Raf-1/MEK/ERK pathway. *AM J PHYSIOL-CELL PH*. 2018 Dec 1;315(6):C863-C872.
79. Zhang H, Nan W, Wang S, Si H, Li G Balance between fibroblast growth factor 10 and secreted frizzled-relate protein-1 controls the development of hair follicle by competitively regulating β -catenin signaling. *Biomed Pharmacother*. 2018 Jul;103:1531-1537.
80. Chen Y, Zhuang H, Chen X, Shi Z, Wang X Spermidine-induced growth inhibition and apoptosis via autophagic activation in cervical cancer. *Oncol Rep*. 2018 Jun;39(6):2845-2854.
81. Cui D, Zhao Y, Xu J Activation of CXCL5-CXCR2 axis promotes proliferation and accelerates G1 to S phase transition of papillary thyroid carcinoma cells and activates JNK and p38 pathways. *Cancer Biol Ther*. 2018 Nov 7:1-9.
82. Luo Z, Xu X, Sho T, Zhang J, Xu W, Yao J, Xu J ROS-induced autophagy regulates porcine trophoblast cell apoptosis, proliferation, and differentiation. *AM J PHYSIOL-CELL PH*. 2019 Feb 1;316(2):C198-C209.
83. Xiong Y, Li T, Assani G, Ling H, Zhou Q, Zeng Y, Zhou F, Zhou Y Ribociclib, a selective cyclin D kinase 4/6 inhibitor, inhibits proliferation and induces apoptosis of human cervical cancer in vitro and in vivo. *Biomed Pharmacother*. . 112:108602.
84. Zhou GZ, Shi YY, Wei LL, Sun GC Autophagy induction and antiproliferative effect of a novel curcumin derivative MOMI-1 on the human lung cancer cells A549. *J BIOCHEM MOL TOXIC*. . 33(4):e22280.
85. Wang W, Wu J, Zhang Q, Li X, Zhu X, Wang Q, Cao S, Du L Mitochondria-mediated apoptosis was induced by oleuropein in H1299 cells involving activation of p38 MAP kinase. *J Cell Biochem*. . 120(4):5480-5494.
86. Tu W, Dong C, Fu J, Pan Y, Kobayashi A, Furusawa Y, Konishi T, Shao C Both

- irradiated and bystander effects link with DNA repair capacity and the linear energy transfer. *Life Sci.* . 222:228-234.
87. Zhang Y, Li Z, Zhao W, Hu H, Zhao L, Zhu Y, Yang X, Gao B, Yang H, Huang Y, Song X WD repeat and SOCS box containing protein 2 in the proliferation, cycle progression, and migration of melanoma cells. *Biomed Pharmacother.* . 116:108974.
 88. Li Z, Chen S, Chen S, Huang D, Ma K, Shao Z Moderate activation of Wnt/ β -catenin signaling promotes the survival of rat nucleus pulposus cells via regulating apoptosis, autophagy, and senescence. *J Cell Biochem.* . 120(8):12519-12533.
 89. Cheng S, Li X, Jia Z, Lin L, Ying J, Wen T, Zhao Y, Guo Z, Zhao X, Li D, Ji W, Wang D, Ruan D The inflammatory cytokine TNF- α regulates the biological behavior of rat nucleus pulposus mesenchymal stem cells through the NF- κ B signaling pathway in vitro. *J Cell Biochem.* . 120(8):13664-13679.
 90. Wang X, Ma S, Liu Y, Lu W, Sun L, Zhao P, Xia Q Transcriptional repression of endogenous genes in BmE cells using CRISPRi system. *INSECT BIOCHEM MOLEC.* . 111:103172.
 91. Gao S, Chen S, Chen L, Zhao Y, Sun L, Cao M, Huang Y, Niu Q, Wang F, Yuan C, Li C, Zhou X Brain-derived neurotrophic factor: A steroidogenic regulator of Leydig cells. *J Cell Physiol.* . 234(8):14058-14067.
 92. Wang J, Lei B, Yan J, Li J, Zhou X, Ren F, Guo H Donkey milk oligosaccharides influence the growth-related characteristics of intestinal cells and induce G2/M growth arrest via the p38 pathway in HT-29 cells. *Food Funct.* . 10(8):4823-4833.
 93. Jia H, Xu M, Bo Y, Li W, Zhou R Ras-ERK1/2 signaling accelerates the progression of colorectal cancer via mediation of H2BK5ac. *Life Sci.* . 230:89-96.
 94. Li XY, Wu SY, Leung PS Human Fetal Bone Marrow-Derived Mesenchymal Stem Cells Promote the Proliferation and Differentiation of Pancreatic Progenitor Cells and the Engraftment Function of Islet-Like Cell Clusters. *Int J Mol Sci.* . 20(17). pii: E4083.
 95. Wang W, Xia S, Zhan W The Long Non-Coding RNA ENST00000489676 Influences Papillary Thyroid Cancer Cell Proliferation and Invasion through Regulating MiR-922. *J Cancer.* . 10(22):5434-5446.
 96. Xu Z, Zheng X, Xia X, Wang X, Luo N, Huang B, Pan X 17 β -estradiol at low concentrations attenuates the efficacy of tamoxifen in breast cancer therapy. *Environ Pollut.* . 255(Pt 2):113228.
 97. Zhang S, Xu J, He Z, Xue F, Jiang T, Xu M Sodium Selenate Ameliorates Cardiac Injury Developed from High-Fat Diet in Mice through Regulation of Autophagy Activity. *SCI REP-UK.* . 9(1):18752.
 98. Wu X, Kong W, Qi X, Wang S, Chen Y, Zhao Z, Wang W, Lin X, Lai J, Yu Z, Lai G Icaritin induces apoptosis of human lung adenocarcinoma cells by activating the mitochondrial apoptotic pathway. *Life Sci.* . 239:116879.
 99. Zhang HT, Sun J, Yan Y, Cui SH, Wang H, Wang CH, Qiu C, Chen X, Ding JS, Qian HG, Wang JC, Zhang Q Encapsulated microRNA by gemcitabine prodrug for cancer treatment. *J Control Release.* . 316:317-330.
 100. Xie K, Xu C, Zhang M, Wang M, Min L, Qian C, Wang Q, Ni Z, Mou S, Dai H, Pang H, Gu L Yes-associated protein regulates podocyte cell cycle re-entry and dedifferentiation in adriamycin-induced nephropathy. *Cell Death Dis.* . 10(12):915.
 101. You Y, Bi FF, Jiang Y, Xu YT, An YY, Li D, Yang Q BRCA1 affects the resistance and stemness of SKOV3-derived ovarian cancer stem cells by regulating autophagy. *CANCER MED-US.* . 8(2):656-668.
 102. Guan L, Chen Y, Wang Y, Zhang H, Fan S, Gao Y, Jiao T, Fu K, Sun J, Yu A, Huang M, Bi H Effects of carnitine palmitoyltransferases on cancer cellular senescence. *J Cell Physiol.* . 234(2):1707-1719.
 103. Bu JY, Lv WZ, Liao YF, Xiao XY, Lv BJ Long non-coding RNA LINC00978 promotes cell proliferation and tumorigenesis via regulating microRNA-497/NTRK3 axis in gastric cancer. *Int J Biol Macromol.* . 123:1106-1114.
 104. Ji AL, Li T, Zu G, Feng DC, Li Y, Wang GZ, Yao JH, Tian XF Ubiquitin-specific protease 22 enhances intestinal cell proliferation and tissue regeneration after intestinal ischemia reperfusion injury. *WORLD J GASTROENTERO.* . 25(7):824-836.
 105. Xie WY, He RH, Zhang J, He YJ, Wan Z, Zhou CF, Tang YJ, Li Z, Mcleod HL, Liu J β -blockers inhibit the viability of breast cancer cells by regulating the ERK/COX-2 signaling pathway and the drug response is affected by ADRB2 single-nucleotide polymorphisms. *Oncol Rep.* . 41(1):341-350.
 106. Li H, Wu X, Yang B, Li J, Xu L, Liu H, Li S, Xu J, Yang M, Wei M Evaluation of biomimetically synthesized mesoporous silica nanoparticles as drug carriers: Structure, wettability, degradation, biocompatibility and brain distribution. *MAT SCI ENG C-MATER.* . 94:453-464.
 107. Hu X, Jiang J, Ni C, Xu Q, Ye S, Wu J, Ge F, Han Y, Mo Y, Huang D, Yang L HBV Integration-mediated Cell Apoptosis in HepG2.2.15. *J Cancer.* . 10(17):4142-4150.
 108. Liu X, Sun J, Yuan P, Shou K, Zhou Y, Gao W, She J, Hu J, Yang J, Yang J Mfn2 inhibits proliferation and cell-cycle in HeLa cells via Ras-NF- κ B signal pathway. *Cancer Cell Int.* . 19:197.
 109. Zhu D, Zhou J, Zhao J, Jiang G, Zhang X, Zhang Y, Dong M ZC3H13 suppresses colorectal cancer proliferation and invasion via inactivating Ras-ERK signaling. *J Cell Physiol.* . 234(6):8899-8907.
 110. Li Y, Ren Y, Wang Y, Tan Y, Wang Q, Cai J, Zhou J, Yang C, Zhao K, Yi K, Jin W, Wang L, Liu M, Yang J, Li M, Kang C A Compound AC1Q3QWB Selectively Disrupts HOTAIR-Mediated Recruitment of PRC2 and Enhances Cancer Therapy of DZNeP. *Theranostics.* . 9(16):4608-4623.
 111. Yu HX, Yao Y, Bu FT, Chen Y, Wu YT, Yang Y, Chen X, Zhu Y, Wang Q, Pan XY, Meng XM, Huang C, Li J Blockade of YAP alleviates hepatic fibrosis through accelerating apoptosis and reversion of activated hepatic stellate cells. *Mol Immunol.* . 107:29-40.
 112. Li L, Zhang C, Chen JL, Hong FF, Chen P, Wang JF Effects of simulated microgravity on the expression profiles of RNA during osteogenic differentiation of human bone marrow mesenchymal stem cells. *CELL PROLIFERAT.* . 52(2):e12539.
 113. Wang Y, Li X, Gong X, Zhao Y, Wu J MicroRNA-322 Regulates Self-renewal of Mouse Spermatogonial Stem Cells through Rassf8. *Int J Biol Sci.* . 15(4):857-869.
 114. Wei X, Guo J, Li Q, Jia Q, Jing Q, Li Y, Zhou B, Chen J, Gao S, Zhang X, Jia M, Niu C, Yang W, Zhi X, Wang X, Yu D, Bai L, Wang L, Na J, Zou Y, Zhang J, Zhang S, Meng D Bcl1 regulates self-renewal and impedes mesodermal differentiation of human embryonic stem cells. *Sci Adv.* . 5(3):eaau7887.
 115. Huang JF, Jiang HY, Cai H, Liu Y, Zhu YQ, Lin SS, Hu TT, Wang TT, Yang WJ, Xiao B, Sun SH, Ma LY, Yin HR, Wang F Genome-wide screening identifies oncofetal lncRNA Ptn-dt promoting the proliferation of hepatocellular carcinoma cells by regulating the Ptn receptor. *Oncogene.* . 38(18):3428-3445.
 116. Liu Q, Yang X, Ji J, Zhang SL, He Y Novel nannocystin A analogues as anticancer therapeutics: Synthesis, biological evaluations and structure-activity relationship studies. *Eur J Med Chem.* . 170:99-111.
 117. Deng Z, Chai J, Zeng Q, Zhang B, Ye T, Chen X, Xu X The anticancer properties and mechanism of action of tablysin-15, the RGD-containing disintegrin, in breast cancer cells. *Int J Biol Macromol.* . 129:1155-1167.
 118. Sun H, Hong M, Yang Q, Li C, Zhang G, Yue Q, Ma Y, Li X, Li CZ Visualizing the down-regulation of hTERT mRNA expression using gold-nanoflare probes and verifying the correlation with cancer cell apoptosis. *Analyst.* . 144(9):2994-3004.
 119. Zhang H, Shi Q, Nan W, Wang Y, Wang S, Yang F, Li G Ginkgolide B and bilobalide promote the growth and increase β -catenin expression in hair follicle dermal papilla cells of American minks. *Biofactors.* . 45(6):950-958.
 120. Zhang M, Jin C, Yang Y, Wang K, Zhou Y, Zhou Y, Wang R, Li T, Hu R AIM2 promotes non-small-cell lung cancer cell growth through inflammasome-dependent pathway. *J Cell Physiol.* . 234(11):20161-20173.
 121. Zhang L, Tao X, Fu Q, Ge C, Li R, Li Z, Zhu Y, Tian H, Li Q, Liu M, Hu H, Zeng B, Lin Z, Li C, Luo R, Song X Curcumin inhibits cell proliferation and migration in NSCLC through a synergistic effect on the TLR4/MyD88 and EGFR pathways. *Oncol Rep.* . 42(5):1843-1855.
 122. Jin L, Wang Q, Chen J, Wang Z, Xin H, Zhang D Efficient Delivery of Therapeutic siRNA by Fe3O4 Magnetic Nanoparticles into Oral Cancer Cells. *Pharmaceutics.* .

- 11(11). pii: E615.
123. Du XS, Li HD, Yang XJ, Li JJ, Xu JJ, Chen Y, Xu QQ, Yang L, He CS, Huang C, Meng XM, Li J Wogonin attenuates liver fibrosis via regulating hepatic stellate cell activation and apoptosis. *Int Immunopharmacol.* . 75:105671.
 124. Zhu HQ, Zhang C, Guo ZY, Yang JM, Guo JH, Chen C, Yao QH, Liu F, Zhang QW, Gao FH Oridonin induces Mdm2-p60 to promote p53-mediated apoptosis and cell cycle arrest in neuroblastoma. *CANCER MED-US.* . 8(11):5313-5326.
 125. Zhang L, Liao C, Yang Y, Wang YZ, Ding K, Huo D, Hou C Response of lipid biosynthesis in *Chlorella pyrenoidosa* to intracellular reactive oxygen species level under stress conditions. *BIORESOURTE TECHNOL.* . 287:121414.
 126. Hou L, Wang R, Wei H, Li S, Liu L, Lu X, Yu H, Liu Z ABT737 enhances ovarian cancer cells sensitivity to cisplatin through regulation of mitochondrial fission via Sirt3 activation. *Life Sci.* . 232:116561.
 127. Zhang Y, Xiao JF, Yang HF, Jiao Y, Cao WW, Shi HM, Cun JF, Tay FR, Ping J, Xiao YH N-Acetyl Cysteine as a Novel Polymethyl Methacrylate Resin Component: Protection against Cell Apoptosis and Genotoxicity. *Oxid Med Cell Longev.* . 2019:1301736.
 128. Zhong B, Shi D, Wu F, Wang S, Hu H, Cheng C, Qing X, Huang X, Luo X, Zhang Z, Shao Z Dynasore suppresses cell proliferation, migration, and invasion and enhances the antitumor capacity of cisplatin via STAT3 pathway in osteosarcoma. *Cell Death Dis.* . 10(10):687.
 129. Jing X, Xu Y, Liu D, Wu Y, Zhou N, Wang D, Yan K, Meng L Intelligent nanoflowers: a full tumor microenvironment-responsive multimodal cancer theranostic nanoplatform. *Nanoscale.* . 11(33):15508-15518.
 130. Wei X, Yang X, Wang B, Yang Y, Fang Z, Yi C, Shi L, Song D LncRNA MBNL1-AS1 represses cell proliferation and enhances cell apoptosis via targeting miR-135a-5p/PHLPP2/FOXO1 axis in bladder cancer. *CANCER MED-US.* . 9(2):724-736.
 131. Chen Y, Zhang X, An Y, Liu B, Lu M LncRNA HCP5 promotes cell proliferation and inhibits apoptosis via miR-27a-3p/IGF-1 axis in human granulosa-like tumor cell line KGN. *Mol Cell Endocrinol.* . 503:110697.
 132. Li Y, Cao Y, Wang J, Fu S, Cheng J, Ma L, Zhang Q, Guo W, Kan X, Liu J Kp-10 promotes bovine mammary epithelial cell proliferation by activating GPR54 and its downstream signaling pathways. *J Cell Physiol.* . 235(5):4481-4493.
 133. Yan Chen, Ning Li, Haijing Wang, Ningning Wang, Hui Peng, Jing Wang, Yihong Li, Mingdi Liu, Hui Li, Yu Zhang, Zhaohui Wang Amentoflavone suppresses cell proliferation and induces cell death through triggering autophagy-dependent ferroptosis in human glioma *Life Sci.* . doi: 10.1016/j.lfs.2020.117425.
 134. Ruiyun Wu, Nan Shang, Meng Gui, Jian Yin, Pinglan Li Sturgeon (Acipenser)-Derived Chondroitin Sulfate Suppresses Human Colon Cancer HCT-116 Both In Vitro and In Vivo by Inhibiting Proliferation and Inducing Apoptosis *Nutrients.* . doi: 10.3390/nu12041130.
 135. Xu Bai, Yu Wang, Bo Hu, Qi Cao, Maochen Xing, Shuliang Song, Aiguo Ji Fucoindan Induces Apoptosis of HT-29 Cells via the Activation of DR4 and Mitochondrial Pathway *Mar Drugs.* . doi: 10.3390/md18040220.
 136. Chen Cai, Wei Dang, Shilei Liu, Ling Huang, Yang Li, Guoqiang Li, Siyuan Yan, Chengkai Jiang, Xiaoling Song, Yunping Hu, Jun Gu Anthrax toxin receptor 1/tumor endothelial marker 8 promotes gastric cancer progression through activation of the PI3K/AKT/mTOR signaling pathway *Cancer Sci.* . doi: 10.1111/cas.14326.
 137. Zhen-Yu She, Yue-Ling Li, Yang Lin, Ming-Hui Lu, Ya-Lan Wei, Kai-Wei Yu, Ning Zhong, Yu Xiao Kinesin-6 family motor KIF20A regulates central spindle assembly and acrosome biogenesis in mouse spermatogenesis *BBA-MOL CELL RES.* . doi: 10.1016/j.bbamcr.2019.118636.
 138. Hui-Hui Zhang, Zhong-Xin Huang, Su-Quan Zhong, Kui-Lin Fei, You-Han Cao miR-21 inhibits autophagy and promotes malignant development in the bladder cancer T24 cell line *Int J Oncol.* . doi: 10.3892/ijo.2020.4984.
 139. Wei-Feng Wu, Jia-Nan Wang, Zeng Li, Biao Wei, Juan Jin, Li Gao, Hai-Di Li, Jun Li, Hai-Yong Chen, Xiao-Ming Meng 7-Hydroxycoumarin protects against cisplatin-induced acute kidney injury by inhibiting necroptosis and promoting Sox9-mediated tubular epithelial cell proliferation *Phytomedicine.* . doi: 10.1016/j.phymed.2020.153202.
 140. Fujun Jin, Yexuan Zhu, Jingyi Chen, Rongze Wang, Yiliang Wang, Yanting Wu, Pengjun Zhou, Xiaowei Song, Zhe Ren, Jun Dong BRE Promotes Esophageal Squamous Cell Carcinoma Growth by Activating AKT Signaling *Front Oncol.* . doi: 10.3389/fonc.2020.01407.
 141. Jing-Yu Liu, Xiao-Xiao Hou, Zhuo-Yu Li, Shu-Hua Shan, Ming-Chang Chang, Cui-Ping Feng, Yin Wei Isolation and structural characterization of a novel polysaccharide from *Hericium erinaceus* fruiting bodies and its arrest of cell cycle at S-phase in colon cancer cells *Int J Biol Macromol.* . doi: 10.1016/j.ijbiomac.2020.04.162.
 142. Suresh Veeraperumal, Hua-Mai Qiu, Shan-Shan Zeng, Wan-Zhi Yao, Bao-Peng Wang, Yang Liu, Kit-Leong Cheong Polysaccharides from *Gracilaria lemaneiformis* promote the HaCaT keratinocytes wound healing by polarised and directional cell migration *CARBOHYD POLYM.* . doi: 10.1016/j.carbpol.2020.116310.
 143. Qingye Zeng, Baishuang Zeng, Jinwei Chai, Jiena Wu, Ruiyin Guo, Yahua Gao, Xiaoyan Han, Jie Yang, Michail Kotsyfakis, Xueqing Xu Antioxidant properties and neuroprotective effects of Esc-1GN through the regulation of MAPK and AKT signaling *Life Sci.* . doi: 10.1016/j.lfs.2020.117753.
 144. Miaomiao Chai, Ce Gu, Qihua Shen, Jiaying Liu, Yi Zhou, Ziyang Jin, Wanli Xiong, Yan Zhou, Wensong Tan Hypoxia alleviates dexamethasone-induced inhibition of angiogenesis in cocultures of HUVECs and rBMSCs via HIF-1 α Stem Cell Res Ther. . doi: 10.1186/s13287-020-01853-x.
 145. Dingsheng Liu, Hong Zhang, Mingming Cui, Chunsheng Chen, Yong Feng Hsa-miR-425-5p promotes tumor growth and metastasis by activating the CTNND1-mediated β -catenin pathway and EMT in colorectal cancer *Cell Cycle.* . doi: 10.1080/15384101.2020.1783058.
 146. Xiaoduan Li, Meiling Tang Exosomes released from M2 macrophages transfer miR-221-3p contributed to EOC progression through targeting CDKN1B *CANCER MED-US.* . doi: 10.1002/cam4.3252.
 147. Jingjing Ma, Yunting Zhang, Huan Deng, Yinghui Liu, Xiaofei Lei, Pengzhan He, Weiguo Dong Thymoquinone inhibits the proliferation and invasion of esophageal cancer cells by disrupting the AKT/GSK-3 β /Wnt signaling pathway via PTEN upregulation *Phytother Res.* . doi: 10.1002/ptr.6795.
 148. Jiayu Gao, Yu Qin, Keke Luo, Xilan Wang, Chun Yu, Aihua Zhang, Xueli Pan Downregulation of miR-4755-5p promotes fluoride-induced osteoblast activation via targeting Cyclin D1 *J TRACE ELEM MED BIO.* . doi: 10.1016/j.jtemb.2020.126626.
 149. Guangyao Zhou, Chunxue Li, Yating Zhan, Rongrong Zhang, Boyu Lv, Wujun Geng, Jianjian Zheng Pinostilbene hydrate suppresses hepatic stellate cell activation via inhibition of miR-17-5p-mediated Wnt/ β -catenin pathway *Phytomedicine.* . doi: 10.1016/j.phymed.2020.153321.
 150. Fang Yuan, Cong Cheng, Feiyan Xiao, Hongcai Liu, Shan Cao, Gan Zhou Inhibition of mTORC1/P70S6K pathway by Metformin synergistically sensitizes Acute Myeloid Leukemia to Ara-C *Life Sci.* . doi: 10.1016/j.lfs.2020.117276.
 151. Yi Qu, Ying He, Yang Yang, Shaoqing Li, Wei An, Zhilin Li, Xue Wang, Zhengxue Han, Lizheng Qin ALDH3A1 acts as a prognostic biomarker and inhibits the epithelial mesenchymal transition of oral squamous cell carcinoma through IL-6/STAT3 signaling pathway *J Cancer.* . doi: 10.7150/jca.40171.
 152. Lin Hao, Kun Pang, Hui Pang, Junjie Zhang, Zhiguo Zhang, Houguang He, Rongsheng Zhou, Zhenduo Shi, Conghui Han Knockdown of P3H4 inhibits proliferation and invasion of bladder cancer *AGING-US.* . doi: 10.18632/aging.102732.
 153. Xiaohui Tian, Yuanyuan Wu, Yating Yang, Jiaying Wang, Menglan Niu, Shanjun Gao, Tao Qin, Dengke Bao Long noncoding RNA LINC00662 promotes M2 macrophage polarization and hepatocellular carcinoma progression via activating Wnt/ β -catenin signaling *Mol Oncol.* . doi: 10.1002/1878-0261.12606.
 154. Ya-Sai Sun, Kiran Thakur, Fei Hu, Jian-Guo Zhang, Zhao-Jun Wei Icariside II inhibits tumorigenesis via inhibiting AKT/Cyclin E/ CDK 2 pathway and activating mitochondria-dependent pathway *Pharmacol Res.* . doi: 10.1016/j.phymed.2020.153202.

- 10.1016/j.phrs.2019.104616.
155. Jun Gao, Chao Dai, Xin Yu, Xiang-Bao Yin, Fan Zhou microRNA-485-5p inhibits the progression of hepatocellular carcinoma through blocking the WBP2/Wnt signaling pathway *Cell Signal.* . doi: 10.1016/j.cellsig.2019.109466.
 156. Shidi Wu, Mengjie Zhao, Yanhong Sun, Meng Xie, Kehao Le, Ming Xu, Changzheng Huang The potential of Diosgenin in treating psoriasis: Studies from HaCaT keratinocytes and imiquimod-induced murine model *Life Sci.* . doi: 10.1016/j.lfs.2019.117115.
 157. Yake Huang, Jiani Huang, Yan Huang, Lei Gan, Ling Long, Aimin Pu, Rongkai Xie TFRC promotes epithelial ovarian cancer cell proliferation and metastasis via up-regulation of AXIN2 expression *Am J Cancer Res.* .
 158. Chuanyong Wu, Xiao-Ting Zhu, Lei Xia, Lin Wang, Wenjun Yu, Qiaomei Guo, Mingna Zhao, Jiatao Lou High Expression of Long Noncoding RNA PCNA-AS1 Promotes Non-Small-Cell Lung Cancer Cell Proliferation and Oncogenic Activity via Upregulating CCND1 *J Cancer.* . doi: 10.7150/jca.39087.
 159. Pei Liu, Chenchen Zhu, Juanjuan Luo, Sheng Lan, Dongsheng Su, Qiongjin Wang, Zhe Wei, Wei Cui, Chuan Xu, Xiaojun Yang Par6 regulates cell cycle progression through enhancement of Akt/PI3K/GSK-3 β signaling pathway activation in glioma *FASEB J.* . doi: 10.1096/fj.201901629RR.
 160. Hua Su, Chen Ye, Chun-Tao Lei, Hui Tang, Jie-Yu Zeng, Fan Yi, Chun Zhang Subcellular trafficking of tubular MDM2 implicates in acute kidney injury to chronic kidney disease transition during multiple low-dose cisplatin exposure *FASEB J.* . doi: 10.1096/fj.201901412R.
 161. Fei Dai, Fei Luo, Rui Zhou, Qiang Zhou, Jianzhong Xu, Zehua Zhang, Jun Xiao, Lei Song Calponin 3 is associated with poor prognosis and regulates proliferation and metastasis in osteosarcoma *AGING-US.* . doi: 10.18632/aging.103224.
 162. Ping Xiang, Kun Wang, Jue Bi, Mengying Li, Rui-Wen He, Daolei Cui, Lena Q Ma Organic extract of indoor dust induces estrogen-like effects in human breast cancer cells *Sci Total Environ.* . doi: 10.1016/j.scitotenv.2020.138505.
 163. Chen Peng, Yansong Yang, Li Ji, Panpan Yang, Xiaoqing Yang, Yuquan Zhang Cyclin H predicts the poor prognosis and promotes the proliferation of ovarian cancer *Cancer Cell Int.* . doi: 10.1186/s12935-020-01406-5.
 164. Peng Xu, Yu Xin, Zheng Zhang, Xiangyu Zou, Ke Xue, Huizhong Zhang, Wenjie Zhang, Kai Liu Extracellular vesicles from adipose-derived stem cells ameliorate ultraviolet B-induced skin photoaging by attenuating reactive oxygen species production and inflammation *Stem Cell Res Ther.* . doi: 10.1186/s13287-020-01777-6.
 165. Zhuo Wang, Weixiao Nan, Huazhe Si, Shiyong Wang, Haihua Zhang, Guangyu Li Pantothenic acid promotes dermal papilla cell proliferation in hair follicles of American minks via inhibitor of DNA Binding 3/Notch signaling pathway *Life Sci.* . doi: 10.1016/j.lfs.2020.117667.
 166. Yuelai Zhou, Xiaolei Li, Huilin Yang LINC00612 functions as a ceRNA for miR-214-5p to promote the proliferation and invasion of osteosarcoma in vitro and in vivo *Exp Cell Res.* . doi: 10.1016/j.yexcr.2020.112012.
 167. Xiangru Zheng, Jinlai Wei, Wenjun Li, Xiaoli Li, Wuyi Wang, Jinbao Guo, Zhongxue Fu PRDX2 removal inhibits the cell cycle and autophagy in colorectal cancer cells *AGING-US.* 2019 Apr. doi: 10.18632/aging.103690.
 168. Kai Zheng, Jingxin Ma, Yifei Wang, Zhendan He, Kejun Deng Sulforaphane Inhibits Autophagy and Induces Exosome-Mediated Paracrine Senescence via Regulating mTOR/TFE3 *Mol Nutr Food Res.* 2019 Apr. doi: 10.1002/mnfr.201901231.
 169. Da Liu, Zixuan Song, Xiaoying Wang, Ling Ouyang Ubiquitin C-Terminal Hydrolase L5 (UCHL5) Accelerates the Growth of Endometrial Cancer via Activating the Wnt/ β -Catenin Signaling Pathway *Front Oncol.* 2019 Apr 1. doi: 10.3389/fonc.2020.00865.
 170. Shangling Zhu, Yuanmei Ye, Yiming Shi, Junlong Dang, Xiaoxue Feng, Yingdi Chen, Fang Liu, Nancy Olsen, Jianlin Huang, Song Guo Zheng Sonic Hedgehog Regulates Proliferation, Migration and Invasion of Synoviocytes in Rheumatoid Arthritis via JNK Signaling *Front Immunol.* 2019 Aug. doi: 10.3389/fimmu.2020.01300.
 171. Shan Li, Fangmin Song, Xu Lei, Jingtao Li, Fang Li, Huabing Tan hsa_circ_0004018 suppresses the progression of liver fibrosis through regulating the hsa-miR-660-3p/TEP1 axis *AGING-US.* 2019 Aug. doi: 10.18632/aging.103257.
 172. Yanhua Yang, Shujun Xia, Lu Zhang, Wenhan Wang, Lin Chen, Weiwei Zhan MiR-324-5p/PTPRD/CBEPD axis promotes papillary thyroid carcinoma progression via microenvironment alteration *Cancer Biol Ther.* 2019 Aug. doi: 10.1080/15384047.2020.1736465.
 173. Shengwu Yuan, Kongrui Zhu, Mei Ma, Xiaoshan Zhu, Kaifeng Rao, Zijian Wang In vitro oxidative stress, mitochondrial impairment and G1 phase cell cycle arrest induced by alkyl-phosphorus-containing flame retardants *Chemosphere.* 2019 Aug. doi: 10.1016/j.chemosphere.2020.126026.
 174. Tong Zhao, Mingyang Cai, Man Liu, Guangsong Su, Daniel An, Byoungsan Moon, Guochang Iyu, Yibo Si, Lingyi Chen, Wange Lu lncRNA 5430416N02Rik Promotes the Proliferation of Mouse Embryonic Stem Cells by Activating Mid1 Expression through 3D Chromatin Architecture *Stem Cell Reports.* 2019 Aug. doi: 10.1016/j.stemcr.2020.02.002.
 175. Xiaosong Wei, Beibei Wang, Qi Wang, Xiaoming Yang, Yang Yang, Zhiwei Fang, Chengzhi Yi, Lei Shi, Xin Fan, Jin Tao, Yufeng Guo, Dongkui Song MiR-362-5p, Which Is Regulated by Long Non-Coding RNA MBNL1-AS1, Promotes the Cell Proliferation and Tumor Growth of Bladder Cancer by Targeting QKI *Front Pharmacol.* 2019 Aug 1. doi: 10.3389/fphar.2020.00164.
 176. Cong Wang, Ke Wang, Shu-Fang Li, Su-Jing Song, Yan Du, Ruo-Wen Niu, Xue-Wen Qian, Xiao-Qing Peng, Fei-Hu Chen 4-Amino-2-trifluoromethyl-phenyl retinate induced differentiation of human myelodysplastic syndromes SKM-1 cell lines by up-regulating DDX23 *Biomed Pharmacother.* 2019 Aug 1. doi: 10.1016/j.biopha.2019.109736.
 177. Tongke Chen, Liqian Zhao, Shinuo Chen, Bin Zheng, Hong Chen, Tianni Zeng, Hanxiao Sun, Sijia Zhong, Wencan Wu, Xiaokun Lin, Lihua Wang The curcumin analogue WZ35 affects glycolysis inhibition of gastric cancer cells through ROS-YAP-JNK pathway *Food Chem Toxicol.* 2019 Aug 21. doi: 10.1016/j.fct.2020.111131.
 178. Gaigai Guo, Zhouyong Gao, Mengsha Tong, Dongdong Zhan, Guangshun Wang, Yi Wang, Jun Qin NQO1 is a determinant for cellular sensitivity to anti-tumor agent Napabucasin *Am J Cancer Res.* 2019 Aug 29.
 179. Jiawen Zhang, Qinyi Zhang, Cong Sun, Yanjuan Huang, Jing Zhang, Qingying Wang Clinical relevance of ARF/ARL family genes and oncogenic function of ARL4C in endometrial cancer *Biomed Pharmacother.* 2019 Dec. doi: 10.1016/j.biopha.2020.110000.
 180. Luzhang Ji, Xiangru Huo, Yuwen Zhang, Zixiang Yan, Qianfeng Wang, Bo Wen TOPORS, a tumor suppressor protein, contributes to the maintenance of higher-order chromatin architecture *BBA-GENE REGUL MECH.* 2019 Dec 10. doi: 10.1016/j.bbagr.2020.194518.
 181. Chunli Xu, Qingke Bai, Chen Wang, Qiuyu Meng, Yuming Gu, Qiwei Wang, Wenjie Xu, Ying Han, Yong Qin, Song Jia, Junfang Zhang, Jie Xu, Jiao Li, Miao Chen, Feng Wang miR-433 Inhibits Neuronal Growth and Promotes Autophagy in Mouse Hippocampal HT-22 Cell Line *Front Pharmacol.* 2019 Dec 15. doi: 10.3389/fphar.2020.536913.
 182. Jing Wang, Yiyi Xiong HSH2D contributes to methotrexate resistance in human T-cell acute lymphoblastic leukaemia *Oncol Rep.* 2019 Dec 28. doi: 10.3892/or.2020.7772.
 183. Zhiren Yao, Yu Li, Zhiwei Wang, Yuanyuan Lan, Tiexin Zeng, Hongfei Gong, Kaimei Zhu, Hui Tang, Shengjiu Gu Research on anti-hepatocellular carcinoma activity and mechanism of Polygala fallax Hems1 *J Ethnopharmacol.* 2019 Dec 30. doi: 10.1016/j.jep.2020.113062.
 184. Huifang Liu, Xiaowen Ma, Xin Wang, Xirui Ma, Ziming Mao, Jing Zhu, Fengling Chen Hsa_circ_0000345 regulates the cellular development of ASMCs in response to oxygenized low-density lipoprotein *J Cell Mol Med.* 2019 Dec 4. doi: 10.1111/jcmm.15801.
 185. Jiaming Gao, Yating Wu, Dajun He, Xiaoqi Zhu, Hongbin Li, Haifeng Liu, Hailiang Liu Anti-aging effects of Ribes meyeri anthocyanins on neural stem cells and aging mice *AGING-US.* 2019 Feb. doi: 10.18632/aging.103955.

186. Qing Zhou, Peiyu Jin, Jieyu Liu, Sihao Li, Weijue Liu, Shuhua Xi Arsenic-induced HER2 promotes proliferation, migration and angiogenesis of bladder epithelial cells via activation of multiple signaling pathways in vitro and in vivo *Sci Total Environ.* 2019 Feb. doi: 10.1016/j.scitotenv.2020.141962.
187. Guanting Li, Xianwei Wu, Peng Sun, Zhiyang Zhang, Enxian Shao, Jianping Mao, Hua Cao, Hongliang Huang Dithiolation indolizine exerts viability suppression effects on A549 cells via triggering intrinsic apoptotic pathways and inducing G2/M phase arrest *Biomed Pharmacother.* 2019 Feb 21. doi: 10.1016/j.biopha.2020.110961.
188. Junlong Zhong, Bin Wang, Biao Wu, Jie Shang, Ning Jiang, Aiyong Cui, Yufeng Chen, Xianghe Xu, Huading Lu Clock knockdown attenuated reactive oxygen species-mediated senescence of chondrocytes through restoring autophagic flux *Life Sci.* 2019 Jan. doi: 10.1016/j.lfs.2021.119036.
189. Haiyun Luo, Wenjing Liu, Yanli Zhang, Yeqing Yang, Xiao Jiang, Shiqing Wu, Longquan Shao METTL3-mediated m⁶A modification regulates cell cycle progression of dental pulp stem cells *Stem Cell Res Ther.* 2019 Jan 1. doi: 10.1186/s13287-021-02223-x.
190. Can Yang, Jian-Fei Hu, Qian Zhan, Zu-Wei Wang, Ge Li, Jing-Jing Pan, Long Huang, Cheng-Yu Liao, Yi Huang, Yi-Feng Tian, Bai-Yong Shen, Jiang-Zhi Chen, Yao-Dong Wang, Shi Chen SHCBP1 interacting with EOGT enhances O-GlcNAcylation of NOTCH1 and promotes the development of pancreatic cancer *Genomics.* 2019 Jul 10. doi: 10.1016/j.ygeno.2021.01.010.
191. Guangyan Kan, Ziyang Wang, Chunjie Sheng, Chen Yao, Yizhi Mao, Shuai Chen Inhibition of DKC1 induces telomere-related senescence and apoptosis in lung adenocarcinoma *J Transl Med.* 2019 Jul 29. doi: 10.1186/s12967-021-02827-0.
192. Ping Wang, Xinwei Hao, Xiaohan Li, Yizhi Yan, Wenxiu Tian, Lin Xiao, Zhenming Wang, Junhong Dong Curcumin inhibits adverse psychological stress-induced proliferation and invasion of glioma cells via down-regulating the ERK/MAPK pathway *J Cell Mol Med.* 2019 Jun. doi: 10.1111/jcmm.16749.
193. Li Liu, Yinliang Qi, Dake Huang, Yechuan Xu, Mingcong Li, Zhou Hong, Yanan Gu, Bo Jia, Xin Luo, Sumei Zhang, Shengquan Zhang Downregulated expression of IL-28RA is involved in the pathogenesis of pancreatic ductal adenocarcinoma *Int J Oncol.* 2019 Jun 24. doi: 10.3892/ijo.2021.5235.
194. Ruihua Jing, Conghui Hu, Tiantian Qi, Jiaqi Yue, Guangyan Wang, Ming Zhang, Chan Wen, Cheng Pei, Bo Ma FILIP1L-mediated cell apoptosis, epithelial-mesenchymal transition and extracellular matrix synthesis aggravate posterior capsular opacification *Life Sci.* 2019 Mar. doi: 10.1016/j.lfs.2021.120061.
195. Chao Wei, Donghua Li, Yu Liu, Wenna Wang, Tiantian Qiu Curdione Induces Antiproliferation Effect on Human Uterine Leiomyosarcoma via Targeting IDO1 *Front Oncol.* 2019 Mar. doi: 10.3389/fonc.2021.637024.
196. Qingying Wang, Zuozeng Cao, Yingze Wei, Jiawen Zhang, Zhongping Cheng Potential Role of SWI/SNF Complex Subunit Actin-Like Protein 6A in Cervical Cancer *Front Oncol.* 2019 Mar 1. doi: 10.3389/fonc.2021.724832.
197. Lixia Zhou, Tenghua Yu, Fei Yang, Jingjing Han, Bin Zuo, Lulu Huang, Xia Bai, Miao Jiang, Depei Wu, Suning Chen, Lijun Xia, Jia Ruan, Changgeng Ruan G Protein-Coupled Estrogen Receptor Agonist G-1 Inhibits Mantle Cell Lymphoma Growth in Preclinical Models *Front Oncol.* 2019 Mar 13. doi: 10.3389/fonc.2021.668617.
198. Ying Li, Jun Ma, Ziteng Song, Yinan Zhao, Han Zhang, Yeling Li, Jing Xu, Yuanqiang Guo The Antitumor Activity and Mechanism of a Natural Diterpenoid From *Casearia graveolens* *Front Oncol.* 2019 May. doi: 10.3389/fonc.2021.688195.
199. Zhongwei Wu, Zhaohui Ding, Bo Cheng, Zongchao Cui The inhibitory effect of human DEFA5 in growth of gastric cancer by targeting BMI1 *Cancer Sci.* 2019 May 15. doi: 10.1111/cas.14827.
200. Lei Yu, Yuxi Chen, Shaohui Yuan, Yang Cao, Zhenggang Bi Peiminine Induces G0/G1-Phase Arrest, Apoptosis, and Autophagy via the ROS/JNK Signaling Pathway in Human Osteosarcoma Cells in Vitro and in Vivo *Front Pharmacol.* 2019 May 7. doi: 10.3389/fphar.2021.770846.
201. Cuixiang Xu, Xiaoyan Huang, Xiaohua Lei, Zhankui Jin, Min Wu, Xiao Liu, Yubin Huang, Xiangrong Zhao, Yue Xiong, Jingying Sun, Xianglong Duan, Jianhua Wang Costunolide-Induced Apoptosis via Promoting the Reactive Oxygen Species and Inhibiting AKT/GSK3 β Pathway and Activating Autophagy in Gastric Cancer *Front Cell Dev Biol.* 2019 Nov. doi: 10.3389/fcell.2021.722734.
202. Simeng Zhang, Zhongyan Hua, Gen Ba, Ning Xu, Jianing Miao, Guifeng Zhao, Wei Gong, Zhihui Liu, Carol J Thiele, Zhijie Li Antitumor effects of the small molecule DMAMCL in neuroblastoma via suppressing aerobic glycolysis and targeting PFKFB3 *Cancer Cell Int.* 2019 Nov. doi: 10.1186/s12935-021-02330-y.
203. Xuan Zhang, Qing Li, Aibe Du, Yifei Li, Qing Shi, Yanrong Chen, Yang Zhao, Bin Wang, Feng Pan Adipocytic Glutamine Synthetase Upregulation via Altered Histone Methylation Promotes 5FU Chemoresistance in Peritoneal Carcinomatosis of Colorectal Cancer *Front Oncol.* 2019 Nov. doi: 10.3389/fonc.2021.748730.
204. Li Cheng, Dongmei Zhang, Wei Yan Ultrasound-targeted microbubble destruction-mediated overexpression of Sirtuin 3 inhibits the progression of ovarian cancer *Oncol Rep.* 2019 Nov 17. doi: 10.3892/or.2021.8171.
205. Jiaqi Chen, Qiaoya Pan, Yang Bai, Xuepeng Chen, Yi Zhou Hydroxychloroquine Induces Apoptosis in Cholangiocarcinoma via Reactive Oxygen Species Accumulation Induced by Autophagy Inhibition *Front Mol Biosci.* 2019 Sep. doi: 10.3389/fmolb.2021.720370.
206. Song Yu, Bingxuan Jia, Na Liu, Dianzhen Yu, Shuo Zhang, Aibo Wu Fumonisin B1 triggers carcinogenesis via HDAC/PI3K/Akt signalling pathway in human esophageal epithelial cells *Sci Total Environ.* 2019 Sep. doi: 10.1016/j.scitotenv.2021.147405.
207. Limei Yan, Zeping He, Wei Li, Ning Liu, Song Gao P76RBE silencing inhibits ovarian cancer cell proliferation, migration, and invasion via suppressing the integrin β 1/NF- κ B pathway *Cell Cycle.* 2019 Sep 1. doi: 10.1080/15384101.2021.1963910.
208. Liguang Xu, Xiuxiu Wang, Weiwei Wang, Maozhong Sun, Won Jin Choi, Ji-Young Kim, Changlong Hao, Si Li, Aihua Qu, Meiru Lu, Xiaoling Wu, Felipe M Colombari, Weverson R Gomes, Asdrubal L Blanco, Andre F de Moura, Xiao Guo, Hua Kuang, Nicholas A Kotov, Chuanlai Xu Enantiomer-dependent immunological response to chiral nanoparticles *Nature.* 2019 Sep 15. doi: 10.1038/s41586-021-04243-2.
209. Guihong Lu, Xiaojun Wang, Feng Li, Shuang Wang, Jiawei Zhao, Jinyi Wang, Jing Liu, Chengliang Lyu, Peng Ye, Hui Tan, Weiping Li, Guanghui Ma, Wei Wei Engineered biomimetic nanoparticles achieve targeted delivery and efficient metabolism-based synergistic therapy against glioblastoma *Nat Commun.* 2019 Sep 18. doi: 10.1038/s41467-022-31799-y.
210. Liang Yong, Yafen Yu, Bao Li, Huiyao Ge, Qi Zhen, Yiwen Mao, Yanxia Yu, Lu Cao, Ruixue Zhang, Zhuo Li, Yirui Wang, Wencheng Fan, Chang Zhang, Daiyue Wang, Sihan Luo, Yuanming Bai, Shirui Chen, Weiwei Chen, Miao Liu, Jijia Shen, Liangdan Sun Calcium/calmodulin-dependent protein kinase IV promotes imiquimod-induced psoriatic inflammation via macrophages and keratinocytes in mice *Nat Commun.* 2019 Sep 7. doi: 10.1038/s41467-022-31935-8.
211. Zhiqing Bai, Zhiying Guo, Jiaxing Liu, Yu-Ann Chen, Qian Lu, Ping Zhang, Lili Hong, Yunfang Wang, Jiahong Dong Lapatinib Suppresses HER2-Overexpressed Cholangiocarcinoma and Overcomes ABCB1-Mediated Gemcitabine Chemoresistance *Front Oncol.* 2020 Jan. doi: 10.3389/fonc.2022.860339.
212. Jun Bai, Jun Li, Ning Liu, Hai Jia, Xuemeng Si, Yusong Zhou, Zhian Zhai, Ying Yang, Fazheng Ren, Zhenlong Wu Glucosamine alleviates zearalenone-induced damage to porcine trophectoderm cells by activating the PI3K/AKT signaling pathway *Food Funct.* 2020 Mar 1. doi: 10.1039/d2fo00928e.
213. Huifeng Jiang, Erjiang Tang, Ying Chen, Hailong Liu, Yun Zhao, Moubin Lin, Luwei He Squalene synthase predicts poor prognosis in stage I-III colon adenocarcinoma and synergizes squalene epoxidase to promote tumor progression *Cancer Sci.* 2020 May. doi: 10.1111/cas.15248.

注：更多使用本产品的文献请参考产品网页。