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Methods	Separation principle	Advantages	Disadvantages
Ultracentrifugation (超速离心法)	Density; Floatation speed (in case of density gradient)	Most widely used as the Gold Standard; Robust method; Possibility to identify subpopulations when density gradient is applied.	Heavy workload; Might lead to ruptured EVs; Low EV yield; Results are difficult to compare between different studies due to different rotor types and centrifugation speeds; Large volume of starting material

			required.
Precipitation (沉淀法)	Solubility	User-friendly method; High EV yield; Many kits commercially available; Relatively cheap method; Possibility for clinical application.	High risk of co-isolating non-vesicular contaminants; Polymers may be problematic for some downstream applications; Might not be suitable for performing functional studies.
Size-exclusion chromatography (尺寸排阻色谱法)	Hydrodynamic radius	Fast and simple; EV integrity is maintained during isolation; May be used in upscaling processes.	Potentially significant contamination with non-EV particles (e.g. lipoproteins, albumin).
Ultrafiltration (超滤法)	Hydrodynamic radius	Can be very well combined with other methods; Easy to perform in every lab.	Filter type and ultrafiltration method influence EV yield and composition considerably; Rather limited resolution.
Field-flow fractionation (场流分离法)	Hydrodynamic radius; Electrophoretic mobility	Can identify smaller EV subsets (e.g. exomeres).	Specific equipment required and the associated expertise; Cannot be used as a stand-alone method.
Immunoaffinity capture (免疫亲和分离法)	Surface markers	Very pure EV isolations; Can be combined with light scattering flow cytometry; Capable of isolating subpopulations.	Currently not suitable for a complete isolation of all EVs; Expensive.
Microfluidics (微流控技术)	Surface markers; Hydrodynamic radius; Density	EV isolation can be combined with EV characterization; Subpopulations can be isolated; Small sample volumes can be used for EV isolation.	Requires a specific level of expertise; Not suitable for preparative purposes (e.g. therapeutic applications).
Membrane affinity methods (膜亲和法)	Membrane affinity	Is able to isolate all EVs instead of only EV subpopulations.	Cannot discriminate between EVs and cell(s) fragments.

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