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
Base pairs	6097bp
Feature Nucleotide	Position
PGK promoter	20-538
luc2 reporter gene	549-2202
Multiple cloning region	2208-2261
SV40 late poly (A) sinal	2268-2512
SV40 early enhancer/promoter	2554-2972
Synthetic neomycin phosphotransferase(Neor) coding region	2998-3792
Synthetic poly(A) signal	3817-3865
Synthetic Beta-lactamase(Ampr) coding region	4980-5840
Synthetic poly(A) signal/transcriptional pause site	5945-6097

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- Eosinophils Recruitment in OVA-Induced Murine Asthma by Targeting Inducible Co-Stimulator. *Int Arch Allergy Immunol.* 2019;180(3):182-194.
3. Yang K, Shi J, Hu Z, Hu X. The deficiency of miR-214-3p exacerbates cardiac fibrosis via miR-214-3p/NLRP5 axis. *CLIN SCI.* 2019 Sep 10;133(17):1845-1856
 4. Liu Y, Chen G, Liu H, Li Z, Yang Q, Gu X, Du Z, Zhang G, Wang J. Integrated bioinformatics analysis of miRNA expression in Ewing sarcoma and potential regulatory effects of miR-21 via targeting ALCAM/CD166. *ARTIF CELL NANOMED B.* 2019 Dec;47(1):2114-2122
 5. Jiekun Guo, Qiang Wu, Xiaoming Peng, Bin Yu. miR-509-5p Inhibits the Proliferation and Invasion of Osteosarcoma by Targeting TRIB2. *Biomed Res Int.* 2019 Dec 16;2019:2523032.
 6. Zhaochuan Yang, Zhenghai Qu, Mingji Yi, Zhidong Lv, Yanxia Wang, Yanchun Shan, Ni Ran, Xinjie Liu. MiR-204-5p Inhibits Transforming Growth Factor- β 1-Induced Proliferation and Extracellular Matrix Production of Airway Smooth Muscle Cells by Regulating Six1 in Asthma. *Int Arch Allergy Immunol.* 2020;181(4):239-248.
 7. Lu K, Feng F, Yang Y, Liu K, Duan J, Liu H, Yang J, Wu M, Liu C, Chang Y. High-throughput screening identified miR-7-2-3p and miR-29c-3p as metastasis suppressors in gallbladder carcinoma. *J Gastroenterol.* 2020 Jan;55(1):51-66
 8. Kai Lu, Feiling Feng, Yingcheng Yang, Kai Liu, Jicheng Duan, Hu Liu, Jiahe Yang, Mengchao Wu, Chen Liu, Yanxin Chang. High-throughput screening identified miR-7-2-3p and miR-29c-3p as metastasis suppressors in gallbladder carcinoma. *J Gastroenterol.* 2020 Jan;55(1):51-66.
 9. Tianyu Feng, Zhongkai Zhu, Yaqian Jin, Hao Wang, Xiaohan Mao, Dan Liu, Yiling Li, Lixia Lu, Guowei Zuo. The microRNA-708-5p/ZEB1/EMT axis mediates the metastatic potential of osteosarcoma. *Oncol Rep.* 2020 Feb;43(2):491-502.
 10. Gaoyang Chen, Qingyu Wang, Zhaoyan Li, Qiwei Yang, Yuzhe Liu, Zhenwu Du, Guizhen Zhang, Yang Song. Circular RNA CDR1as promotes adipogenic and suppresses osteogenic differentiation of BMSCs in steroid-induced osteonecrosis of the femoral head. *Bone.* 2020 Apr;133:115258.
 11. Xuehan Fan, Jingya Wang, Tingting Qin, Yujia Zhang, Wenting Liu, Kaiting Jiang, Dingzhi Huang. Exosome miR-27a-3p secreted from adipocytes targets ICOS to promote antitumor immunity in lung adenocarcinoma. *Thorac Cancer.* 2020 Jun;11(6):1453-1464.
 12. Junmei Niu, Bianfang Zhang, Kuili Cui, Yuan Gao, Zhenkui Li, Zhibin Qian. Suppression of miR-147b contributed to H37Rv-infected macrophage viability and migration in tuberculosis in vitro. *MICROB PATHOGENESIS.* 2020 Jul;144:104125.
 13. Weitie Wang, Yong Wang, Hulin Piao, Bo Li, Zhicheng Zhu, Dan Li, Tiance Wang, Kexiang Liu. Bioinformatics Analysis Reveals MicroRNA-193a-3p Regulates ACTG2 to Control Phenotype Switch in Human Vascular Smooth Muscle Cells. *Front Genet.* 2021 Jan 12;11:572707
 14. Fei Xue, Jian Wu, Wei Feng, Ting Hao, Yuan Liu, Wenbo Wang. MicroRNA-141 inhibits the differentiation of bone marrow-derived mesenchymal stem cells in steroid-induced osteonecrosis via E2F3. *Mol Med Rep.* 2022 Jul;26(1):234.
 15. Junzhe Zhang, Kaini Yang, Junfeng Bu, Jiayan Yan, Xiaoqiang Hu, Ke Liu, Si Gao, Shuibin Tang, Lili Gao, Wei Chen. IGF2BP3 promotes progression of gallbladder carcinoma by stabilizing KLK5 mRNA in N6-methyladenosine-dependent binding. *Front Oncol.* 2022 Oct 13;12:1035871.

Version 2024.03.12