碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-168-3301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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1. Zheng S, Wang G, Lin S. Heat shock effects and population survival in the polar dinoflagellate *Polarella glacialis* Journal of Experimental Marine Biology and Ecology. 2012 Dec;138:100-8.
2. Zhang Y, Xu Y, Zhu Q, Zhao F, Luo J, Zhang X, Wang X. Upregulation of dysbindin in temporal lobe epileptic foci of human and experimental animals. Synapse. 2012 Jul;66(7):622-9.
3. Zhu W, Mao Z, Zhu C, Li M, Cao C, Guan Y, Yuan J, Xie G, Guan X. Adolescent exposure to cocaine increases anxiety-like behavior and induces morphologic and neurochemical changes in the hippocampus of adult rats. Neuroscience. 2016 Jan 28;313:174-83.
4. Tang Y, Lv P, Sun Z, Han L, Zhou W. 14-3-3 $\beta$  Promotes Migration and Invasion of Human Hepatocellular Carcinoma Cells by Modulating Expression of MMP2 and MMP9 through PI3K/Akt/NF- $\kappa$ B Pathway. PLoS One. 2016 Jan 5;11(1):e0146070.
5. Lv L, Liu HG, Dong SY, Yang F, Wang QX, Guo GL, Pan YF, Zhang XH. Upregulation of CD44v6 contributes to acquired chemoresistance via the modulation of autophagy in coloncancer SW480 cells. Tumour Biol. 2016 Jul;37(7):8811-24.

6. Shi B, Li Y, Wang X, Yang Y, Li D, Liu X, Yang X . Silencing of hypoxia inducible factor-1 $\alpha$  by RNA interference inhibits growth of SK-NEP-1 Wilms tumour cells in vitro, and suppresses tumourigenesis and angiogenesis in vivo. *Clin Exp Pharmacol Physiol*. 2016 Jun;43(6):626-33.
7. Liu C, Huang J, Li H, Yang Z, Zeng Y, Liu J, Hao Y, Li R . Ghrelin accelerates wound healing through GHS-R1a-mediated MAPK-NF- $\kappa$ B/GR signaling pathways in combined radiation and burn injury in rats. *Sci Rep*. 2016 Jun 7;6:27499.
8. Fang K, Dong H, Jiang S, Li F, Wang D, Yang D, Gong J, Huang W, Lu F . Diosgenin and 5-Methoxypsoralen Ameliorate Insulin Resistance through ER- $\alpha$ /PI3K/Akt-Signaling Pathways in HepG2 Cells. *Evid Based Complement Alternat Med*. 2016;2016:7493694.
9. Cai M, Wang H, Li JJ, Zhang YL, Xin L, Li F, Lou SJ . The signaling mechanisms of hippocampal endoplasmic reticulum stress affecting neuronal plasticity-related protein levels in high fat diet-induced obese rats and the regulation of aerobic exercise. *Brain Behav Immun*. 2016 Oct;57:347-59.
10. Xu Y, Liu J, He M, Liu R, Belegu V, Dai P, Liu W, Wang W, Xia QJ, Shang FF, Luo CZ, Zhou X, Liu S, McDonald J, Liu J, Zuo YX, Liu F, Wang TH . Mechanisms of PDGF siRNA-mediated inhibition of bone cancer pain in the spinal cord. *Sci Rep*. 2016 Jun 10;6:27512.
11. Pang M, Teng Y, Huang J, Yuan Y, Lin F, Xiong C . Substrate stiffness promotes latent TGF- $\beta$ 1 activation in hepatocellular carcinoma. *Biochem Biophys Res Commun*. 2017 Jan 29;483(1):553-558.
12. Zhao Y, Yu Y, Zhang Y, He L, Qiu L, Zhao J, Liu M, Zhang J . Letrozole regulates actin cytoskeleton polymerization dynamics in a SRC-1 dependent manner in the hippocampus of mice. *J Steroid Biochem Mol Biol*. 2017 Mar;167:86-97.

Version 2017.03.05