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
- Cai H, Yin D, Zhang L, Yang X, Xu X, Liu W, Zheng X, Zhang H, Wang J, Xu Y, Cheng D, Zheng M, Han Y, Wu M, Wang Y. Preparation and biological evaluation of 2-amino-6-[18F]fluoro-9-(4-hydroxy-3-hydroxy-methylbutyl) purine (6-[18F]FPCV) as a novel PET probe for imaging HSV1-tk reporter gene expression. *Nucl Med Biol.* 2007 Aug; 34(6):717-25.
- He M, Tang WX, Jiang R, Zhang N, Cha L, Fan J, Wang ZW. Differential Expression of DNMTs Between Gastric Tissues with and without H. Pylori Infection. *Chin J Biochem Mol Biol.* 2011;27(11):1032-8.
- Xie P, Pan L, Xu W, Yue F. The roles of serine protease, intracellular and extracellular phenoloxidase in activation of prophenoloxidase system, and characterization of phenoloxidase from shrimp haemocytes induced by lipopolysaccharide or dopamine. *Chinese Journal of Oceanology and Limnology.* 2013 Sep; 31(5):1018-1027.
- Chen X, Wang H, Ou-Yang XN, Xie FW, Wu JJ. Research on drug resistance mechanism of trastuzumab caused by activation of the PI3K/Akt signaling pathway. *Contemp Oncol (Pozn).* 2013;17(4):363-9.
- Wang HJ, Tan YZ. Methods for assessing effects of Wnt/ β -catenin signaling in senescence of mesenchymal stem cells. *Methods Mol Biol.* 2013;976:111-30.
- Liao W, Sun L, Wang C, Huang H, Liu J, Liao W, Shi M. LC3A-positive "stone-like" structures predict an adverse prognosis of gastric cancer. *Anat Rec (Hoboken).* 2014 Apr;297(4):653-62.
- Cai C, Zhao D, Ma C, Zhang Y, Wu X, Wei G, He D. Connexin 43 expression in Sprague-Dawley rat seminiferous epithelium after in utero exposure to flutamide. *Syst Biol Reprod Med.* 2014 Oct;60(5):257-62.
- Liu G, Song Y, Cui L, Wen Z, Lu X. Inositol hexaphosphate suppresses growth and induces apoptosis in HT-29 colorectal cancer cells in culture: PI3K/Akt pathway as a potential target. *Int J Clin Exp Pathol.* 2015 Feb 1;8(2):1402-10.
- Zhu MY, Lin D, Liu J, Mo HY. Artesunate interfere in modulation of Foxp3 expression in synovial cells in collagen-induced arthritis rats. *Chin J Integr Med.* 2016 Jun 28:1-7.
- Gong T, Wei Q, Mao D, Shi F. Expression patterns of taste receptor type 1 subunit 3 and α -gustducin in the mouse testis during development. *Acta Histochem.* 2016 Jan;118(1):20-30.
- Li L, Li YM, Liu ZL, Zhang JG, Liu Q, Yi LT. The renal protective effects of Anoectochilus roxburghii polysaccharose on diabetic mice induced by high-fat diet and streptozotocin. *J Ethnopharmacol.* 2016 Feb 3;178:58-65.
- Sui Y, Fu X, Wang Y, Hu W, Zhang T, Liu W, Jiang L, Xing S, Fu X, Xu X. Expression, purification and characterization of a catalytic domain of human protein tyrosine phosphatase non-receptor 12 (PTPN12) in Escherichia coli with FKBP-type PPIase as a chaperon. *Protein Expr Purif.* 2018 Feb;142:45-52.
- Cao W, Zhang J, Wang G, Lu J, Wang T, Chen X. Reducing-Autophagy Derived Mitochondrial Dysfunction during Resveratrol Promotes Fibroblast-Like Synovial Cell Apoptosis. *Anat Rec (Hoboken).* 2018 Feb 20.
- Wu G, Song D, Wei Q, Xing J, Shi X, Shi F. Melatonin mitigates bisphenol A-induced estradiol production and proliferation by porcine ovarian granulosa cells in vitro. *Anim Reprod Sci.* 2018 May;192:91-98.
- Chen C, Chen P, Fang X. The FLC dimer with lambda type may false-migrate to the position of "albumin" band by urine protein electrophoresis on Sebia agarose gel-based detection system. *J Clin Lab Anal.* 2019 Feb;33(2):e22658.
- Geng F, Xie Y, Wang J, Li S, Jin Y, Ma M. Large-scale purification of ovalbumin using polyethylene glycol precipitation and isoelectric precipitation. *J POULT SCI.* 2019 Mar 1;98(3):1545-1550.
- Bai Y, Wang J, He Z, Yang M, Li L, Jiang H. Mesenchymal Stem Cells Reverse Diabetic Nephropathy Disease via Lipoxin A4 by Targeting Transforming Growth Factor β (TGF- β)/smad Pathway and Pro-Inflammatory Cytokines. *MED SCI MONITOR.* 2019 Apr 26;25:3069-3076.
- Zhang W, Zheng J, Tian X, Tang Y, Ding G, Yang Z, Jin H. Pepsin-Soluble Collagen from the Skin of Lophius litulo: A Preliminary Study Evaluating Physicochemical, Antioxidant, and Wound Healing Properties. *Mar Drugs.* 2019 Dec 16;17(12). pii: E708.
- Shuang Wu, Tiantian Yang, Kai Cen, Yihuai Zou, Xiaowei Shi, Dongrui Zhou, Yonghong Gao, Limin Chai, Yizhou Zhao, Yikun Sun, Lingqun Zhu. In Vitro Evaluation of the Neuroprotective Effect of Panax notoginseng Saponins by Activating the EGFR/PI3K/AKT Pathway. *Evid Based Complement Alternat Med.* 2020 Jul 31;2020:1403572.
- Xinran Dong, Lin Yang, Kaiyi Liu, Xiaoli Ji, Chuanqing Tang, Wanxing Li, Ling Ma, Yuting Mei, Ting Peng, Ban Feng, Ziyang Wu, Qingyuan Tang, Yanyan Gao, Kai Yan, Wenhao Zhou, Man Xiong. Transcriptional networks identify synaptotagmin-like 3 as a regulator of cortical neuronal migration during early neurodevelopment. *Cell Rep.* 2021 Mar 2;34(9):108802.
- Fengzhen Liu, Xia Zhang, You Peng, Liping Zhang, Yao Yu, Peng Hua, Peiying Zhu, Xinyu Yan, Yin Li, Liang Zhang. miR-24 controls the regenerative competence of hair follicle progenitors by targeting Plk3. *Cell Rep.* 2021 Jun 8;35(10):109225.

Version 2023.11.01