碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Wang F, Liang YJ, Wu XP, Su XD, Fu LW. . ABCG2-overexpressing S1-M1-80 cell xenografts in nude mice keep original biochemistry and cell biological properties. *Chin J Cancer*. 2012 Mar;31(3):150-8.
2. Liu X, Zhang Z, Li P, Zhu L, Wang Y, Wang C. . Polypeptide from *Chlamys farreri* modulates UVB-induced activation of NF-kappaB signaling pathway and protection HaCaT cells from apoptosis. *REGUL PEPTIDES*. 2009 Feb 25;153(1-3):49-55.
3. Liu HT, Li WM, Li XY, Xu QS, Liu QS, Bai XF, Yu C, Du YG. . Chitosan oligosaccharides inhibit the expression of interleukin-6 in lipopolysaccharide-induced human umbilical vein endothelial cells through p38 and ERK1/2 protein kinases. *BASIC CLIN PHARMACOL*. 2010 May;106(5):362-71.
4. Wang XH, Fu ZG, Zhao Y, Shen W, Wu XY, Wang CY. . Survivin: an overexpression protein with notable cellular localization and multiple roles in colon cancer. 《中德临床肿瘤学杂志》. 2010 Sep;9(9):519-23.
5. Wang X, Fu Z, Zhao Y, Wu X, Shen W. . Profile of protein expression of the colon cancer cell line SW480 with survivin/shRNA. *Eur J Cancer Prev*. 2011 May;20(3):190-8.
6. Wang X, Zhao Y, Fu Z, He Y, Xiang D, Zhang L. . Preliminary autogenic endothelial cells in allogeneic vessels inhibits thrombosis and intimal hyperplasia: an efficacy study in dogs. *J Surg Res*. 2011 Jul;169(1):148-55.
7. Zhou L, Xue H, Wang Z, Ni J, Yao T, Huang Y, Yu C, Lu L. . Angiotensin-(1-7) attenuates high glucose-induced proximal tubular epithelial-to-mesenchymal transition via inhibiting ERK1/2 and p38 phosphorylation. *Life Sci*. 2012 Mar 10;90(11-12):454-62.
8. Nie J, Cao X, Zhou M, Zhang X, Zhang R, Niu L, Xia Y, Hou L, Wang C. . *Musca domestica* pupae lectin induces apoptosis in HepG2 cells through a NF-kB/p65-mediated caspase pathway. *B CANCER*. 2012 Apr 1;99(4):E49-54.
9. Meng X, Li F, Chen S, Tang C, Zhang W, Wang Z, Zhao S. . Cloning and expression of neuron-specific enolase in the corpus luteum of dairy goats. *Gene*. 2012 Jul 25;503(2):222-8.
10. Yuan Z, Li Y, Hu Y, You J, Higashisaka K, Nagano K, Tsutsumi Y, Gao J. . Chitosan nanoparticles and their Tween 80 modified counterparts disrupt the developmental profile of zebrafish embryos. *INT J PHARMACOL*. 2016 Dec 30;515(1-2):644-656.
11. Dong XL, Liu TH, Wang W, Pan CX, Du GY, Wu YF, Pan MH, Lu C. . BmNHR96 participate BV entry of BmN-SWU1 cells via affecting the cellular cholesterol level. *BIOCHEM BIOPH RES CO*. 2017 Jan 22;482(4):1484-1490.
12. Liao S, Gan L, Qin W, Liu C, Mei Z. . Inhibition of GSK3 and MEK induced cancer stem cell generation via the Wnt and MEK signaling pathways. *Oncol Rep*. 2018 Oct;40(4):2005-2013.
13. Shen J, Tu L, Chen D, Tan T, Wang Y, Wang S. . TRPV4 channels stimulate Ca²⁺-induced Ca²⁺ release in mouse neurons and trigger endoplasmic reticulum stress after intracerebral hemorrhage. *Brain Res Bull*. 2019 Mar 146:143-152.
14. Yang K, Yao Y. . Mechanism of GPER promoting proliferation, migration and invasion of triple-negative breast cancer cells through CAF. *Am J Transl Res*. 2019 Sep 15 11(9):5858-5868. eCollection 2019.
15. Xie W, Zhu T, Dong X, Nan F, Meng X, Zhou P, Sun G, Sun X. . HMGB1-triggered inflammation inhibition of notoginseng leaf triterpenes against cerebral ischemia and reperfusion injury via MAPK and NF-kB signaling pathways. *Biomolecules*. 2019 Sep 20 9(10). pii: E512.
16. Jie Cheng, Meng Zhang, Shaoli Cheng, Fan Li, Bingyi Zhang, Xiaoming Sun, Huijuan Hu, Lina Chen, Zhenghang Zhao, Hao Hu, Zhanqin Zhang. . Low-dose alcohol ameliorates high fat diet-induced anxiety-related behavior via enhancing adiponectin expression and activating the Nrf2 pathway. *Food Funct*. 2021 Jan 7;12(1):241-251.
17. Huiling Zhang, Xuedi Zhang, Chunxiu Ling, Change Liu, Shaodong Hua, Ziyang Xiong, Huimin Zhao, Qiushuang Feng, Zhichun Feng, Jing Tang, Xiaoyang Hong. . EGFR-TNFR1 pathway in endothelial cell facilitates acute lung injury by NF-kB/MAPK-mediated inflammation and RIP3-dependent necroptosis. *Int Immunopharmacol*. 2023 Apr;117:109902.
18. Jian Cheng, Guizhong Zhang, Tian Deng, Zhao Liu, Mengqi Zhang, Pengchao Zhang, Funmilayo O Adeshakin, Xiangyun Niu, Dehong Yan, Xiaochun Wan, Guang Yu. . CD317 maintains proteostasis and cell survival in response to proteasome inhibitors by targeting calnexin for RACK1-mediated autophagic degradation. *Cell Death Dis*. 2023 May 20;14(5):333.
19. Lijuan Chen, Jing Li, Xinglian Liu, Zhiwei Zhao, Yan Jin, Yikun Fu, Aiqin Zhou, Chengqun Wang, Yan Zhou. . Vitamin B6 Deficiency Induces Autism-Like Behaviors in Rats by Regulating mTOR-Mediated Autophagy in the Hippocampus. *Behav Neurol*. 2023 May 9;2023:6991826.

Version 2025.03.12