碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Md Sodrul Islam, Lingyan Miao, Hui Yu, Ziyi Han, Hongxiang Sun. Ethanol Extract of Illicium henryi Attenuates LPS-Induced Acute Kidney Injury in Mice via Regulating Inflammation and Oxidative Stress. *Nutrients*. 2019 Jun 23;11(6):1412.
2. Ting Liu, Jing-Zhuo Huang, Ze-Yuan Lei, Rong-Shuai Yan, Dong-Li Fan. FoxO3a depletion accelerates cutaneous wound healing by regulating epithelial-mesenchymal transition through β -catenin activation. *Mol Med Rep*. 2020 Mar;21(3):1224-1232.
3. Jian-Hua Zhou, Zhi-Xian Yao, Zhong Zheng, Jun Yang, Rui Wang, Shi-Jie Fu, Xu-Feng Pan, Zhi-Hong Liu, Ke Wu. G-MDSCs-Derived Exosomal miRNA-143-3p Promotes Proliferation via Targeting of ITM2B in Lung Cancer. *Onco Targets Ther*. 2020 Sep 30;13:9701-9719.
4. Bin Zhan, Linjin Huang, Yachun Chen, Wen Ye, Jingkun Li, Jianhui Chen, Sheng Yang, Wei Jiang. miR-196a-mediated downregulation of p27 kip1 protein promotes prostate cancer proliferation and relates to biochemical recurrence after radical prostatectomy. *Prostate*. 2020 Sep;80(12):1024-1037.
5. Lin Li, Haiyan Gan, Huaqian Jin, Yan Fang, Yan Yang, Jianping Zhang, Xiaowei Hu, Lisheng Chu. Astragaloside IV promotes microglia/macrophages M2 polarization and enhances neurogenesis and angiogenesis through PPAR γ pathway after cerebral ischemia/reperfusion injury in rats. *Int*

Immunopharmacol. 2021 Mar;92:107335.

6. Shuai Ren, Yiqing Cai, Shunfeng Hu, Jiarui Liu, Yi Zhao, Mengfei Ding, Xiaomin Chen, Linqun Zhan, Xiangxiang Zhou, Xin Wang. Berberine exerts anti-tumor activity in diffuse large B-cell lymphoma by modulating c-myc/CD47 axis. *Biochem Pharmacol.* 2021 Jun;188:114576.
7. Siliang Wang, Cheng Tang, Heng Zhao, Peiliang Shen, Chao Lin, Yun Zhu, Dan Han. Network Pharmacological Analysis and Experimental Validation of the Mechanisms of Action of Si-Ni-San Against Liver Fibrosis. *Front Pharmacol.* 2021 Jul 1;12:656115.
8. Zhihui Guan, Lingmin Zhou, Yu Zhang, Yu Zhang, Hongping Chen, Feifei Shao. Sulforaphane Ameliorates the Liver Injury of Traumatic Hemorrhagic Shock Rats. *J Surg Res.* 2021 Nov;267:293-301.
9. Shijing Sun, Qingyu Kong, Zhifeng Cai, Minmin Wang, Haizhao Zhao, Cuifen Zhao. circ-Grm1 promotes pulmonary artery smooth muscle cell proliferation and migration via suppression of GRM1 expression by FUS. *Int J Mol Med.* 2021 Nov;48(5):202.
10. Peng Li, Li-Yan Yao, Yu-Jia Jiang, Dan-Dan Wang, Ting Wang, Yan-Ping Wu, Bai-Xiang Li, Xue-Ting Li. Soybean isoflavones protect SH-SY5Y neurons from atrazine-induced toxicity by activating mitophagy through stimulation of the BEX2/BNIP3/NIX pathway. *Ecotoxicol Environ Saf.* 2021 Dec 20;227:112886.
11. Shuang Liu, Wen Yan, Yanbing Hu, Huiying Wu. Shikonin Alleviates Endothelial Cell Injury Induced by ox-LDL via AMPK/Nrf2/HO-1 Signaling Pathway. *Evid Based Complement Alternat Med.* 2021 Dec 6;2021:5881321.
12. Yu-Qin Di, Xiao-Lin Han, Xin-Le Kang, Di Wang, Cai-Hua Chen, Jin-Xing Wang, Xiao-Fan Zhao. Autophagy triggers CTSD (cathepsin D) maturation and localization inside cells to promote apoptosis. *Autophagy.* 2021 May;17(5):1170-1192.

Version 2023.11.13