Heparin Agarose (肝素琼脂糖凝胶)

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Characteristics	Description
Product content	50% settled gel in specific protective buffer
Agarose structure	6% cross-linked agarose
Average particle size	45-165 μ m
Coupled ligand	Heparin
Ligand concentration	~4mg Heparin per ml settled gel
Binding capacity	Per ml settled gel: \geq 2mg antithrombin III (ATIII)
Specificity	Electrostatic effect
Application	Separation of biomolecules with an affinity for heparin, including antithrombin III, coagulation factors and other plasma proteins, DNA binding proteins, lipoproteins, protein synthesis factors, enzymes that act on nucleic acids, and steroid receptors, viruses et al.

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产品编号	产品名称	包装
P2015-2ml	Protein A Agarose (Fast Flow, 抗体纯化用)	2ml
P2015-10ml	Protein A Agarose (Fast Flow, 抗体纯化用)	10ml
P2015-50ml	Protein A Agarose (Fast Flow, 抗体纯化用)	50ml
P2015-200ml	Protein A Agarose (Fast Flow, 抗体纯化用)	200ml
P2017-2ml	Protein G Agarose (Fast Flow, 抗体纯化用)	2ml
P2017-10ml	Protein G Agarose (Fast Flow, 抗体纯化用)	10ml
P2017-50ml	Protein G Agarose (Fast Flow, 抗体纯化用)	50ml
P2017-200ml	Protein G Agarose (Fast Flow, 抗体纯化用)	200ml
P2019-2ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	2ml

P2019-10ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	10ml
P2019-50ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	50ml
P2019-200ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	200ml
P2157-1ml	Biotin Agarose (生物素琼脂糖凝胶)	1ml
P2157-5ml	Biotin Agarose (生物素琼脂糖凝胶)	5ml
P2157-20ml	Biotin Agarose (生物素琼脂糖凝胶)	20ml
P2159-1ml	Streptavidin Agarose (链霉亲和素琼脂糖凝胶)	1ml
P2159-5ml	Streptavidin Agarose (链霉亲和素琼脂糖凝胶)	5ml
P2159-20ml	Streptavidin Agarose (链霉亲和素琼脂糖凝胶)	20ml
P2165-1ml	Heparin Agarose (肝素琼脂糖凝胶)	1ml
P2165-5ml	Heparin Agarose (肝素琼脂糖凝胶)	5ml
P2165-20ml	Heparin Agarose (肝素琼脂糖凝胶)	20ml

Version 2024.05.05