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Characteristics	Description
Product content	50% settled gel in 50% glycerol with 10 mM phosphate-buffered saline and preservative (pH7.4)
Matrix	4% agarose
Average bead size	90µm
Antibody	Mouse monoclonal antibody against HA-tag
Isotype	IgG1
M.W. of antibody	Approximately 150kDa
Antibody concentration	Approximately 7.5mg HA antibody per ml settled gel
Binding capacity	Approximately 1.5mg HA-tagged protein per ml settled gel
Elution method	Acid, alkaline, neutral, peptide competitive or SDS-PAGE loading buffer elution. Note: If elute with SDS-PAGE loading buffer, the light (~25kDa) and heavy (~50kDa) chain of antibody will be denatured and release from the gel.
Reagent compatibility	Chaotropic reagents will denature the target HA-tagged protein. Do not exceed 0.3M GuHCl or 1.5M Urea.
Application	Suitable for IP、Co-IP and protein purification.
Storage	-20°C

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Steps	Solution required	Volume per assay
Cell lysis and sample preparation	Lysis Buffer with Protease Inhibitor Cocktail	100µl
Preparation of HA Peptide elution buffer and agarose gel	TBS	~1.6ml
Immunoprecipitation	Agarose gel	20µl
Wash (3 times)	Lysis Buffer with Protease Inhibitor Cocktail	500µl each time
Peptide competitive elution (optional)	Peptide solution (1X)	100µl

Acid elution and neutralization (optional)	Acid Elution Buffer	100µl
	Neutralization Buffer	10µl
SDS-PAGE sample loading buffer elution (optional)	SDS-PAGE Sample Loading Buffer (2X)	20µl

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Problem	Possible Causes	Solution
Very few or no tagged protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the tagged protein by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
	Detection system is inadequate.	If Western blot detection is used:

		<ol style="list-style-type: none"> 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the reactivity was adequate by using prestained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Background is too high.	Proteins bind nonspecifically to the monoclonal antibody, insufficient washing on agarose gel, or the microcentrifuge tubes.	<ol style="list-style-type: none"> 1. Pre-clear lysate with Mouse IgG Agarose (P2265) to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before centrifugation.
	Washes are insufficient.	<ol style="list-style-type: none"> 1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Choose other wash buffers. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non-specific binding.	<ol style="list-style-type: none"> 1. Prepare cell lysate again. 2. Add additional wash steps.

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1. Lee C. Methods Mol Biol. 2007. 362:401-406.
2. Gevaert K, Vandekerckhove J. Electrophoresis. 2000. 21:1145-1154.
3. Toshiyuki U, David J F W, Pier-Luc T, Kelly P N, Joy E W, et al. ACS Synth Biol. 2019. 8(8):1809-1817.
4. Bargonetti J, Reynisdóttir I, Friedman PN, Prives C. Genes Dev. 1992. 6(10):1886-98.

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