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Characteristics	Description
Product Content	25% slurry in specific protective buffer
Beads Size	30~150μm
Magnetization	Superparamagnetic
Coupled Ligand	Iminodiacetic acid (IDA)
Chelating Metalion	Co ²⁺
Metal ion Density	30-50μmol/ml beads
Binding Capacity	20-30mg His-tagged fusion protein per ml beads (100%)
Specificity	Met-His-tag-Protein, His-tag-Protein, Protein-His-tag
Application	Protein purification, IP, Co-IP

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Buffer	Components
Binding/Wash Buffer	10mM Tris, 500mM NaCl, pH7.4
Elution Buffer	10mM Tris, 500mM NaCl, 500mM Imidazole, pH7.4

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Problem	Possible Causes	Solution
Very few or no His-tagged protein exists in the eluate.	Protein is not completely eluted.	1. Increase the Imidazole concentration of elution buffer. 2. Increase the elution or incubation time.
	No target protein expressed.	Make sure the protein of interest contains the His-tag by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using pre-stained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
The purity of the target protein is low.	There are other polyhistidine-rich proteins in sample	Try a pH step gradient elution or an imidazole step gradient elution.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Add 5-20mM imidazole in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
The binding capacity of the magnetic beads has declined.	Protein or other substances have aggregated on the beads	Wash the beads with NaOH
	The beads have been reused too many times.	Use new beads.

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- Bromberg R, Cai K, Guo Y, Plymire D, Emde T, et al. Front Mol Biosci. 2022. 9:912072.
- Zeng K, Sun EJ, Liu ZW, Guo J, Yuan C, et al. RSC Adv. 2020. 10(19):11524-11534.

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产品编号	产品名称	包装
P2135	BeyoMag™ Anti-His Magnetic Beads (Anti-His磁珠)	0.5/2ml
P2210	BeyoGold™ His-tag Purification Resin (耐还原螯合型)	10ml
P2218	BeyoGold™ His-tag Purification Resin (耐还原螯合型)	100ml
P2220	BeyoGold™ His-tag Purification Resin (耐还原螯合型)	1000ml
P2221	BeyoGold™ His-tag Purification Resin (Fast Flow, 耐还原螯合型)	10/100/1000ml
P2222	BeyoGold™ His-tag Purification Resin (NTA-Co)	10/100/1000ml
P2226	His标签蛋白纯化试剂盒(耐还原螯合型)	10ml

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