碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- Wei Lv, Wei Jiang, Hongmei Luo, Qian Tong, Xiaoyu Niu, Xiao Liu, Yang Miao, Jingnan Wang, Yiwen Guo, Jianan Li, Xizhen Zhan, Yunqing Hou, Yaxin Peng, Jian Wang, Shuhong Zhao, Zaiyan Xu, Bo Zuo. Long noncoding RNA lncMREF promotes myogenic differentiation and muscle regeneration by interacting with the Smarca5/p300 complex. *Nucleic Acids Res.* 2022 Oct 14;50(18):10733-10755.
- Anming Ding, Xianfeng Tang, Dahai Yang, Meng Wang, Angyan Ren, Zongchang Xu, Ruibo Hu, Gongke Zhou, Malcolm O'Neill, Yingzhen Kong. ERF4 and MYB52 transcription factors play antagonistic roles in regulating homogalacturonan de-methylesterification in Arabidopsis seed coat mucilage. *Plant Cell.* 2021 Apr 17;33(2):381-403.
- Lin Zhao, Xiangying Deng, Yueshuo Li, Jianmin Hu, Longlong Xie, Feng Shi, Min Tang, Ann M Bode, Xin Zhang, Weihua Liao, Ya Cao. Conformational change of adenine nucleotide translocase-1 mediates cisplatin resistance induced by EBV-LMP1. *EMBO Mol Med.* 2021 Dec 7;13(12):e14072.
- Tao Li, Bin Fu, Xin Zhang, Yunjiang Zhou, Mengdi Yang, Mengran Cao, Yaxin Chen, Yingying Tan, Rong Hu. Overproduction of Gastrointestinal 5-HT Promotes Colitis-Associated Colorectal Cancer Progression via Enhancing NLRP3 Inflammasome Activation. *Cancer Immunol Res.* 2021 Sep;9(9):1008-1023.
- Qidong Li, Hanhui Yao, Yichun Wang, Yang Wu, Rick F Thorne, Youming Zhu, Mian Wu, Lianxin Liu. circPRKAA1 activates a Ku80/Ku70/SREBP-1 axis driving de novo fatty acid synthesis in cancer cells. *Cell Rep.* 2022 Nov 22;41(8):111707.
- Wuyi Wang, Jinlai Wei, Hairong Zhang, Xiangru Zheng, He Zhou, Yajun Luo, Jianguo Yang, Qican Deng, Siqi Huang, Zhongxue Fu. PRDX2 promotes the proliferation of colorectal cancer cells by increasing the ubiquitinated degradation of p53. *Cell Death Dis.* 2021 Jun 11;12(6):605.
- Yun-Long Wang, Wan-Wen Zhao, Shao-Mei Bai, Yan Ma, Xin-Ke Yin, Li-Li Feng, Guang-Dong Zeng, Fang Wang, Wei-Xing Feng, Jian Zheng, Ying-Nai Wang, Bing Zeng, Quentin Liu, Mien-Chie Hung, Xiang-Bo Wan. DNA damage-induced paraspeckle formation enhances DNA repair and tumor radioresistance by recruiting ribosomal protein P0. *Cell Death Dis.* 2022 Aug 16;13(8):709.
- Hairong Shi, Hongshun Wang, Cheng Zhang, Yajie Lu, Jun Yao, Zhibin Chen, Guangqian Xing, Qinjun Wei, Xin Cao. Mutations in OSBPL2 cause hearing loss associated with primary cilia defects via sonic hedgehog signaling. *JCI Insight.* 2022 Feb 22;7(4):e149626.
- Yanli Yang, Hui-Guang Li, Jie Wang, Hou-Ling Wang, Fang He, Yanyan Su, Ying Zhang, Cong-Hua Feng, Mengxue Niu, Zhonghai Li, Chao Liu, Weilun Yin, Xinli Xia. ABF3 enhances drought tolerance via

- promoting ABA-induced stomatal closure by directly regulating ADF5 in *Populus euphratica*. *J Exp Bot.* 2020 Dec 31;71(22):7270-7285.
10. Chen Xie, Liu Yang, Guixian Jia, Kang Yan, Shizhong Zhang, Guodong Yang, Changai Wu, Yingping Gai, Chengchao Zheng, Jinguang Huang. Maize HEAT UP-REGULATED GENE 1 plays vital roles in heat stress tolerance. *J Exp Bot.* 2022 Oct 18;73(18):6417-6433.
 11. Xin-Juan Fan, Yun-Long Wang, Wan-Wen Zhao, Shao-Mei Bai, Yan Ma, Xin-Ke Yin, Li-Li Feng, Wei-Xing Feng, Ying-Nai Wang, Quentin Liu, Mien-Chie Hung, Xiang-Bo Wan. NONO phase separation enhances DNA damage repair by accelerating nuclear EGFR-induced DNA-PK activation. *Am J Cancer Res.* 2021 Jun 15;11(6):2838-2852.
 12. Xiao N, Li H, Yu W, Gu C, Fang H, Peng Y, Mao H, Fang Y, Ni W, Yao M. SUMO-specific protease 2 (SEN2) suppresses keratinocyte migration by targeting NDR1 for de-SUMOylation. *FASEB J.* 2019 Jan;33(1):163-174.
 13. Zhou X, Liu Y, Huang J, Liu Q, Sun J, Cai X, Tang P, Liu W, Miao W. High temperatures affect the hypersensitive reaction, disease resistance and gene expression induced by a novel harpin HpaG-Xcm. *SCI REP-UK.* 2019 Jan 30;9(1):990.
 14. Niu J, Huang Y, Li Y, Wang Z, Tang J, Wang B, Lu Y, Cai J, Jian J. Characterization of a tandem-repeat galectin-9 from Nile tilapia (*Oreochromis niloticus*) involved in the immune response against bacterial infection. *FISH SHELLFISH IMMUN.* 2019 Sep;92:216-223.
 15. Ni W, Li Y, Cai L, Dong C, Fang H, Chen Y, Li H, Yao M, Xiao N. SUMOylation is required for PIPK1 γ -driven keratinocyte migration and growth. *FEBS J.* 2019 Dec;286(23):4709-4720.
 16. Xiong Y, Yi Y, Wang Y, Yang N, Rudd CE, Liu H. Ubc9 Interacts with and SUMOylates the TCR Adaptor SLP-76 for NFAT Transcription in T Cells. *J Immunol.* 2019 Dec 1;203(11):3023-3036.
 17. Wen Zhang, Shou-Song Tao, Ting Wang, Jie Zhang, Xian Liu, Ya-Ting Li, Hui Chen, Yi-Qun Zhan, Miao Yu, Chang-Hui Ge, Chang-Yan Li, Guang-Ming Ren, Xiao-Ming Yang, Rong-Hua Yin. ABRO1 stabilizes the deubiquitinase BRCC3 through inhibiting its degradation mediated by the E3 ubiquitin ligase WWP2. *FEBS Lett.* 2021 Jan;595(2):169-182.
 18. Ning Zhang, Lingran Zhang, Linjie Li, Junyou Geng, Lei Zhao, Yan Ren, Zhongdong Dong, Feng Chen. Global Profiling of 2-hydroxyisobutyrylome in Common Wheat. *Genomics Proteomics Bioinformatics.* 2021 Feb 10;S1672-0229(21)00007-3.
 19. Yueshuo Li, Feng Shi, Jianmin Hu, Longlong Xie, Lin Zhao, Min Tang, Xiangjian Luo, Mao Ye, Hui Zheng, Min Zhou, Na Liu, Ann M Bode, Jia Fan, Jian Zhou, Qiang Gao, Shuangjian Qiu, Weizhong Wu, Xin Zhang, Weihua Liao, Ya Cao. Stabilization of p18 by deubiquitylase CYLD is pivotal for cell cycle progression and viral replication. *NPJ Precis Oncol.* 2021 Mar 2;5(1):14.
 20. Songbiao Chen, Bichen Miao, Nannan Chen, Caiyi Chen, Ting Shao, Xuezi Zhang, Lingling Chang, Xiujuan Zhang, Qian Du, Yong Huang, Dewen Tong. SYNCRIP facilitates porcine parvovirus viral DNA replication through the alternative splicing of NS1 mRNA to promote NS2 mRNA formation. *Vet Res.* 2021 May 25;52(1):73.
 21. Xuehui Zhao, Binbin Wen, Chen Li, Li Liu, Xiude Chen, Dongmei Li, Ling Li, Xiling Fu. PpEBB1 directly binds to the GCC box-like element of auxin biosynthesis related genes. *Plant Sci.* 2021 May;306:110874.
 22. Chao Li, Xue-Jiao Liu, Yan Yan, Mohammad Shah Alam, Zhen Liu, Zhen-Kun Yang, Ruo-Fu Tao, Er-Kui Yue, Ming-Hua Duan, Jian-Hong Xu. OsLHY is involved in regulating flowering through the Hd1- and Ehd1- mediated pathways in rice (*Oryza sativa* L.). *Plant Sci.* 2022 Feb;315:111145.
 23. Yuan Gao, Juan Xiong, Qing-Zhu Chu, Wei-Ke Ji. PDZD8-mediated lipid transfer at contacts between the ER and late endosomes/lysosomes is required for neurite outgrowth. *J Cell Sci.* 2022 Mar 1;135(5):jcs255026.
 24. Yan Hu, Ya-Ping Liang, Hang-Yu Tian, Cui-Xian Xu, Dan Yu, Pan Zhang, Hui Ye, Ming Li. Microplitis bicoloratus bracovirus regulates cyclophilin A-apoptosis-inducing factor interaction to induce cell apoptosis in the insect immunosuppressive process. *Arch Insect Biochem Physiol.* 2022 May;110(1):e21877.
 25. Guangyong Yang, Jie Yang, Qiwei Zhang, Wenfeng Wang, Liping Feng, Li Zhao, Bang An, Qiannan Wang, Chaozu He, Hongli Luo. The Effector Protein CgNLP1 of *Colletotrichum gloeosporioides* Affects Invasion and Disrupts Nuclear Localization of Necrosis-Induced Transcription Factor HbMYB8-Like to Suppress Plant Defense Signaling. *Front Microbiol.* 2022 Jun 13;13:911479.
 26. Yang Wu, Yichun Wang, Hanhui Yao, Heng Li, Fanzheng Meng, Qidong Li, Xiansheng Lin, Lianxin Liu. MNX1-AS1, a c-Myc induced lncRNA, promotes the Warburg effect by regulating PKM2 nuclear translocation. *J Exp Clin Cancer Res.* 2022 Dec 7;41(1):337.
 27. Qiuchen Li, Zhimin Jiang, Shuning Ren, Hui Guo, Zhimin Song, Saini Chen, Xintao Gao, Fanfeng Meng, Junda Zhu, Litao Liu, Qi Tong, Honglei Sun, Yipeng Sun, Juan Pu, Kin-Chow Chang, Jinhua Liu. SRSF5-Mediated Alternative Splicing of M Gene is Essential for Influenza A Virus Replication: A Host-Directed Target Against Influenza Virus. *Adv Sci (Weinh).* 2022 Dec;9(34):e2203088.
 28. Qianghua Zhou, Xu Chen, Kai Yao, Yangjie Zhang, Haixia He, Hao Huang, Hao Chen, Shengmeng Peng, Ming Huang, Liang Cheng, Qiang Zhang, Ruihui Xie, Kaiwen Li, Tianxin Lin, Hai Huang. TSPAN18 facilitates bone metastasis of prostate cancer by protecting STIM1 from TRIM32-mediated ubiquitination. *J Exp Clin Cancer Res.* 2023 Aug 5;42(1):195.

Version 2023.12.08