碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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- Xing Li, Manyu Chen, Peilei Wang, Ya Yao, Xiaowen Han, Jie Liang, Qing Jiang, Yong Sun, Yujiang Fan, Xingdong Zhang . A highly interweaved HA-SS-nHAp/collagen hybrid fibering hydrogel enhances osteoinductivity and mineralization Nanoscale. 2020 Jun 28;12(24):12869-12882.
- Fang Geng, Yunxiao Xie, Yi Wang, Jinqiu Wang . Depolymerization of chicken egg yolk granules induced by high-intensity ultrasound Food Chem. 2021 Aug 30;354:129580.
- Shuai Zhang, Shuwei Hu, Wanting Dong, Songqiang Huang, Zhexiao Jiao, Zewen Hu, Shiyun Dai, Yiwen Yi, Xiaohan Gong, Ke Li, Hui Wang, Dan Xu . Prenatal dexamethasone exposure induces anxiety- and depressive-like behavior of male offspring rats through intrauterine programming of the activation of NRG1-ErbB4 signaling in hippocampal PV intern Cell Biol Toxicol. 2021 Jun 29.
- Lulu Zhao, Shengnan Wang, Hong Liu, Xiaoli Du, Ren Bu, Bing Li, Ruilan Han, Jie Gao, Yang Liu, Jian Hao, Jianrong Zhao, Yan Meng, Gang Li . The Pharmacological Effect and Mechanism of Lanthanum Hydroxide on Vascular Calcification Caused by Chronic Renal Failure Hyperphosphatemia Front Cell Dev Biol. 2021 Apr 13;9:639127.
- Yunyun Cheng, Ting Chen, Jie Song, Qien Qi, Chunli Wang, Qianyun Xi, Songcai Liu, Linlin Hao, Yongliang Zhang . miR-709 inhibits GHRP6 induced GH synthesis by targeting PRKCA in pituitary Mol Cell Endocrinol. 2020 Apr 15;506:110763.
- Jiixin Li, Qian Zhou, Xin Zhou, Baodong Wei, Yingbo Zhao, Shujuan Ji . Calcium Treatment Alleviates Pericarp Browning of 'Nanguo' Pears by Regulating the GABA Shunt After Cold Storage Front Plant Sci. 2020 Sep 15;11:580986.
- Lei Cui, Qi Zhou, Xiufeng Zheng, Bowen Sun, Shilei Zhao . Mitoquinone attenuates vascular calcification by suppressing oxidative stress and reducing apoptosis of vascular smooth muscle cells via the Keap1/Nrf2 pathway Free Radic Biol Med. 2020 Dec;161:23-31.
- Hongzhou Shen, Jun Shi, Yin Zhi, Xiaoyan Yang, Yuan Yuan, Jiawen Si, Steve G F Shen . Improved BMP2-CPC-stimulated osteogenesis in vitro and in vivo via modulation of macrophage polarization Mater Sci Eng C Mater Biol Appl. 2021 Jan;118:111471.
- Songbin Wu, Shaomin Yang, Mingxi Ou, Jiamin Chen, Jiabing Huang, Donglin Xiong, Wuping Sun, Lizu Xiao . Transcriptome Analysis Reveals the Role of Cellular Calcium Disorder in Varicella Zoster Virus-Induced Post-Herpetic Neuralgia Front Mol Neurosci. 2021 May 17;14:665931.
- Zichao Xiang, Qionghui Wu, Yu Wang, Peng Wang, Yingyou He, Jihua Li . eIF2 $\alpha$ -ATF4 Pathway Activated by a Change in the Calcium Environment Participates in BCP-Mediated Bone Regeneration ACS Biomater Sci Eng. 2021 Jul 12;7(7):3256-3268.

11. Lanlan Li, Wenting Zhou, Sasa Miao, Xinyang Dong, Xiaoting Zou . Effects of Zinc on Cell Proliferation, Zinc Transport, and Calcium Deposition in Primary Endometrial Epithelial Cells of Laying Hens In Vitro Biol Trace Elem Res. 2021 Nov;199(11):4251-4259.
12. Xintao Wang, Bonan Tong, Rongji Hui, Congcong Hou, Zilu Zhang, Ludi Zhang, Bing Xie, Zhiyu Ni, Bin Cong, Chunling Ma, Di Wen . The Role of Hyperthermia in Methamphetamine-Induced Depression-Like Behaviors: Protective Effects of Coral Calcium Hydride Front Mol Neurosci. 2022 Jan 4;14:808807.

Version 2024.03.12