



碧云天生物技术/Beyotime Biotechnology
订热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

U0126 (MEK1/2抑制剂)

产品编号	产品名称	包装
S1901-1mg	U0126 (MEK1/2抑制剂)	10mg/ml×0.1ml
S1901-5mg	U0126 (MEK1/2抑制剂)	5mg
S1901-25mg	U0126 (MEK1/2抑制剂)	25mg
S1901-100mg	U0126 (MEK1/2抑制剂)	100mg

产品简介:

- U0126即1,4-Diamino-2,3-dicyano-1,4-bis(o-aminophenylmercapto)butadiene, 是一种常用的MAPKK抑制剂, 即MAPK kinase抑制剂或MEK1/2抑制剂。U0126可以通透细胞, 选择性抑制MEK1/2, 从而抑制MAP kinase (Erk1/2即p44/42 MAPK)的磷酸化和激活。同时可以抑制AP-1依赖的基因转录激活。
- U0126分子量为380.50, 分子式为C₁₈H₁₆N₆S₂, CAS Number: 109511-58-2。本产品纯度大于98%。
- 本U0126为进口分装, 其中10mg/ml包装用DMSO配制, 共0.1ml。5mg, 25mg和100mg为粉末装。

包装清单:

产品编号	产品名称	包装
S1901-1mg	U0126 (MEK1/2抑制剂)	10mg/ml×0.1ml
S1901-5mg	U0126 (MEK1/2抑制剂)	5mg
S1901-25mg	U0126 (MEK1/2抑制剂)	25mg
S1901-100mg	U0126 (MEK1/2抑制剂)	100mg
—	说明书	1份

保存条件:

-20°C避光保存, 一年有效。

注意事项:

- 本产品对人体有害, 操作时请小心, 并注意有效防护以避免直接接触人体或吸入体内。
- 本U0126在4°C、冰浴等较低温度情况下会凝固而粘在离心管管底、管壁或管盖内, 可以20-25°C水浴温育片刻至全部融解后使用。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. U0126常见使用浓度范围为1-20μM。具体的最佳工作浓度请参考相关文献, 或根据实验目的, 以及所培养的特定细胞和组织, 通过实验进行摸索和优化。

使用本产品的文献:

1. Pu J, Peng G, Li L, Na H, Liu Y, Liu P. Palmitic acid acutely stimulates glucose uptake via activation of Akt and ERK1/2 in skeletal muscle cells. *J Lipid Res.* 2011 Jul;52(7):1319-27.
2. Liang X, So YH, Cui J, Ma K, Xu X, Zhao Y, Cai L, Li W The low-dose ionizing radiation stimulates cell proliferation via activation of the MAPK/ERK pathway in rat cultured mesenchymal stem cells. *J Radiat Res.* 2011;52(3):380-6.
3. Xie H, Hua C, Sun L, Zhao X, Fan H, Dou H, Sun L, Hou Y. 17β-estradiol induces CD40 expression in dendritic cells via MAPK signaling pathways in aminochrome maintenance protein 6-dependent manner. *Arthritis Rheum.* 2011 Aug;63(8):2425-35.
4. Fang S, Jin Y, Zheng H, Yan J, Cui Y, Bi H, Jia H, Zhang H, Wang Y, Na L, Gao X, Zhou H High glucose condition upregulated Txnip expression level in rat mesangial cells through ROS/MEK/MAPK pathway. *Mol Cell Biochem.* 2011 Jan;347(1-2):175-82.
5. Ma J, Zhang L, Han W, Shen T, Ma C, Liu Y, Nie X, Liu M, Ran Y, Zhu D. Activation of JNK/c-Jun is required for the proliferation, survival, and angiogenesis induced by EET in pulmonary artery endothelial cells. *J Lipid Res.* 2012 Jun;53(6):1093-105.
6. Zhu F, Shen F, Fan Y, Xie Y, Xia Y, Kong Y. Osteopontin increases the expression of β1, 4-Galactosyltransferase-I and promotes adhesion in human RL95-2 cells. *Glycoconj J.* 2012 Aug;29(5-6):347-56.
7. Zhao X, Liu L, Liu D, Fan H, Wang Y, Hu Y, Hou Y. Progesterone Enhances Immunoregulatory Activity of Human Mesenchymal Stem Cells Via PGE2 and IL-6. *Am J Reprod Immunol.* 2012 Oct; 68(4):290-300.
8. Zhou F, Chen S, Xiong J, Li Y, Qu L. Luteolin Reduces Zinc-Induced Tau Phosphorylation at Ser262/356 in an ROS-Dependent Manner in SH-SY5Y Cells. *Biol Trace Elem Res.* 2012 Nov;149(2):273-9.
9. Hu Q, Li B, Xu R, Chen D, Mu C, Fei E, Wang G. The protease Omi cleaves the mitogen-activated protein kinase kinase MEK1 to inhibit microglial activation. *Sci Signal.* 2012 Aug 21;5(238):ra61.
10. Zhou R, Yao X, Xu X, Wang G, Zhu Z, Chen J, Chen L, Shen X. Blockage of progesterone receptor effectively protects pancreatic islet beta cell viability. *Steroids.* 2013 Jul 1. pii: S0039-128X(13)00146-3.

11. Li L, Han M, Li S, Wang L, Xu Y. Cyclic Tensile Stress During Physiological Occlusal Force Enhances Osteogenic Differentiation of Human Periodontal Ligament Cells via ERK1/2-Elk1 MAPK Pathway. *DNA Cell Biol.* 2013 Oct;78(10):987-95.
12. You M, Dong G, Li F, Ma F, Ren J, Xu Y, Yue H, Tang R, Ren D, Hou Y. Ligation of CD180 inhibits IFN- α signaling in a Lyn-PI3K-BTK-dependent manner in B cells. *Cell Mol Immunol.* 2013 Feb;27(2):251-7.
13. Ye L, Shi Z, Liu H, Yang X, Wang K. GdCl₃ induced Hep G2 cell death through mitochondrial and external death pathways without significant elevation of ROS generation. *Biol Trace Elem Res.* 2013 Jan;151(1):148-55.
14. Li HJ, Guo LM, Yang LL, Zhou YC, Zhang YJ, Guo J, Xie XJ, Guo GZ. Electromagnetic-pulse-induced activation of p38 MAPK pathway and disruption of blood-retinal barrier. *Toxicol Lett.* 2013 Jun 20;220(1):35-43.
15. Wang S, Zhao Z, Bi Y, Sun L, Liu X, Liu W. Tyrosine 132 phosphorylation of influenza A virus M1 protein is crucial for virus replication by controlling the nuclear import of M1. *J Virol.* 2013 Jun;87(11):6182-91.
16. Zhang Z, Miao L, Lv C, Sun H, Wei S, Wang B, Huang C, Jiao B. Wentilactone B induces G2/M phase arrest and apoptosis via the Ras/Raf/MAPK signaling pathway in human hepatoma SMMC-7721 cells. *Cell Death Dis.* 2013 Jun 6;4:e657.
17. Liu P, Wang YY, Qi X, Gu Q, Geng M, Li J. Undecylprodigiosin induced apoptosis in p38^{delta} cancer cells is associated with its binding to ribosome. *PLoS One.* 2013 Jun 14;8(6):e65381.
18. Zhang X, Ding R, Zhou Y, Zhu R, Liu W, Jin L, Yao W, Gao X. Toll-like receptor 2 and Toll-like receptor 4-dependent activation of B cells by a polysaccharide from marine fungus *Phoma herbarum* YS4108. *PLoS One.* 2013;8(3):e60781.
19. Liu T, Yi W, Feng B, Zhou Z, Xiao G. IGF-1-induced enhancement of PRNP expression depends on the negative regulation of transcription factor FOXO3a. *PLoS One.* 2013 Aug 14;8(8):e71896.
20. Shen Y, Zhao Y, Zheng D, Chang X, Ju S, Guo L. Effects of orexin A on GLUT4 expression and lipid content via MAPK signaling in 3T3-L1 adipocytes. *J Steroid Biochem Mol Biol.* 2013 Nov;138:376-83.
21. Xue P, Wu X, Zhou L, Ma H, Wang Y, Liu Y, Ma J, Li Y. IGF1 promotes osteogenic differentiation of mesenchymal stem cells derived from rat bone marrow by increasing TAZ expression. *Biochem Biophys Res Commun.* 2013 Apr 5;433(2):226-31.
22. Shi R, Lin J, Guo Y, Gong YP. The MEK1/2 inhibitor U0126 reverses imatinib resistance through down-regulating activation of Lyn/ERK signaling pathway in imatinib-resistant K562R leukemia cells. *Pharmazie.* 2014 May;69(5):346-52.
23. Shen M, Liu Z, Li B, Teng Y, Zhang J, Tang Y, Sun SC, Liu H. Involvement of FoxO1 in the effects of follicle-stimulating hormone on inhibition of apoptosis in mouse granulosa cells. *Cell Death Dis.* 2014 Oct 16;5:e1475.
24. Song MY, Yu JZ, Zhao DM, Wei S, Liu Y, Hu YM, Zhao WC, Yang Y, Wu H. The time-dependent manner of sinusoidal electromagnetic fields on rat bone marrow mesenchymal stem cells proliferation, differentiation, and mineralization. *Cell Biochem Biophys.* 2014 May; 69(1):47-54.
25. Lin B, Xu Y, Zhang B, He Y, Yan Y, He MC. MEK inhibition reduces glial scar formation and promotes the recovery of sensorimotor function in rats following spinal cord injury. *Exp Ther Med.* 2014 Jan;7(1):66-72.
26. Lv S1, Yi PF, Shen HQ, Zhang LY, Dong HB, Wu SC, Xia F, Guo X, Wei XB, Fu BD. Ginsenoside Rh2-B1 stimulates cell proliferation and IFN- γ production by activating the p38 MAPK and ERK-dependent signaling pathways in CTLL-2 cells. *Immunopharmacol Immunotoxicol.* 2014 Feb;36(1):43-51.
27. Liu Y, Mei C, Liu H, Wang H, Zeng G, Lin J, Xu M. Modulation of cytokine expression in human macrophages by endocrine-disrupting chemical Bisphenol-A. *Biochem Biophys Res Commun.* 2014 Sep 5; 451(4):592-8.
28. Jiao S, Wu M, Ye F, Tang H, Xie X, Xie X. BikDDA, a mutant of Bik with longer half-life expression protein, can be a novel therapeutic gene for triple-negative breast cancer. *PLoS One.* 2014 Mar 17;9(3):e92172.
29. Wang Y, Liu D, Zhao H, Jiang H, Luo C, Wang M, Yin H. Cordyceps sinensis polysaccharide CPS-2 protects human mesangial cells from PDGF-BB-induced proliferation through the PDGF/ERK and TGF- β 1/Smad pathways. *Mol Cell Endocrinol.* 2014 Feb 15; 382(2):979-88.
30. Yuan H, Chen R, Wu L, Chen Q, Hu A, Zhang T, Wang Z, Zhu X. The Regulatory Mechanism of Neurogenesis by IGF-1 in Adult Mice. *Mol Neurobiol.* 2015 Apr;51(2):512-22.
31. Liu C, Li L, Ha M, Qi S, Duan P, Yang K. The PI3K/Akt and ERK pathways elevate thyroid hormone receptor β 1 and TRH receptor to decrease thyroid hormones after exposure to PCB153 and p,p'-DDE. *Chemosphere.* 2015 Jan;118:229-38.
32. Liu P, Kong F, Wang J, Lu Q, Xu H, Qi T, Meng J. Involvement of IGF-1 and MEOX2 in PI3K/Akt1/2 and ERK1/2 pathways mediated proliferation and differentiation of perivascular adipocytes. *Exp Cell Res.* 2015 Feb 1;331(1):82-96.
33. Gu YY, Yang M, Zhao M, Luo Q, Yang L, Peng H, Wang J, Huang SK, Zheng ZX, Yuan XH, Liu P, Huang CZ. The de-ubiquitinase UCHL1 promotes gastric cancer metastasis via the Akt and Erk1/2 pathways. *Tumour Biol.* 2015 Nov;36(11):8379-87.
34. Guan J, Du S, Lv T, Qu S, Fu Q, Yuan Y. Oxygen-glucose Deprivation Preconditioning Protects Neurons against Oxygen-glucose Deprivation/ reperfusion Induced Injury via Bone Morphogenetic Protein-7 Mediated ERK, p38 and Smad Signaling Pathways. *Clin Exp Pharmacol Physiol.* 2015 Sep 19;43(1):1440-681.
35. Xia Q, Hu Q, Wang H, Yang H, Gao F, Ren H, Chen D, Fu C, Zheng L, Zhen X, Ying Z, Wang G. Induction of COX-2-PGE2 synthesis by activation of the MAPK/ERK pathway contributes to neuronal death triggered by TDP-43-depleted microglia. *Cell Death Dis.* 2015 Mar 26;6:e1702.
36. Chang X, Zhang J, Huang C, Pang X, Luo Q, Zhang H, Zhang S. Sex-determining region Y-related high mobility group box (SOX)-2 is overexpressed in cervical squamous cell carcinoma and contributes cervical cancer cell migration and invasion in vitro. *Tumour Biol.* 2015 Sep;36(10):7725-33.
37. Li T, Tang ZH, Xu WS, Wu GS, Wang YF, Chang LL, Zhu H, Chen XP, Wang YT, Chen Y, Lu JJ. Platycodin D triggers autophagy through activation of extracellular signal-regulated kinase in hepatocellular carcinoma HepG2 cells. *Eur J Pharmacol.* 2015 Feb 15;749:81-8.
38. Wang Y, Zhao S, Wu Z, Feng Y, Zhao C, Zhang C. Oxytocin in the regulation of social behaviours in medial amygdala-lesioned mice via the inhibition of the extracellular signal-regulated kinase signalling pathway. *Clin Exp Pharmacol Physiol.* 2015 May;42(5):465-74.
39. Shi F, Len Y, Gong Y, Shi R, Yang X, Naren D, Yan T. Ribavirin Inhibits the Activity of mTOR/eIF4E, ERK/Mnk1/eIF4E Signaling Pathway and Synergizes with Tyrosine Kinase Inhibitor Imatinib to Impair Bcr-Abl Mediated Proliferation and Apoptosis in Ph+ Leukemia. *PLoS One.* 2015 Aug 28;10(8):e0136746.
40. Chen Y, Yan H, Song Z, Chen F, Wang H, Niu J, Shi X, Zhang D, Zhang N, Zhai Z, Zhong B, Cheng L, Qian T, Hao F. Downregulation of TNIP1 Expression Leads to Increased Proliferation of Human Keratinocytes and Severe Psoriasis-Like Conditions in an Imiquimod-Induced Mouse Model of Dermatitis. *PLoS One.* 2015 Jun 5;10(6):e0127957.
41. Fan X, Xiong H, Wei J, Gao X, Feng Y, Liu X, Zhang G, He QY, Xu J, Liu L. Cytoplasmic hnRNPk interacts with GSK3 β and is essential for the osteoclast differentiation. *Sci Rep.* 2015 Dec 7;5:17732.
42. Li M, Li S, Li Y. Liraglutide Promotes Cortical Neurite Outgrowth via the MEK-ERK Pathway. *Cell Mol Neurobiol.* 2015 Oct;35(7):987-93.
43. Shang L, Jia SS, Jiang HM, Wang H, Xu WH, Lv CJ. Simvastatin downregulates expression of TGF- β RII and inhibits proliferation of A549 cells via ERK. 2015 Jun;36(6):4819-24.
44. Li W, Li G, Zhang Y, Wei S, Song M, Wang W, Yuan X, Wu H, Yang Y. Role of P2 \times 7 receptor in the differentiation of bone marrow stromal cells into osteoblasts and adipocytes. *Exp Cell Res.* 2015 Dec 10; 339(2):367-79.
45. Li T, Xu XH, Tang ZH, Wang YF, Leung CH, Ma DL, Chen XP, Wang YT, Chen Y, Lu JJ. Platycodin D induces apoptosis and triggers ERK- and JNK-mediated autophagy in human hepatocellular carcinoma BEL-7402 cells. *Acta Pharmacol Sin.* 2015 Dec;36(12):1503-13.
46. Wen J, Zhao Y, Guo L. Orexin A induces autophagy in HCT-116 human colon cancer cells through the ERK signaling pathway. *Int J Mol Med.* 2016 Jan;37(1):126-32.
47. Xiang T, Fei R, Wang Z, Shen Z, Qian J, Chen W. Nicotine enhances invasion and metastasis of human colorectal cancer cells through the nicotinic acetylcholine receptor downstream p38 MAPK signaling pathway. *Oncol Rep.* 2016 Jan;35(1):205-10.
48. Yu M, Chen R, Jia Z, Chen J, Lou J, Tang S, Zhang X. MWCNTs Induce

- ROS Generation, ERK Phosphorylation, and SOD-2 Expression in Human Mesothelial Cells. *Int J Toxicol*. 2016 Jan-Feb;35(1):17-26.
49. Hu A, Yuan H, Wu L, Chen R, Chen Q, Zhang T, Wang Z, Liu P, Zhu X. The effect of constitutive over-expression of insulin-like growth factor 1 on the cognitive function in aged mice. *Brain Res*. 2016 Jan; 1631(15):204-13.
 50. Guan J, Du S, Lv T, Qu S, Fu Q, Yuan Y. Oxygen-glucose deprivation preconditioning protects neurons against oxygen-glucose deprivation/reperfusion induced injury via bone morphogenetic protein-7 mediated ERK, p38 and Smad signalling pathways. *CLIN EXP PHARMACOL P*. 2016 Jan;43(1):125-34.
 51. Yu M, Chen R, Jia Z, Chen J, Lou J, Tang S, Zhang X. MWCNTs Induce ROS Generation, ERK Phosphorylation, and SOD-2 Expression in Human Mesothelial Cells. *Int J Toxicol*. 2016 Jan-Feb;35(1):17-26.
 52. Xiang T, Fei R, Wang Z, Shen Z, Qian J, Chen W. Nicotine enhances invasion and metastasis of human colorectal cancer cells through the nicotinic acetylcholine receptor downstream p38 MAPK signaling pathway. *Oncol Rep*. 2016 Jan;35(1):205-10.
 53. Li Z, Wang S, Zhao W, Sun Z, Yan H, Zhu J. Oxidized low-density lipoprotein upregulates microRNA-146a via JNK and NF- κ B signaling. *Mol Med Rep*. 2016 Feb;13(2):1709-16.
 54. Han Y, Wu Z, Wu T, Huang Y, Cheng Z, Li X, Sun T, Xie X, Zhou Y, Du Z. Tumor-suppressive function of long noncoding RNA MALAT1 in glioma cells by downregulation of MMP2 and inactivation of ERK/MAPK signaling. *Cell Death Dis*. 2016 Mar 3;7:e2123.
 55. Gao M, Kong Q, Hua H, Yin Y, Wang J, Luo T, Jiang Y. AMPK-mediated up-regulation of mTORC2 and MCL-1 compromises the anti-cancer effects of aspirin. *ONCOTARGET*. 2016 Mar 29;7(13):16349-61.
 56. Lv S, Xu QY, Sun EC, Zhang JK, Wu DL. Dissection and integration of the autophagy signaling network initiated by bluetongue virus infection: crucial candidates ERK1/2, Akt and AMPK. *SCI REP-UK*. 2016 Mar 15;6:23130.
 57. Chen Q, Zhu L, Zheng B, Wang J, Song X, Zheng W, Wang L, Yang D, Wang J. Effect of AQP9 Expression in Androgen-Independent Prostate Cancer Cell PC3. *Int J Mol Sci*. 2016 May 14;17(5): pii: E738.
 58. Wang Y, Cao J, Fan Y, Xie Y, Xu Z, Yin Z, Gao L, Wang C. Artemisinin inhibits monocyte adhesion to HUVECs through the NF- κ B and MAPK pathways in vitro. *Int J Mol Med*. 2016 Jun;37(6):1567-75.
 59. Wu W, Sun Y, Zhao C, Zhao C, Chen X, Wang G, Pang W, Yang G. Lipogenesis in myoblasts and its regulation of CTRP6 by AdipoR1/Erk/PPAR γ signaling pathway. *ACTA BIOCH BIPH SIN*. 2016 Jun;48(6):509-19.
 60. Wang K, Zhang C, Bao J, Jia X, Liang Y, Wang X, Chen M, Su H, Li P, Wan JB, He C. Synergistic chemopreventive effects of curcumin and berberine on human breast cancer cells through induction of apoptosis and autophagic cell death. *SCI REP-UK*. 2016 Jun 6;6:26064.
 61. Wu D, Cheng J, Sun G, Wu S, Li M, Gao Z, Zhai S, Li P, Su D, Wang X. p70S6K promotes IL-6-induced epithelial-mesenchymal transition and metastasis of head and neck squamous cell carcinoma. *ONCOTARGET*. 2016 Jun 14;7(24):36539-36550.
 62. Li P, Xu Y, Gan Y, Song L, Zhang C, Wang L, Zhou Q. Role of the ERK1/2 Signaling Pathway in Osteogenesis of Rat Tendon-Derived Stem Cells in Normoxic and Hypoxic Cultures. *Int J Med Sci*. 2016 Jul 18;13(8):629-37.
 63. Wang J, Chen S, Liao Y, Zhang E, Feng S, Yu S, Li LF, He WR, Li Y, Luo Y, Sun Y, Zhou M, Wang X, Munir M, Li S, Qiu HJ. Mitogen-activated Protein Kinase Kinase 2 (MEK2), a Novel E2-interacting Protein, Promotes the Growth of Classical Swine Fever Virus via Attenuation of the JAK-STAT Signaling Pathway. *J Virol*. 2016 Sep 7. pii: JVI.01407-16.
 64. Ma K, Huang MY, Guo YX, Hu GQ. Matrine-induced autophagy counteracts cell apoptosis via the ERK signaling pathway in osteosarcoma cells. *Oncol Lett*. 2016 Sep;12(3):1854-1860.
 65. Ji M, Lu Y, Zhao C, Gao W, He F, Zhang J, Zhao D, Qiu W, Wang Y. C5a Induces the Synthesis of IL-6 and TNF- α in Rat Glomerular Mesangial Cells through MAPK Signaling Pathways. *PLoS One*. 2016 Sep 1;11(9):e0161867.
 66. Zhang M, Song S, Yi Z, Zhao X, Fu L, Wang L, Ma C, Mao M, Xing Y, Zhu D. Human biliverdin reductase promotes EMT through the ERK1/2 signal pathway in breast cancer. *Eur J Pharmacol*. 2016 Oct 5;788:45-53.
 67. Han Y, Zhou L, Wu T, Huang Y, Cheng Z, Li X, Sun T, Zhou Y, Du Z. Downregulation of lncRNA-MALAT1 Affects Proliferation and the Expression of Stemness Markers in Glioma Stem Cell Line SHG139S. *Cell Mol Neurobiol*. 2016 Oct;36(7):1097-107.
 68. Zhou M, Fu J, Xiao L, Yang S, Song Y, Zhang X, Feng X, Sun H, Xu W, Huang W. miR-196a overexpression activates the MEK/ERK signal and represses the progesterone receptor and decidualization in eutopic endometrium from women with endometriosis. *Hum Reprod*. 2016 Nov;31(11):2598-2608.
 69. Lin X, Wu S, Wang Q, Shi Y, Liu G, Zhi J, Wang F. Saikosaponin-D Reduces H2O2-Induced PC12 Cell Apoptosis by Removing ROS and Blocking MAPK-Dependent Oxidative Damage. *Cell Mol Neurobiol*. 2016 Nov;36(8):1365-1375.
 70. Liu Q, Li H, Wang Q, Zhang Y, Wang W, Dou S, Xiao W. Increased expression of TROP2 in airway basal cells potentially contributes to airway remodeling in chronic obstructive pulmonary disease. *RESP RES*. 2016 Nov 25;17(1):159.
 71. Shen J, Jiang F, Yang Y, Huang G, Pu F, Liu Q, Chen L, Ju L, Lu M, Zhou F, Zhang C, Luo X, Yang X, Jiao C, Li X, Li Z, Li Y, Zhang J. 14-3-3 η is a novel growth-promoting and angiogenic factor in hepatocellular carcinoma. *J Hepatol*. 2016 Nov;65(5):953-962.
 72. Jiang Z, Hua Y. Hydrogen sulfide promotes osteogenic differentiation of human periodontal ligament cells via p38-MAPK signaling pathway under proper tension stimulation. *Arch Oral Biol*. 2016 Dec;72:8-13.
 73. Liu X, Ren K, Suo R, Xiong SL, Zhang QH, Mo ZC, Tang ZL, Jiang Y, Peng XS, Yi GH. ApoA-I induces S1P release from endothelial cells through ABCA1 and SR-BI in a positive feedback manner. *J Physiol Biochem*. 2016 Dec;72(4):657-667.
 74. Yu HP, Zhang N, Zhang T, Wang ZL, Li N, Tang HH, Zhang R, Zhang MN, Xu B, Fang Q, Wang R. Activation of NPFF2 receptor stimulates neurite outgrowth in Neuro 2A cells through activation of ERK signaling pathway. *Peptides*. 2016 Dec;86:24-32.
 75. He LX, Tong X, Zeng J, Tu Y, Wu S, Li M, Deng H, Zhu M, Li X, Nie H, Yang L, Huang F. Paeonol Suppresses Neuroinflammatory Responses in LPS-Activated Microglia Cells. *Inflammation*. 2016 Dec;39(6):1904-1917.
 76. Fan B, Wang BF, Che L, Sun YJ, Liu SY, Li GY. Blockage of NOX2/MAPK/NF- κ B Pathway Protects Photoreceptors against Glucose Deprivation-Induced Cell Death. *Oxid Med Cell Longev*. 2017;2017:5093473.
 77. Qi R, Liu H, Wang Q, Wang J, Yang F, Long D, Huang J. Expressions and Regulatory Effects of P38/ERK/JNK Maps in the Adipogenic Trans-Differentiation of C2C12 Myoblasts. *CELL PHYSIOL BIOCHEM*. 2017;44(6):2467-2475.
 78. Yang M, Wang X, Wang L, Wang X, Li J, Yang Z. IL-1 α Up-Regulates IL-6 Expression in Bovine Granulosa Cells via MAPKs and NF- κ B Signaling Pathways. *CELL PHYSIOL BIOCHEM*. 2017;41(1):265-273.
 79. Geng H, Zhao L, Liang Z, Zhang Z, Xie D, Bi L, Wang Y, Zhang T, Cheng L, Yu D, Zhong C. Cigarette smoke extract-induced proliferation of normal human urothelial cells via the MAPK/AP-1 pathway. *Oncol Lett*. 2017 Jan;13(1):469-475.
 80. Xu Y, Zhao Y, Xu Y, Guan Y, Zhang X, Chen Y, Wu Q, Zhu G, Chen Y, Sun F, Wang J, Yu Y. Blocking inhibition to YAP by Actinomycin D enhances anti-tumor efficacy of Corosolic acid in treating liver cancer. *Cell Signal*. 2017 Jan;29:209-217.
 81. Yu K, Ma Y, Li X, Wu X, Liu W, Li X, Shen J, Wang H. Lipopolysaccharide increases IL-6 secretion via activation of the ERK1/2 signaling pathway to up-regulate RANKL gene expression in MLO-Y4 cells. *Cell Biol Int*. 2017 Jan;41(1):84-92.
 82. Xu Y, Zhao Y, Xu Y, Guan Y, Zhang X, Chen Y, Wu Q, Zhu G, Chen Y, Sun F, Wang J, Yu Y. Blocking inhibition to YAP by Actinomycin D enhances anti-tumor efficacy of Corosolic acid in treating liver cancer. *Cell Signal*. 2017 Jan;29:209-217.
 83. Yu K, Ma Y, Li X, Wu X, Liu W, Li X, Shen J, Wang H. Lipopolysaccharide increases IL-6 secretion via activation of the ERK1/2 signaling pathway to up-regulate RANKL gene expression in MLO-Y4 cells. *Cell Biol Int*. 2017 Jan;41(1):84-92.
 84. Yu D, Geng H, Liu Z, Zhao L, Liang Z, Zhang Z, Xie D, Wang Y, Zhang T, Min J, Zhong C. Cigarette smoke induced urocytic epithelial mesenchymal transition via MAPK pathways. *ONCOTARGET*. 2017 Jan 31;8(5):8791-8800.
 85. You N, Li J, Gong Z, Huang X, Wang W, Wang L, Wu K, Zheng L. COMMD7 functions as molecular target in pancreatic ductal adenocarcinoma. *MOL CARCINOGEN*. 2017 Feb;56(2):607-624.
 86. Wang H, Li C, Li J, Zhu Y, Jia Y, Zhang Y, Zhang X, Li W, Cui L, Li W, Liu Y. Naringin enhances osteogenic differentiation through the

- activation of ERK signaling in humanbone marrow mesenchymal stem cells. *Iran J Basic Med Sci* . 2017 Apr;20(4):408-414.
87. Zhao K, Zhang S, Song X, Yao Y, Zhou Y, You Q, Guo Q, Lu N. Gambogic acid suppresses cancer invasion and migration by inhibiting TGF β 1-induced epithelial-to-mesenchymal transition. *ONCOTARGET* . 2017 Apr 18;8(16):27120-27136.
 88. Li Y, Wu Z, Yuan J, Sun L, Lin L, Huang N, Bin J, Liao Y, Liao W. Long non-coding RNA MALAT1 promotes gastric cancer tumorigenicity and metastasis by regulating vasculogenic mimicry and angiogenesis. *Cancer Lett* . 2017 Jun 1;395:31-44.
 89. Dai XQ, Cai WT, Wu X, Chen Y, Han FM. Protocatechuic acid inhibits hepatitis B virus replication by activating ERK1/2 pathway and down-regulating HNF4 α and HNF1 α in vitro. *Life Sci* . 2017 Jul 1;180:68-74.
 90. Huang WS, Xu FM, Zeng QZ, Liu XH, Gao XJ, Liu LX. ERK1/2-mediated Cytoplasmic Accumulation of hnRNPK Antagonizes TRAIL-induced Apoptosis through Upregulation of XIAP in H1299 Cells. *Biomed Environ Sci* . 2017 Jul;30(7):473-481.
 91. Zhang C, Yu P, Zhu L, Zhao Q, Lu X, Bo S. Blockade of α 7 nicotinic acetylcholine receptors inhibit nicotine-induced tumor growth and vimentin expression in non-small cell lung cancer through MEK/ERK signaling way. *Oncol Rep* . 2017 Dec;38(6):3309-3318.
 92. Zhang M, Xin W, Yi Z, Li Y, Liu Y, Zhang H, Chen H, Chen X, Tan S, Zhu D. Human biliverdin reductase regulates the molecular mechanism underlying cancer development. *J Cell Biochem* . 2018 Feb;119(2):1337-1345.
 93. Liu Q, Song J, Li H, Dong L, Dai S. Schizandrin B inhibits the cis DDP induced apoptosis of HK 2 cells by activating ERK/NF κ B signaling to regulate the expression of survivin. *Int J Mol Med* . 2018 Apr;41(4):2108-2116.
 94. Wang T, Cheng J, Wang S, Wang X, Jiang H, Yang Y, Wang Y, Zhang C, Liang W, Feng H. α -Lipoic acid attenuates oxidative stress and neurotoxicity via the ERK/Akt-dependent pathway in the mutant hSOD1 related Drosophila model and the NSC34 cell line of amyotrophic lateral sclerosis. *Brain Res Bull* . 2018 Jun;140:299-310.
 95. Geng B, Liang M, Qin L, Zhao W, Wang H, Wang L, Pan X, Chen X. An TRIM59-CDK6 axis regulates growth and metastasis of lung cancer. *J Cell Mol Med*. 2019 Feb;23(2):1458-1469.
 96. Lin Y, Cong H, Liu K, Jiao Y, Yuan Y, Tang G, Chen Y, Zheng Y, Xiao J, Li C, Chen Z, Cao P. Microbicidal Phagocytosis of Nucleus Pulposus Cells Against Staphylococcus aureus via the TLR2/MAPKs Signaling Pathway. *Front Immunol*. 2019 May 22;10:1132.
 97. Wang H, Ning T, Song C, Luo X, Xu S, Zhang X, Deng Z, Ma D, Wu B. Priming integrin α 5 promotes human dental pulp stem cells odontogenic differentiation due to extracellular matrix deposition and amplified extracellular matrix-receptor activity. *J Cell Physiol*. 2019 Aug;234(8):12897-12909.
 98. Zhang C, Xie H, Yang Q, Yang Y, Li W, Tian H, Lu L, Wang F, Xu JY, Gao F, Wang J, Jin C, Xu G, Xu GT, Zhang J. Erythropoietin protects outer blood-retinal barrier in experimental diabetic retinopathy by up-regulating ZO-1 and occludin. *Clin Exp Ophthalmol*. 2019 Dec;47(9):1182-1197.

Version 2021.12.02