## IGEPAL® CA-630 (Nonidet P-40 Substitute, Reagent grade)


| 产品编号         | 产品名称  | 包装    |
|--------------|---|-------|
| ST2045-100ml | IGEPAL® CA-630 (Nonidet P-40 Substitute, Reagent grade) | 100ml |
| ST2045-500ml | IGEPAL® CA-630 (Nonidet P-40 Substitute, Reagent grade) | 500ml |

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| CAS Number | Chemical Formula            | Molecular Weight | Purity | Grade         |
|------------|-----------------------------|------------------|--------|---------------|
| 9002-93-1  | $(C_2H_4O)_n C_{14}H_{22}O$ | ~603             | -      | Reagent grade |

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|  |  |
|--|--|
|  | <p>intracellular tubulin polymerization assays for screening of microtubule targeting agents' (Biochem. Cell Biol. 95: 379–384 (2017)), IGEPAL® CA-630是Nonidet P-40最适合的替代物; 同时, 如果需要和Nonidet P-40达到一样的效果, 本产品的使用浓度可能需要根据实验进行一定的调整。</p> <p>IGEPAL, with an HLB number of 13, is the most similar substitute to the Shell Nonidet P-40. Tergitol and the other two detergents consist of nonyl-phenyl polyethylene glycol with an HLB number of 17.8. The HLB number relates to the lipid solubility of the compound, with a lower HLB indicating the more lipophilic substance (HLB = 0 indicating 100% lipophilic and HLB = 20 indicating 100% hydrophilic). All of the 4 substitute detergents, at concentrations around 0.1%–0.15%, demonstrated good, consistent preservation of polymerized tubulin in the hypotonic lysate before separation into soluble and polymerized tubulin fractions.</p> <p>The results of this study indicate that the substitute detergents are approximately 10-fold more potent than the original Shell Nonidet P-40 detergent that was generally used at a 1% working concentration. Given the importance and ubiquitous use of the Shell Nonidet P-40 in various biochemical buffers and assays, it was concluded that all of the laboratories using the new, substitute detergents in their research should in the first instance reduce the concentration of the buffer by about 10-fold to obtain the same result as with the discontinued Shell Nonidet P-40.</p> |
|--|--|

➤ **性质(Properties):**

|                           |  |
|---------------------------|--|
| Description               | non-ionic                              |
| Appearance (Form)         | viscous liquid                         |
| Appearance (Turbidity)    | clear to slightly hazy                 |
| Appearance (Color)        | colorless to faint yellow              |
| mol wt                    | ~603                                   |
| CMC                       | 0.08 mM (20-25°C)                      |
| Viscosity                 | 250-350 cP, neat (25 °C)(lit.)         |
| Brookfield Viscosity      | 1.06 g/cm <sup>3</sup> at 25 °C (lit.) |
| HLB                       | 13                                     |
| Water (by Karl Fischer)   | ≤0.5 %                                 |
| Cloud Point (1% Solution) | 63.0 - 67.0 °C                         |
| Neutralization Number     | ≤0.5 mg KOH/g (pH7)                    |

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