



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

SOD (抗氧化酶)

产品编号	产品名称	包装
S0086	SOD (抗氧化酶)	3KU

产品简介:

- SOD即Superoxide dismutase, 中文名为超氧化物歧化酶。SOD可以催化超氧自由基转化为过氧化氢和分子氧, 在抵御氧自由基导致的细胞损伤中起关键作用。SOD可以抑制一些细胞的凋亡。有报道SOD可以促进一氧化氮的作用。
- 本SOD为进口分装, 纯化自牛红细胞, 酶活力为2500-7000UN/mg蛋白。一个活力单位的SOD, 在25°C, pH7.8, xanthine oxidase 耦联体系存在的情况下, 可以抑制细胞色素C还原的50%。
- SOD分子量31.2kD, 纯度>99%。

包装清单:

产品编号	产品名称	包装
S0086	SOD (抗氧化酶)	3KU/100微升
—	说明书	1份

保存条件:

-20°C保存, 一年有效。

注意事项:

- 尽量避免多次反复冻融, 宜适当分装后-20°C保存。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用本产品的文献:

1. Hu Y, Cai K, Luo Z, Xu D, Xie D, Huang Y, Yang W, Liu P. . TiO₂ nanotubes as drug nanoreservoirs for the regulation of mobility and differentiation of mesenchymal stem cells. *Acta Biomater.* 2012 Jan;8(1):439-48.
2. Zhang Y, Xu Z, Ding J, Tan C, Hu W, Li Y, Huang W, Xu Y. . HZ08 suppresses RelB-activated MnSOD expression and enhances Radiosensitivity of prostate Cancer cells. *J EXP CLIN CANC RES.* 2018 Jul 27;37(1):174.
3. Zheng X, Wu J, Shao Q, Li X, Kou J, Zhu X, Zhong Z, Jiang Y, Liu Z, Li H, Tian Y, Yang L. . Apoptosis of THP-1 Macrophages Induced by Pseudohypericin-Mediated Sonodynamic Therapy Through the Mitochondria-Caspase Pathway. *CELL PHYSIOL BIOCHEM.* 2016;38(2):545-57.
4. Zhang X, Wang YN, Zhu JJ, Liu XX, You H, Gong MY, Zou M, Cheng WH, Zhu JH. . N-acetylcysteine negatively regulates Notch3 and its malignant signaling. *ONCOTARGET.* 2016 May 24;7(21):30855-66.
5. Xin Zhang, Nianfeng Zhang, Juan Kan, Rui Sun, Sixue Tang, Zhihao Wang, Mengfei Chen, Jun Liu, Changhai Jin . Anti-inflammatory activity of alkali-soluble polysaccharides from *Arctium lappa* L. and its effect on gut microbiota of mice with inflammation *Int J Biol Macromol.* 2020 Jul 1;154:773-787.
6. Xie Y, Liu D, Cai C, Chen X, Zhou Y, Wu L, Sun Y, Dai H, Kong X, Liu P. . Size-dependent cytotoxicity of Fe₃O₄ nanoparticles induced by biphasic regulation of oxidative stress indifferent human hepatoma cells. *INT J NANOMED.* 2016 Jul 29;11:3557-70.
7. Li NS, Luo XJ, Dai Z, Liu B, Zhang YS, Yang ZC, Peng J. . Beneficial effects of capsate on ethanol-induced mucosal injury in rats are related to stimulation of calcitonin gene-related Peptide release. *Planta Med.* 2012 Jan;78(1):24-30.
8. Zhang X, Song Y, Han X, Feng L, Wang R, Zhang M, Zhu M, Jia X, Hu S. . Liquiritin attenuates advanced glycation end products-induced endothelial dysfunction via RAGE/NF-κB pathway in human umbilical vein endothelial cells. *Mol Cell Biochem.* 2013 Feb;374(1-2):191-201.
9. Kong X, Ma MZ, Qin L, Zhang Y, Li XY, Wang GD, Su Q, Zhang DY. . Pioglitazone enhances the blood pressure-lowering effect of losartan via synergistic attenuation of angiotensin II-induced vasoconstriction. *J RENIN-ANGIO-ALDO S.* 2014 Sep;15(3):259-70.
10. Cai G, Yang X, Lai Q, Yu X, Zhang H, Li Y, Chen Z, Lei X, Zheng W, Xu H, Zheng T. . Lysing bloom-causing alga *Phaeocystis globosa* with microbial algicide: An efficient process that decreases the toxicity of algal exudates. *SCI REP-UK.* 2016 Feb 5;6:20081.
11. Zhang X, Wang L, Wang R, Luo X, Li Y, Chen Z. . Protective effects of rice dreg protein hydrolysates against hydrogen peroxide-induced oxidative stress in HepG-2 cells. *Food Funct.* 2016 Mar;7(3):1429-37.
12. Wang J, Wang G, Sun GW. . Role of PPARα in down-regulating AGE-induced TGF-β and MMP-9 expressions in chondrocytes. *GENET MOL RES.* 2016 May 9;15(2).
13. Fu P, Hu Q. . 3,4-Dihydroxyphenylethanol alleviates early brain injury by modulating oxidative stress and Akt and nuclear factor-κB pathways in a rat

- model of subarachnoid hemorrhage. *Exp Ther Med*. 2016 May;11(5):1999-2004.
14. Yu J, Liu C, Li Z, Zhang C, Wang Z, Liu X. . Inhibitory effects and mechanism of 25-OH-PPD on glomerular mesangial cell proliferation induced by highglucose. *ENVIRON TOXICOL CHEM*. 2016 Jun;44:93-8.
 15. Jiang Q, Hao R, Wang W, Gao H, Wang C. . SIRT1/Atg5/autophagy are involved in the antiatherosclerosis effects of ursolic acid. *Mol Cell Biochem*. 2016 Sep;420(1-2):171-84.
 16. Gong W, Chen C, Xiong F, Yang Z, Wang Y, Huang J, Liu P, Huang H. . CKIP-1 ameliorates high glucose-induced expression of fibronectin and intercellular cell adhesion molecule-1 by activating the Nrf2/ARE pathway in glomerular mesangial cells. *Biochem Pharmacol*. 2016 Sep 15;116:140-52.
 17. Liu Y, Zhi D, Li M, Liu D, Wang X, Wu Z, Zhang Z, Fei D, Li Y, Zhu H, Xie Q, Yang H, Li H. . Shengmai Formula suppressed over-activated Ras/MAPK pathway in *C. elegans* by opening mitochondrialpermeability transition pore via regulating cyclophilin D. *SCI REP-UK*. 2016 Dec 16;6:38934.
 18. Zhang W, Li X, Liu Y, Chen H, Gong J. . Activation of imidazoline II receptor by moxonidine regulates the progression of liver fibrosis in the Nrf2-dependent pathway. *Biomed Pharmacother*. 2017 Jun;90:821-834.
 19. Shi H, Wang Q, Yang L, Xie S, Zhu H. . IMM-H007, a new therapeutic candidate for nonalcoholic fatty liver disease, improves hepaticsteatosis in hamsters fed a high-fat diet. *FEBS Open Bio*. 2017 Aug 29;7(9):1379-1391.
 20. Qiao H, Ren H, Du H, Zhang M, Xiong X, Lv R. . Liraglutide repairs the infarcted heart: The role of the SIRT1/Parkin/mitophagy pathway. *Mol Med Rep*. 2018 Mar;17(3):3722-3734.
 21. Li Y, Kong X, Chen J, Liu H, Zhang H. . Characteristics of the Copper,Zinc Superoxide Dismutase of a Hadal Sea Cucumber (*Paelopatides* sp.) from the Mariana Trench. *Mar Drugs*. 2018 May 18;16(5). pii: E169.
 22. Gao LX, Bian C, Wu Y, Nisar MF, Chen S, Li CM, Yu L, Ji P, Huang E, Zhong JL. . Label-free electrochemical sensor to investigate the effect of tocopherol on generation of superoxide ions following UV irradiation. *J Biol Eng*. 2018 Sep 12;12:17.
 23. Yu C, Wang W, Jin X. . Hirudin Protects Ang II-Induced Myocardial Fibroblasts Fibrosis by Inhibiting the Extracellular Signal-Regulated Kinase1/2 (ERK1/2) Pathway. *MED SCI MONITOR*. 2018 Sep 8;24:6264-6272.
 24. Ding X, Sun W, Chen J. . IL-2 augments the sorafenib-induced apoptosis in liver cancer by promoting mitochondrial fission and activating the JNK/TAZ pathway. *Cancer Cell Int*. 2018 Nov 9;18:176.
 25. Wu Y, Li Y, Jiang T, Yuan Y, Li R, Xu Z, Zhong X, Jia G, Liu Y, Xie L, Xu K, Zhang H, Li X, Xiao J. . Reduction of cellular stress is essential for Fibroblast growth factor 1 treatment for diabetic nephropathy. *J Cell Mol Med*. 2018 Dec;22(12):6294-6303.
 26. Zhang L, Li S, Wang R, Chen C, Ma W, Cai H. . Cytokine augments the sorafenib-induced apoptosis in Huh7 liver cancer cell by inducing mitochondrial fragmentation and activating MAPK-JNK signalling pathway. *Biomed Pharmacother*. 2019 Feb 110:213-223.
 27. Haijun Mao, Guangyue Xu. . Protective effect and mechanism of microRNA-146a on ankle fracture *Exp Ther Med*. 2020 Nov;20(5):3.
 28. Li-Zhe Xu, Kai-Xiang He, Jin-Zhuo Ning, Fan Cheng. . Oleuropein attenuates testicular ischemia-reperfusion by inhibiting apoptosis and inflammation *Tissue Cell*. 2022 Oct;78:101876.

Version 2024.03.12